


## Custom product key features

The complete range of data sheets is available on [www.phonak.com](http://www.phonak.com)


Custom products	Power		Power		Power	
	Certéna CIC/MC Petite	Certéna CIC/MC P Petite	Certéna ITC/HS Petite	Certéna ITC/HS P Petite	Certéna FS	Certéna FS P
Performance profile						
Maximum Power Output (dB SPL)						
2 cc coupler	113	115	113	115	115	122
Ear Simulator	123	125	123	126	125	130
Maximum gain (dB)						
2 cc coupler	35	55	48	55	50	55
Ear Simulator	45	65	59	66	60	63
Frequency range – 2 cc coupler (Hz)	<100 – 7000	<100 – 6200	<100 – 6600	<100 – 6700	<100 – 6400	<100 – 5500
Battery size	10	10	312	312	13	13
Working current (mA)	0.9	1.1	1.2	1.1	1.2	1.2

## Fitting range



## Technical description

- 1 Microphone input (2 inputs on directional models)
- 2 TacTronic Switch
- 3 Battery compartment and programming socket
- 4 Volume control
- 5 Sound outlet
- 6 Custom made shell
- 7 Venting



## Product information



Certéna, based on the proven wireless chip technology platform CORE, offers a complete product line including a wide range of custom products as well as micro Petite BTE, micro BTE and traditional BTE, accommodating mild to profound hearing losses. All products offer the sophisticated WhistleBlockTechnology that suppresses feedback successfully without gain loss or artifacts.

The automatic SoundFlow Standard instantly optimizes the hearing instrument setting to the changing environment. Certéna offers a wide range of accessories and FM possibilities, giving the convenience of wireless connectivity, or the ease of a remote control etc.

## Key features\*

### Features

- Fully programmable 6 channel digital hearing instrument
- SoundFlow Standard – a completely new approach to self-adapting automatic functionality based on 5 fine tunable configuration points
- Digital Audio Zoom – Single channel adaptive directionality
- WhistleBlock Technology – the new benchmark in feedback elimination
- Wind Et Weather Protector – protects the microphones
- NoiseBlock Processing – 6 channel SNR dependent noise canceller
- AOV – Acoustically Optimized Vent – delivers the right sized custom vent solution for all Certéna custom products and ear pieces
- EasyPhone – automatic activation of telephone program
- EasyFM – automatic activation of FM program
- EasyAudio – automatic program for binaural stereo streaming of audio signal
- EasyBluetooth – automatic program for handsfree binaural mobile phone use
- Acoustic confirmation of manual interventions and end of battery life
- DataLogging with SmartTips
- Up to four manually accessible programs (two acoustic, two non-acoustic)
- BassBoost
- QuickSync – binaural volume and program control with one-touch simplicity

### Connectivity

- iCom – easy, binaural handsfree connection to communication and entertainment devices
- CableFree Fitting – completely cableless, smooth and reliable fitting with iCube

### Accessories

- myPilot – a command center with full-color display and unique access to hearing and instrument control and status readout
- KeyPilot2 / WatchPilot2 remote controls
- iView – status viewer
- Click'nTalk – easy to use, intuitive, Plug-in device for the mobile phone. Completely cable and headset free mobile phone use.

### Fitting

- iPFG version 2.4 or higher, NOAH compatible
- PC (IBM compatible) with NOAHlink and/or HI-PRO
- Optional: iCube device for CableFree Fitting

\*please note that some features are model dependent.



PHONAK life is on

## Programming

### CableFree Fitting with iCube\*



1. Switch iCube on and check the Bluetooth connection to the computer.
2. Put iCube around the patient's neck and connect the neckloop to the iCube.

\*Certéna Petite instruments are not compatible with CableFree Fitting.

### Programming a Certéna BTE via cable



1. Flip open the cover for the programming socket.
2. Connect the programming plug into the socket. Use programming cable CS-44 A (Phonak version) only.

### Programming a Certéna micro and micro Petite via cable



1. Open the battery compartment.
2. Connect the programming plug into the socket. Use programming cable CS-44 A (Phonak version) only.

### Programming a Certéna custom product via cable




1. Open the battery compartment and lift the cover out of the hinge, using the battery compartment tool (twisting motion).
2. Connect the cable. Attach the cables with the blue and red markers to the left and the right hearing instruments, respectively.

## BTE key features


The complete range of data sheets is available on [www.phonak.com](http://www.phonak.com)

BTE	Power		SuperPower		
	Certéna micro Petite	Certéna micro	Certéna M	Certéna P	Certéna SP
<b>Performance profile</b>	micro tube/HE9 680	micro tube/HE9 680	HE7/HE7 680	HE7/HE7 680	HE7/HE7 680
Maximum Power Output (dB SPL)					
2 cc coupler	123/125	124/125	128/125	135/130	136/131
Ear Simulator	128/131	129/131	135/131	140/138	141/139
Maximum gain (dB)					
2 cc coupler	46/52	53/55	58/53	68/63	75/70
Ear Simulator	58/60	61/63	61/58	71/68	80/75
Frequency range – 2 cc coupler (Hz)	<100 – 7600	<100 – 7100	<100 – 6800	<100 – 6100	<100 – 6000
Battery size	10/312	13	13	13	13
Working current (mA)	1.0	1.3	1.1	1.2	1.2


## Fitting range




Mild to moderately-severe hearing loss, all audiometric configurations




Mild to moderately-severe hearing loss, all audiometric configurations



Mild to moderately-severe hearing loss, all audiometric configurations




Moderate to severe hearing loss, all audiometric configurations



Moderately-severe to profound hearing loss, all audiometric configurations

## Technical description

- 1 Transparent hook (or SlimTube)
- 2 Sound outlet
- 3 Miniature microphone with protector
- 4 Program switch
- 5 Access to programming socket
- 6 Digital volume control
- 7 Serial number
- 8 Battery compartment with ON/OFF switch
- 9 Dome or SlimTip (small image)
- 10 Retention aid
- 11 micro tube (or transparent hook)



## Custom product options

	Certéna CIC/MC Petite	Certéna CIC/MC P Petite	Certéna ITC/HS Petite	Certéna ITC/HS P Petite	Certéna ITC/HS	Certéna ITC/HS P	Certéna FS	Certéna FS P
Volume control	0	0	0	0	0	0	0	0
Program switch TacTronic	0	0	0	0	0	0	0	0
EasyPhone			0	0	0	0	0	0
Telecoil			0	0	0	0	0	0
SmartGuard	0	0	S	S	S	S	S	S

S = Standard O = Optional Blank = not available

## BTE options

	Certéna micro Petite	Certéna micro	Certéna M	Certéna P	Certéna SP
Hook	HE7 680 HE7		S	S	S
Hook HE7/HE7 680 mini hook HE7 680 mini hook HE7 680			0	0	0
HE9 680 mini hook HE9 680 HE9 HE 9	S 0	S 0			
Volume control				S 0	S 0
Version 1 Version 2 Cover				S 0	S 0
Program switch cover		0			
EasyPhone	S	S	S	S	S
Telecoil	S	S	S	S	S
SlimTube (open/close) with SlimTip/Dome		0			
micro tube (open/close) with SlimTip/Dome	0	0			
Tamper proof battery compartment		0	0	0	0

S = Standard O = Optional Blank = not available

## Color options

<b>Custom products</b>	<b>Faceplate:</b> Cocoa (22), Pink (26), Brown (28), Tan (14) <b>Shell:</b> Cocoa (22), Pink (26), Brown (28), Tan (14), Red transparent (36), Blue transparent (37), White (19), Transparent (21)
<b>micro Petite/micro</b>	Beige (01), Light gray (04), Pure transparent (13)*, Gray (23), Taupe (35), Palladium (40), Pearl (41), Sahara (64), Cobalt (65), Graphite (66), Cinnamon (67), Light pink (33), Light blue (34)
<b>BTE</b>	Beige (01), Brown (03), Taupe (35), Light gray (04), Black (06), High-Tech gray (12), Gray (23), Blonde (49), Chestnut brown (50), Pure transparent (13), Red transparent (36), Blue transparent (37), Purple transparent (38), Light pink (33), Light blue (34)

\* not available for Certéna micro Petite  
Note: Some options are country dependent. Please refer to your local order form.

## Certéna features and accessories overview

Features Sound Recover	Custom Products								BTE	micro Petite	Micro	M	P	SP	
	CIC/MC Petite	CIC/MC P Petite	ITC/HS Petite	ITC/HS P Petite	ITC/HS	ITC/HS P	FS	FS P							
VoiceZoom															
SoundFlow Standard															
WhistleBlock Technology															
NoiseBlock Processing															
*QuickSync															
Acoustically Optimized Vent (AOV)															
<b>Accessories</b>															
myPilot															
iCom with EasyAudio and EasyBluetooth															
iView															
KeyPilot2															
WatchPilot2															
Spectacles adaptor															
Click'nTalk															
<b>FM receivers</b>															
MyLink															
MicroMLxS or MLxi															
ML9S or ML9i															
ML12i															
Audio shoe AS9-MLxS															
CROSLink															
<b>FM transmitters</b>															
inspiro, SmartLink+, ZoomLink+ or EasyLink+															

\* requires binaural fitting  
\*\* requires remote control option  
\*\*\* requires Telecoil option  
\*\*\*\* requires iCom