

Your Roger On™ welcome guide.

Congratulations on starting your Roger On trial!

Roger™ is your third ear that makes sure you don't miss life's beautiful moments. It allows you to follow conversations, stay close to friends and family, and enjoy an active social life.

Ready to hear what you've been missing?

This guide will walk you through everything you need to know to get the most from your Roger On microphone.

Getting started with your Roger On

1 Set up RogerDirect™

This is a one-time step to connect Roger with your hearing aids. Your audiologist may have already done this during your last clinic visit.

If not, refer to the step-by-step instructions in the RogerDirect Setup Guide included in the box or watch the tutorial video on the Phonak YouTube channel.



Scan here to watch the tutorial video.

2 Get familiar with Roger On

As you unbox your new microphone, take a moment to go through the Quick setup guide provided in the box. It covers what's inside, explains the icons that appear on the screen, and gives you examples of everyday usage.

3 Charge before use

It takes 2 hours to fully charge. A full battery can provide up to 8 hours of continuous use.

4 Download the myRogerMic app

The app lets you adjust Roger's settings easily. It's optional but a great way to enhance your experience.



Google Play



App Store

And that's it - you're all set and ready to bring Roger along wherever you go.



Put Roger On to the test

With Roger on your side, you're equipped with a powerful and versatile hearing solution. Try it in different settings and note down whether it helped (Yes/ No). When your trial period ends, we recommend sharing your findings with your audiologist to get the most from your experience.

At social gatherings	Recommended use	Yes	No
At the coffee shop/ a noisy restaurant, dining with friends and family	Place Roger On in the center of the table. You can select the microphone direction	<input type="checkbox"/>	<input type="checkbox"/>
At a social event, with loud music and people talking	Hold it in your hand and point it towards the person speaking	<input type="checkbox"/>	<input type="checkbox"/>
At the shopping mall/ grocery store	Hold it in your hand and point it towards the person speaking	<input type="checkbox"/>	<input type="checkbox"/>
Sports and other activities			
At the gym or training outside	Give Roger On to your coach/ teammate and ask them to clip it on their shirt. Alternatively, they can wear it around their neck by using the lanyard	<input type="checkbox"/>	<input type="checkbox"/>
Hiking or biking	Give it to a fellow hiker/ biker to clip on their shirt while you're on the go	<input type="checkbox"/>	<input type="checkbox"/>
At home			
Conversation with your partner or a friend when you're not face-to-face (e.g. in a different room)	Give Roger On to your partner/ friend to clip it on their shirt	<input type="checkbox"/>	<input type="checkbox"/>
Watching TV	Use the docking station to connect Roger On to your TV	<input type="checkbox"/>	<input type="checkbox"/>
When travelling			
Have a conversation with your travel companion, while in the car or in public transportation	Give Roger On to your travel companion to clip it on their shirt	<input type="checkbox"/>	<input type="checkbox"/>
Business events			
Lectures and presentations	Give Roger On to the presenter and sit comfortably among the audience	<input type="checkbox"/>	<input type="checkbox"/>
Panel discussion	Place it on the stage, ideally on a table, where the panel discussion takes place. Open the mobile app and choose the Wide pointing mode. You can now sit comfortably among the audience and listen to the discussion	<input type="checkbox"/>	<input type="checkbox"/>
At work			
Team workshops	Place Roger On in the middle of the table when having a team discussion. Give it to the instructor to clip it on their shirt	<input type="checkbox"/>	<input type="checkbox"/>
Standing meetings and hallway conversations	Hold it in your hand and point it towards the person speaking	<input type="checkbox"/>	<input type="checkbox"/>
Small meetings	Place it in the middle of the meeting room table	<input type="checkbox"/>	<input type="checkbox"/>
Online calls	You can also use Roger On as a headset. First, use the cable to connect it to your laptop. Then clip it on your shirt. It will pick up your voice and send it to the computer	<input type="checkbox"/>	<input type="checkbox"/>