# Acoustic Leakage Tester ALT

## User Guide

Easy Line™  
Accessories

# Acoustic Leakage Tester ALT

---

English	2
Deutsch	52
Spanish	102
Français	152
Dutch	202

---

# Contents

1	Welcome	p. 3
2	Introduction	p. 4
3	Language	p. 5
4	What's included	p. 6
5	ANSI table	p. 8
6	First use	p. 10
7	Settings	p. 12
8	Calibrate	p. 14
9	Measure with client report	p. 18
10	Batch report	p. 30
11	Quick measure	p. 37
12	General remarks	p. 44
13	Problem solving	p. 47
14	Safety and precautions	p. 48
15.	Information and explanation of symbols	p. 50

# 1. Welcome

Thank you for buying the Acoustic Leakage Tester (ALT).

ALT is a portable leakage tester for use with both custom and universal earpieces.

ALT delivers single-button measurement and an automated calibration routine, ensuring accurate results every time. Comfort and correct fitting is paramount for the client.

## Features:

- Compatible with Apple mobile devices (iPhone, iPad, iPod)
- Acoustic Sensor Unit replaces filter during testing
- Manufactured for use with Sonova hearing protection
- Leakage testing of both custom and universal earpieces
- Acoustic measurement
- One-button calibration
- One-button test
- Batch report

This User Guide will help you to connect and use the ALT correctly.

## 2. Introduction

The Acoustic Leakage Tester (ALT) has been developed to support supervisors and professionals to ensure that personnel using custom hearing protection are appropriately protected.

It can also be used with universal hearing protection to confirm whether the user is using the correct size. ALT is used to measure the quality of the acoustic seal of the earpiece with the ear canal. The seal is very important for the hearing protection to function correctly.

The earpiece should be able to maintain a good seal by exerting some pressure on the ear canal wall. In order to achieve a seal, the outer surface of the earpiece should match the shape of the ear canal as closely as possible.

ALT checks the quality of the acoustic seal by generating a low level continuous sound in the ear and measuring the response.

### 3. Language

The App is available in five languages: EN, DE, FR, NL, ES. If one of these languages is the default language of the used device, the App will automatically run in that language.

If the device uses another language, the App will automatically run in English.

You can change the App's language to any of the five supported languages at any time. To do this, open the device's Settings section, then scroll down to the Easy Line ALT App. Tap on the App, and then scroll to the Language section – here you can choose from any of the five available languages.

## 4. What's included

The ALT kit consists of a probe that can be connected to both custom and universal earpieces, a calibrator, 3 soft adapters, a quick set-up guide and a carrying case.

To make ALT kit complete, you should download the Easy Line ALT App on the device you will use with the kit. You can find the App by searching in the App store for:


The ALT kit contains the following items:

extraction tool for  
silicone shells

Calibrator

Probe (Acoustic  
Sensor Unit)

Soft adaptor for mounting  
the probe on the calibrator  
and acrylic shells

Not included

USB-C jack adaptor  
(dongle)

Lightning to 3.5mm jack  
adaptor (dongle)

Easy Line ALT

The Apple Store logo is a trademark of Apple Inc.

## 5. ANSI table

### ALT 2 performance declaration

#### Specification


A	Manufacturer
B	FAES brand name
C	System requirements
D	HPDs supported
E	FAES test method
F	FAES type
G	Recommended regular physical calibration
H	Recommended periodic verification procedure
I	Maximum permissible ambient noise
J	Maximum and minimum measurable attenuation values
K	Fit-test output
L	FAES measurement uncertainty
M	Listener requirements

Illustrative values and/or possible descriptors	References to clauses of ANSI S12.71-2018
Sonova Communications AG	
Easyline	
Apple device running iOS 13 or higher.	
The device requires a jack port or an official Apple lightning to jack adapter or USB-C to jack adapter	
Include list of earplug products. (Sonova, Audionova, Phonak, canal type, filter type)	
Attenuation derived from an acoustic leakage measurement	
Physical, using standard HPD	Clause 4
Every two years	Clause 5.1
Daily on start up	Clause 5.1
60 dBA and 65 dBC	Clause 5.2
Able to measure PARs ranging from	Clause 5.5
0dB to 42dB	
Binaural PARx value on report or Fail.	Clause 8.1
3,77 dB for all earmoulds	Clause 8.2.1
Not applicable	Clause 11

## 6. First use

Prior to the first use of the ALT system, please take the following steps:


- Switch on the iPhone®, iPad® or iPod®.
- Create an Apple ID if you do not have one yet.
- Open the App store, search for the Easy Line ALT App, install and start the application.


Easy Line ALT

Incoming calls and messages may disturb the measurement.


- Be sure to put your device on “**do not disturb**”.

## 7. Settings

Press '**Settings**' to enter the Operator menu.


Fill out the ALT User (Operator) information.

Define the PAR (Personal Attenuation Rating) calculation setting:

- ANSI (US Standard)
- CE (European Standard)

For more information click on the 'i' information button.

Portfolio choice:  
Choose your hearing protection brand.

Press '**Save**'.

## 8. Calibrate

In order to make accurate measurements, the system must always be calibrated first.



Press '**Calibrate**' in the Welcome screen to start the calibration process, then follow the on-screen instructions:

### 1. Connect the probe

- Plug the 3.5 mm jack connector of the probe into the 3.5 mm headphone socket of the device.
- If you use a device without a 3.5 mm headphone socket, you should use the included dongle (Lightning to headphone jack adapter).

## 2. Connect the calibrator


- Gently assemble the probe, the soft adapter and the calibrator together.
- Push and rotate the calibrator until the soft adapter and probe are sitting flush within the calibrator, with no gap between the calibrator and soft adapter, or the soft adapter and the probe.


### 3. Calibrate

- Depending on your device, you will be asked to select the audio hardware input source: Jack or Dongle.  
If you use the Dongle you should give the device some time to recognize that the probe is connected (about 10 seconds).
- The software will now calibrate the system making it ready for use. This will take about 10 seconds.
- Press '**Calibrate**' on the bottom of the screen.


Always allow access to microphone

When calibration is complete, one of the following messages will be displayed:

You can continue


Try again: check the connection of the jack connector in the device and press '**Calibrate**' again. See the '**Problem solving**' section below.


## 9. Measure with Client report

(for custom mold only)

Press '**Measure with client report**' to perform a measurement and generate a report.

- Before the first measurement you must perform a calibration to ensure that the system works properly (Please refer to the '**Calibrate**' section above).


Enter the requested information:

- Company name
- Company E-Mail address (the measurement report PDF will be sent to this address)
- Report language

Press '**Next**'.

For other product types:

- Select the '**Mold type**' and the '**Filter type**' of the hearing protection to be tested.


Name Surname	
S/N	2222C007
S/N	2222C008
Filter Phonak KI20	

For Phonak Communications products:

- Select the Mold type.
- Select '**Scan product type**' and take a picture of the sticker with end-customer information (originally on the packaging sleeve).

The information will be automatically added to the client information.


Press '**Next**'.

Enter the requested information:

- Client name
- Right and left SN (serial number of the earpiece)
- Comment (any helpful or additional remark)


Press '**Next**'.

Insert the probe into the right earpiece by removing the filter as follow:

For acrylic earpieces:

- Remove the filter by pulling up the silicon ring support.


For silicon and universal earpieces:

- Remove the filter by pressing the earpiece under the filter or gently insert the extraction tool into the filter slot.


Place the earpiece with the probe attached into the right ear.

Press '**Next**'.

At this point you can:


- Start sound (to start the measurement process).
- Calibrate (to redo the calibration please refer to the '**Calibrate**' section above).

## Press '**Start sound**'

The system will generate a continuous sound at a very low level (which is not dangerous to the user), during this time you can evaluate the acoustic seal of the earpiece.

- When the needle is in the green area it means that the acoustic seal of the tested earpiece is sufficient.
- When the needle is in the red area, it means the acoustic seal of the earpiece is insufficient and too poor to be used as hearing protection.


The needle can also be in the red area when:


- the probe is not correctly placed in the earpiece.
- the earpiece is not correctly placed in the ear canal.

When the needle drops below the red area you get the message '**Signal Too Weak!**'.


Press '**Ok**', then check if everything is correctly placed and press '**Start sound**' to measure again.


- Press '**Stop sound**' if you want to check everything is correctly placed.
- Press '**Measure result**' to start the measurement: the system will make a five second measurement and give the result.


Press '**Ok**'.

- If the measurement result is insufficient you can remeasure the right earpiece by pressing '**Remeasure**'.


- When you finished the measurement of the right earpiece, press '**Measure other ear**'.


- Disconnect the probe from the right earpiece and place it into the left earpiece.
- Place the left earpiece with the probe attached into the left ear.


Press '**Next**', then follow the same steps as for the right ear.

Press '**Next**' to access the Report screen.

- On the Report screen you can double-check the E-Mail recipients or add an additional E-Mail recipient.
- Press '**View**' to see the detailed report.
- Press '**Send Report**'


Press the blue harrow to send the report to the chosen E-Mail addresses.

## 10. Batch Report

You can use the batch report functionality if you want to send several reports at the same time.


On the '**Home**' menu, press '**Batch Report**'

Enter the requested information:

- Company name
- Company E-Mail address  
(the measurement report PDF will be sent to this address)
- Report language

For other product types:

- Select the '**Mold type**' and the '**Filter type**' of the measured hearing protection.

Press '**Next**'.

Name Surname  
S/N 2222C007  
S/N 2222C008  
Filter Phonak KI20

For Phonak Communications products:

- Select the Mold type.
- Select '**Scan product type**' and take a picture of the sticker with end-customer information (originally on the packaging sleeve).

The information will be automatically added to the client information.


Press '**Next**'.

Enter the requested information:

- Client name
- Right and left SN (serial number of the earpiece)
- Comment (any helpful or additional remark)


Press '**Next**'.

Insert the probe into the right earpiece by removing the filter as follow:

For acrylic earpieces:

- Remove the filter by pulling up the silicon ring support.

For silicon and universal earpieces:

- Remove the filter by pressing the earpiece under the filter or gently insert the extraction tool into the filter slot.


Place the earpiece with the probe attached into the right ear.

Press '**Next**'.


At this point you can:

- Start sound (to start the measurement Process).
- Redo the same procedure as done for '**Measure Client**' – see page ....
- Calibrate (to redo the calibration – please refer to the '**Calibrate**' section above).

When you have finished the measurement, press '**Next**'.


- You will see the value of the Personal Attenuation Rating (PAR).
- Press '**Measure Next**' to measure the next person.
- Press '**Finish**' to access the report.

On the Report screen you can:

- View the report
- Double-check the E-Mail recipients or add an additional E-Mail recipient.

Press '**Send Report**'

Press the blue harrow to send the report to the chosen E-Mail addresses.


# 11. Quick Measure

(for custom and universal mold)

Quick Measure allows you to take a measurement without generating a report.


- Press '**Quick Measure**' to perform a measurement without generating a report.
- Before the first measurement you have to perform a calibration to ensure that the system works properly (Please refer to the '**Calibrate**' section above).

Insert the probe into the right earpiece by removing the filter as follow:

For acrylic earpieces:

- Remove the filter by pulling up the silicon ring support.

For silicon and universal earpieces:

- Remove the filter by pressing the earpiece under the filter or gently insert the extraction tool into the filter slot.


Place the earpiece with the probe attached into the right ear.

Press '**Next**'.

At this point you can:

- Start sound (to start the measurement Process).
- Calibrate (to redo the calibration – please refer to the '**Calibrate**' section above).


### Press '**Start sound**'


- The system will generate a continuous sound at a very low level (which is not dangerous to the user), during this time you can evaluate the acoustic seal of the earpiece.
- When the needle is in the green area it means that the acoustic seal of the tested earpiece is sufficient.
- When the needle is in the red area, it means the acoustic seal of the earpiece is insufficient and too poor to be used as hearing protection.

The needle can also be in the red area when:


- the probe is not correctly placed in the earpiece.
- the earpiece is not correctly placed in the ear canal.


- When the needle drops below the red area you get the message '**Signal Too Weak!**'.


- Press '**Ok**', then check if everything is correctly placed and press 'Start sound' to measure again.
- Press '**Stop sound**' if you want to check everything is correctly placed.
- Press '**Measure result**' to start the measurement: the system will make a five second measurement and give the result.

Press 'Ok'.

If the measurement result is insufficient you can remeasure the earpiece by pressing 'Remeasure'.

- When you have finished the measurement of the earpiece, press '**Finish**' to go back to the '**Home**' screen.


Make sure to re-insert the filters into the earpieces after the measurement.

## 12. General remarks

For best results, please follow the recommendations below:

- The person under test should keep their head as stable as possible. Jaw movements, swallowing, speaking and taking a breath may influence the measurement.
- The cable of the probe should be kept still during the measurement process, as microphonics produced by the cable rubbing on clothing etc. can, in some cases, influence results.
- It is recommended to carry out testing in a quiet environment, as loud noises (> 70 dB(A)) may influence the measurement.
- Cerumen build up in the ear canal can affect results due to poor sealing. Please ensure the canal is clean prior to testing.
- Incorrect positioning of the earpiece in the ear canal may show an insufficient acoustic sealing and result in a failed test. You can try to reposition the earpiece and measure again.

- Changing the volume of the device has no effect on the sound levels of the leakage tester tones, as the volume is fixed by the ALT App.
- Please ensure other apps are not running in the background when using the ALT App, as we are unable to guarantee the compatibility of apps from third parties.
- Switching off the Volume limit in Music settings might be a solution when the ALT App doesn't pass the Calibration phase. Go to the settings App, select 'Music', select 'Volume Limit' and switch off 'EU-volume limit'.
- In general, you should calibrate the system before every use. If more than 15 minutes have passed since the last calibration, the following message will be displayed:

- Press 'Ok' and redo the calibration (please refer to the 'Calibrate' section on page x).
- With some people the noise produced by blood flow in the ear can make the measurement less stable. It could cause the measurement result to fluctuate around a certain average based on the frequency of the heartbeat.
- If any problem occurs that cannot be solved using the recommendations above, please visit the following link:

[www.sonova-communications.com/alt-2-0-software/](http://www.sonova-communications.com/alt-2-0-software/)

## 13. Problem solving

If measurements made through the ALT App are unstable, please follow the recommendations in 'General remarks' above.

If measurement values continue to be unstable please carry out the following steps in the order shown below:

- shut down other apps that use audio input or output
- restart the ALT App
- restart the device.

If the measurement value does not stabilize after carrying out these recommendations, or another problem occurs that cannot be solved, please visit the following link:

[www.sonova-communications.com/alt-2-0-software/](http://www.sonova-communications.com/alt-2-0-software/)

## 14. Safety and Warnings

- Please treat all parts with appropriate care. Protect the device from excessive shock and vibration.  
Do not use excessive force when connecting your device to its different cables.
- When inserting or removing the probe from a soft silicone custom earpiece, care should be taken as vigorous twisting can cause damage to the earpiece material.
- Please keep the probe and calibrator clean. Prevent contact with moisture, dust and grease. Clean using a damp cloth. Never use household cleaning products (washing powder, soap, etc) or alcohol to clean the product.
- Do not immerse any part of the ALT in water.
- Keep this product and its accessories out of reach of children below 36 months. Contains small parts and cords that may pose risk of choking or strangulation.

- When you are not using your product, store it in a safe, dry and clean place.
- Do not use your device in areas where electronic equipment is prohibited.

iPhone, iPad and iPod are trademarks of Apple Inc. registered in the U.S. and other countries.

## 15. Information and explanation of symbols

---

Indicates the device manufacturer.

The CE symbol is confirmation by Sonova Communications AG that this product is in compliance with the essential requirements and other relevant provisions of the Electromagnetic Compatibility (EMC) Directive 2014/30/EU

---

The UKCA symbol is confirmation by Sonova Communications AG that this product meets the applicable regulatory requirements for Great Britain.

---

This symbol confers that the electromagnetic interference from the device is under limits approved by the US Federal Communications Commission.

---

Australian EMC and Radiocommunications compliance label."

---

The symbol with the crossed-out garbage bin is to make you aware that this device may not be thrown away as normal household waste.

Please dispose of old or unused device, at waste disposal sites intended for electronic waste, or give your device to your hearing care professional for disposal. Proper disposal protects the environment and health.

---

This symbol indicates that it is important for the user to read and take into account the relevant information in this user guide.

---

Temperature transport and storage: -20° to 60° Celsius  
(-4° to +140° Fahrenheit) Temperature  
operation: 0° to +45° Celsius (+32° to +113°  
Fahrenheit).

---

Keep dry

---

Humidity transportation: up to 90%  
(noncondensing). Humidity storage: 0% to  
70%, if not in use. Humidity operation: <95%  
(noncondensing).

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iPhone, iPad en iPod zijn handelsmerken van Apple Inc. geregistreerd in de VS en andere landen.
Lees de informatie op de volgende pagina's voordat u dit product van Sonova Communications AG gebruikt.