

Phonak Target 11

Phonak Target/ALPS User Guide



This guide contains details on using Phonak Target/ALPS available in Phonak Target fitting software. Phonak Target/ALPS allows you to enter data via Phonak Target directly into ALPS. This one step process avoids duplicate entry and allows you to manage all relevant Lyric fitting and subscription information in an easy way.

To use Phonak Target/ALPS, the computer you are using for fitting must have an internet connection.

For every Lyric fitting, Phonak Target will automatically connect you to ALPS. A wizard will guide you through all the necessary steps to complete a successful Lyric fitting. Please note that the Lyric device activation code must be entered the first time Phonak Target is launched.

More information is also be available under Web Help.


Content

Create a new client in Phonak Target	2
Access an existing client in your ALPS account	6
Device exchange (including ear canal observation)	7
Tourist case: Device exchange for a tourist client	10
Tourist case: Your client has been refit by a different Authorized Lyric Partner	11
Renew a subscription	12
Cancel of an active subscription	13
Offline fitting and synchronization	14
ALPS home and reports	16
Enhanced log files for support	20
Information and description of symbols and System Requirements	20

Create a new client in Phonak Target

To create a new client in Phonak Target:

1. Click on [Instruments] and select Lyric devices in [Hearing aids]
2. Click on [Lyric configuration] and click on [Select device] to open the Phonak Target/ALPS wizard.




The ALPS button indicates your current connection status:

- Green - connected
- Black – disconnected




To connect into ALPS, enter your user name and password to log into your Lyric account.


In case you have forgotten your password, click on the link [forgot password] and change it in www.phonakpro.com directly.



You will be asked to [Add new client to your ALPS account]. Your account information is shown on this screen, that is very helpful for providers working with multiple logins.



In case you don't intend to login with the last login shown or you want to remove the automatic login checkbox, choose [Switch login] under the ALPS button and login with another username.



You can now see the client’s first and last name, date of birth and gender as entered in Phonak Target.

The additional information (date of birth and gender) is optional. You can uncheck [Submit additional information] if you do not wish to transmit this data.

For new clients, you are then asked to confirm that the client understands and gave consent that his/her data will be transferred to Sonova for processing. If the box is not ticked, the fitting cannot continue.

The screenshot shows the 'Enter client details' form in the Lyric software. It includes input fields for 'First name' and 'Last name', both with placeholder text 'Please enter...'. The 'Date of birth' is set to '01.01.1948' and 'Gender' is set to '-'. There is a checked checkbox for 'Submit additional information'. At the bottom, there is a confirmation statement: 'I hereby confirm that I fully informed my client that his personal data will be transferred to the hearing aid manufacturer for processing purpose. My client consented to such transfer.' Below this statement are two buttons: 'Create Lyric ID' and 'Cancel'.

To select the appropriate Lyric device, choose the measured seal size and then match the serial number of the device you intend to use from your stock with the matching serial number in the system. Devices about to expire within 30 days are marked red.

Please note that expired devices can no longer be fit. To return an expired device, use the [Device Return process] in ALPS, select the serial number and the reason for device return as “expired” and send it back to Phonak.

The screenshot shows the 'Select Lyric device from consignment stock' screen. It displays a table of devices with columns for 'Serial#', 'Model', 'Seal size', 'XXL', and 'Insert by'. The 'Insert by' column shows dates, some of which are in red, indicating expiration. A search bar at the top left is labeled 'Serial...'. The 'Former device seal size: XXL' is noted at the top. The 'Insertion date' is '14.03.2023'. At the bottom right, there are 'Select' and 'Cancel' buttons.

Serial#	Model	Seal size	XXL	Insert by
2307NYSEH	Lyric4	XXL		30.03.2023
2307NYSEK	Lyric4	XXL		05.04.2023
2307NYSEL	Lyric4	XXL		05.04.2023
2307NYSEJ	Lyric4	XXL		09.04.2023
2307NY5F2	Lyric4	XXL		13.04.2023
2307NY5EP	Lyric4	XXL		17.02.2024
2307NY5EW	Lyric4	XXL		17.02.2024
2307NY5EV	Lyric4	XXL		17.02.2024
2307NY5EU	Lyric4	XXL		17.02.2024


Select the measured insertion depth.

The screenshot shows the 'Select measured insertion depth' screen. It features a dropdown menu labeled 'Initial device insertion for this ear.' with a list of insertion depths from 3 mm to 15 mm. The '7 mm' option is currently selected and highlighted. At the bottom right, there are 'Continue' and 'Cancel' buttons.

Choose the subscription duration you want to create for your client. Make sure to click the check box, in case the customer is using Lyric for the treatment of tinnitus.


The screenshot shows the 'Select a subscription' screen. It includes a 'Subscription type' dropdown set to '1-year subscription', a 'Payment type' dropdown set to 'Yearly upfront', and a 'Start date' field set to '11.11.2019'. There is a checkbox labeled 'Lyric is used for tinnitus treatment' which is currently unchecked. At the bottom right, there are 'Continue' and 'Cancel' buttons.

Phonak Target confirms the device selection. You can get [More information] about the subscription in this window or choose [Close]. For binaural fittings, you are asked to [Proceed with other side] which leads you directly to selecting another device from your consignment stock.




Below the start and end date of the subscription and actual days remaining to the trial end date are shown.

It is now possible to extend a trial once for an amount of 10 days. The new end date will then be shown in Target/ALPS.



In case a wrong serial number has been chosen, it can be returned to the stock by clicking “return device to stock”.




You can now move to programming the Lyric devices on the fitting tab.


You have the possibility to choose between Gain (Real ear and Coupler 0.4cc) and Output display (SPL real ear and coupler 0.4cc) when fitting Lyric. Click on the [Output] menu to select your desired output curve display. When selecting the output curve display, you may also choose to show or hide the hearing level and UCL.




In the event the achieved insertion depth varies from the measured one, both values can be entered in [Fitting]. A warning will be shown if the achieved insertion depth is higher than the measured one.




By clicking [Save to database] at the end of the session, your fitting data will be synchronized with ALPS.




In case you need to have more notes stored about the patient, fill them in on the new notes section.



A new note can be added and Target allows a “remind me again in next session” function.




All note are then visible under [Client, Notes].



The following reports can be printed from the Print menu:


For Professional:

- Fitting Report
- Lyric Datasheet



For Client:

- Hearing Aid Instructions (SoundLync)
- Communication tips
- Diary of observation
- Understanding hearing loss



The curves reports can be printed as follows:

Choose [Setup], [Reports], [Curve type] and choose between Gain or Output. The selected curve type will then be visible on the curve report.

If the client exists in ALPS but not in Target (or a first access within a satellite office), you can [Access existing client in your ALPS account].

Enter the correct Lyric ID and/or the first and last names as stored in ALPS.

The Lyric ID has to be entered if you are connecting an existing client for the first time. When opening the session again, you will automatically be connected to the ALPS data.

By selecting the search button the P# will populate, if the P# does not populate, this means the patient is not an existing patient or the name does not match the ALPS information.

If the data of the client was found, the client subscription information screen will show all the relevant information.

Lyric

Client has no Lyric ID assigned in Phonak Target

Account 0020035121 LYRIC PORTAL - TARGET TEAM

☐ Add new client to your ALPS account

☒ Access existing client in your ALPS account

☐ Access existing client in a different ALPS account (Tourist case)

Continue

Cancel

Lyric

Lyric client

First name

Last name

Lyric ID

Date of birth

-

Gender

-

Search




Lyric

Client subscription information

First name John

Last name Doe

Lyric ID P000087661

Information	 Right	 Left
Subscription	Trial  (20 days left)	No subscription
Start date	06.06.2014	
End date	05.06.2015	
Remaining devices	11	
Device	Lyric3	
Serial number	1406NY0JT	

Continue

Cancel

Device exchange (including ear canal observation)

To exchange a Lyric device, please log into Phonak Target/ALPS, go to [Instruments], [Lyric configuration] and click [Exchange/remove device].

Subscription

Status: Trial (27 days left)

Start date: 07.11.2017

End date: 07.12.2017

Trial end date: 07.12.2017

Billing date: 07.12.2017

Remaining devices: 0

Lyric3

Serial number: 1733NY3HT (Code: E325)

Seal size: S

Insertion date: 09.11.2017

Insertion depth: Measured 10 mm

Exchange/remove device

Select if the provider or the patient has removed the device and specify the removal date. Choose the corresponding removal reason. The removal reasons are separated into four categories (Discomfort, Device related, Fitting related, Proactive removal and Cancelled Subscription).

Remove Lyric device

Serial number: 1849NY3G3

Removed by: ☒ Provider ☐ Patient self-removed

Removal date: ☒ Today ☐ Previously

Removal reason: --select--

- Discomfort
- Device related
- Fitting related
- Proactive removal
- Cancelled subscription

Remove Cancel

Hovering over individual removal reasons activates the “tooltips”, which gives more information to help determine which reason to select.

Remove Lyric device

Serial number: 1849NY3G3

Removed by: ☒ Provider ☐ Patient self-removed

Removal date: ☒ Today ☐ Previously

Removal reason: --select--

- Discomfort
- Device related
- Fitting related
 - Occlusion
 - Feedback
 - Insufficient gain
 - Sound quality
 - Migrating
- Proactive removal
- Cancelled subscription

Device is not loud enough for patient hearing loss

Remove Cancel

Phonak Target will ask you for the status of the ear canal. You can choose between: [OK for immediate refit], [OK for refit after rest] and [Requires medical referral]. When choosing [OK for immediate refit], the device exchange can happen immediately. If the ear canal was irritated and needed rest or a medical referral, the system will inform you accordingly.*

The comments for removal reason and ear canal observation are both optional.

When [Requires medical referral] is chosen, additional details about the state of the ear canal must be chosen.

Additionally, if “other reason” is chosen, a comment to give more explanation is required.

To remove a second device, click on [Proceed with other side]. There are some reasons (scheduled removal, near end of device life and clinical proactive removal) that will automatically populate the second side as a default but it can be changed if the removal reason is different.

*If you have chosen [OK for refit after rest] or [Requires medical referral], Phonak Target will inform you on the main screen.

Lyric

Remove Lyric device

Serial number

1633H00KR

Removed by

☐ Provider

☒ Patient self-removed

Removal date

☒ Today

☐ Previously

Removal reason

Discomfort

Additional details (optional)...

Ear canal observation

OK for immediate refit

Additional ear canal observation (optional)...

Remove

Cancel

Lyric

Remove Lyric device

Serial number

1849NY3G3

Removed by

☒ Provider

☐ Patient self-removed

Removal date

☒ Today

☐ Previously

Removal reason

Discomfort

Additional details (optional)...

Ear canal observation

Requires medical referral

Additional details

--select--

Redness of tissue

Excess fluid collection on tissue

Swelling of clotted blood below tissue (bruise)

Sore or ulceration of tissue

Medial bulge/growth

Blood/Bleeding

Cerumen management

Other reason

Remove

Cancel

Lyric

Successfully removed

Device successfully removed.

Print return form

Select new device

Close

Select device

OK for refit after rest (13.11.2017)

Redness of tissue - Mild

Before you can exchange a device, it is mandatory to confirm that the [Ear is ok for immediate fit (today)] by checking the box at the bottom of the screen. Adding the comments about the patient's ear canal is optional.

Lyric

Ⓡ Remove Lyric device

Serial number: 1736NVHCG

Removed by: ☒ Provider ☐ Patient self-removed

Removal date: ☒ Today ☐ Previously

Removal reason: Discomfort ▼

Additional details (optional)...

Ear canal observation: OK for immediate refit ▼

Additional ear canal observation (optional)...

Remove Cancel

After confirming the ear is OK for immediate fit, choose another Lyric device to exchange.

Lyric

Ⓡ Select Lyric device from consignment stock

Former device seal size: M

Serial... Insertion date: 16.06.2014

Serial#	Model	Seal size	M ▼	Insert by
1406NY0JV	Lyric3	M		08.05.2015
1406NY0JN	Lyric3	M		08.05.2015
1406NY0K0	Lyric3	M		08.05.2015
1406NY0IX	Lyric3	M		08.05.2015
1406NY0JJ	Lyric3	M		08.05.2015
1406NY0JH	Lyric3	M		08.05.2015
1406NY0JM	Lyric3	M		08.05.2015
1406NY0JK	Lyric3	M		08.05.2015

Select Cancel

The comments you entered are saved under [Client, Notes].

Fitting

Notes

Enter new session note...

16.06.2014 13:57: ECO Right: OK for refit after rest.
ECO Right: OK for refit after rest.
Ear canal observation after device removal.
Ear canal is a bit irritated.

16.06.2014 13:56: ECO Left: OK for immediate refit.

16.06.2014 13:55: ECO Right: OK for immediate refit.

Please note that if the device has been worn longer than 30 days, the [Easy Refit] is enabled and will show default values:

- Removed by provider
- Removed today
- Removal reason: scheduled removal
- Ear canal observation: "ok for immediate refit"

Click [Remove] and then choose a new device accordingly.
When the criteria doesn't match the default value, please change it and continue.

Lyric

Ⓛ Remove Lyric device

Serial number: 1836NV212

Removed by: ☒ Provider ☐ Patient self-removed

Removal date: ☒ Today ☐ Previously

Removal reason: Scheduled removal ▼

Additional details (optional)...

Ear canal observation: OK for immediate refit ▼

Additional ear canal observation (optional)...


Remove Cancel

Tourist case: Device exchange for a tourist client


If you need to exchange a Lyric device for a client who was originally fit by another Authorized Lyric Partner, you need to obtain access to the client's ALPS fitting data. To do so you need the client's individual Lyric ID.

To start the tourist case, proceed as you would with a new client. Log into your ALPS account and select [Access existing client in a different ALPS account] (Tourist case)].

Enter the client's Lyric ID, and check if the name is displayed correctly (as printed on the Lyric ID card). Proceed with the exchange of the device, and adjust the fitting parameters if necessary. If the patient does not know their Lyric ID, it can be obtained by calling Lyric customer service.



The screenshot shows the Lyric ALPS Portal interface. At the top, it says "Lyric". Below that, a message states: "Client has no Lyric ID assigned in Phonak Target". The account information is "Account 0020035121 LYRIC PORTAL - TARGET TEAM". There are three radio button options: "Add new client to your ALPS account", "Access existing client in your ALPS account", and "Access existing client in a different ALPS account (Tourist case)". The third option is selected. A blue information icon with a text box below it says: "Session data in Phonak Target will be updated from ALPS." At the bottom right, there are "Continue" and "Cancel" buttons.



The screenshot shows the Lyric ALPS Portal interface for entering client information. At the top, it says "Lyric". Below that, it says "Lyric client". There are input fields for "First name" (John), "Last name" (Doe), "Lyric ID" (P00008766), "Date of birth" (-), and "Gender" (-). A blue information icon with a text box below it says: "Please contact Lyric customer support if the client's Lyric ID is unknown." At the bottom right, there are "Continue" and "Cancel" buttons.

Tourist case: Your client has been refit by a different Authorized Lyric Partner

If your client has been refit by a different Lyric Provider, you will be alerted that there is a mismatch between ALPS (tourist fitter) and Phonak Target (your fitting). This indicates that the serial number has been changed.

By clicking [Continue] the system will synchronize the newest serial number.

Lyric

Data conflict between ALPS and Phonak Target

First nameJohn

Last nameDoe

Lyric IDP000087661

Device data in ALPS and Phonak Target do not match:

Information

Right

Left

Subscription

Trial (20 days left)

Trial (30 days left)

ALPS

1406NY0YG (Tourist exchange)

No serial number

Target

No serial number

No serial number

Action

Serial number will be updated in Phonak Target

No action required

Press "Cancel" to resolve the conflicts manually.

Continue

Cancel

If the tourist Lyric fitter also changed the fitting parameters, a fitting parameter mismatch will be shown. You can solve the conflict by choosing the preferred option.

Please note that the ALPS data is the most up to date information from the tourist fitter and Phonak Target data is the data you entered before the tourist fitting.

Lyric

Fitting parameter mismatch!

Fitting parameters stored in ALPS and Phonak Target do not match.

Right

Left

Resolve conflict

Use Phonak Target data

Last updated 16.06.2014 13:57

Use ALPS data (newer)

Last updated 16.06.2014 14:05

Updated by tourist fitter.

Fitting parameters are identical

Clicking "Continue" will resolve the conflict.

Press "Cancel" to resolve the conflict manually.

Continue

Cancel


Renew a subscription

As soon as a subscription is active (no longer in trial) it can be renewed by clicking on [Renew subscription].

When the subscription is three months from expiration, an orange [Renew subscription] button will appear.

Define which type of subscription length you want to create in ALPS for your client and follow the instruction wizard.

After renewing a cancelled/expired subscription (trial or active), a new trial will only be permitted after 60 days. A subscription renewed within 60 days of the expiration date will automatically be an active subscription.




Hearing aids		Lyric configuration
Subscription		
Status	Active (61 days left)	Cancel subscription Renew subscription ALPS links
Start date	01.08.2016	
End date	01.08.2016	
Trial end date	01.08.2016	
Remaining devices	7 (in subscription)	

Cancel of an active subscription

If there is a need to cancel an active subscription, you can do so by clicking [Cancel subscription] on the Lyric configuration screen.

The system will make you aware that you are about to cancel an active subscription. To proceed, choose the corresponding cancellation reason and click Continue. The cancellation request will be sent to Phonak.



Offline fitting and synchronization

There might be occasions where you have a client set up via Phonak Target/ALPS, however, due to various reasons (such as temporarily no internet connection or computer problems) you will have one or more sessions where you have to work offline. This means that the new data you have in Phonak Target is not synchronized with ALPS.

Removing a device in offline mode works the same way as in removing a device in online mode.

You are now asked to select a new device.

Enter the serial number as well as the code that is shown on the Lyric package. Enter the measured insertion depth. If needed, the fitting parameters can be changed.

Lyric configuration

Remove Lyric device

Serial number1633H00L5

Removed by

Provider

Patient self-removed

Removal date

Today

Previously

Removal reason

Discomfort

Additional details (optional)...

Ear canal observation

OK for immediate refit

Additional ear canal observation (optional)...

Not connected to ALPS. All changes will be synchronized later.

Remove

Cancel

Lyric configuration

Successfully removed

Device successfully removed.

Select new device

Close

Lyric configuration

Enter device information

Serial number1633H00L8

CodeE396

ModelLyric3


Seal sizeXL

Not connected to ALPS. All changes will be synchronized later.


Select

Cancel

Once the internet connection problem is resolved and you are able to use Phonak Target/ALPS online again, you will be notified by an orange exclamation point [Fitting not saved to ALPS] message for each patient that was fit offline.




As soon as Phonak Target has started connecting into ALPS again, it will automatically synchronize the offline sessions.



A green check mark will indicate that the synchronization has been successful.


After an offline fitting, logs on your desktop are available.



ALPS home and reports

There are some actions that require a direct link to ALPS:


- Open ALPS
- Client details
- Device details
- Print Lyric ID card
- Reports



Alternatively, ALPS can also be opened via Phonakpro.com.

On the [Home Screen] you have an overview of:

- Inventory Audit
- Trials about to expire within 7 days
- Subscriptions about to expire within 90 days
- Devices to expire within 120 days



- Print Lyric ID
- Subscription of that client
- Renew Subscription(s)
- Remove Devices (needed when cancelling the subscription)
- Client details (such as fitting parameter)

File name ↑	Last name	Lyric ID	Patent Status	Patent account number
05-Sep-2021	05-Sep-2021	P000221618	1 year	9020071645
1407	TC	P000101808	4 year	9020071645
1407 changed	TC	P000101808	4 year	9020071645
7.1.9	Bessaid	P000221648	1 year	9020071645
ss	ss	P000221629	1 year	9020071645
sss	sss	P000221686	1 year	9020071645
John	John	P000200596	2 year	9020071645
ALPS	Tes	P000221620	1 year	9020071645
Alps	Ljivi	P000221643	1 year	9020071645
Alps	Ljivi	P000221617	2 year	9020071645
alps	alps	P000221642	1 year	9020071645
ANDI	ANDI	P000221696	1 year	9020071645
Amayazawa39750101	Amayazawa39750101	P000221619	1 year	9020071645

- Terms and Conditions
- Cancel Subscription
- Device History

Subscription

2000568483

Subscription no

2000568483

Client

CHN

Terms and Conditions

Subscription plan

1-yr subscription

Payment method

Trade payment

Billing Day

Contracture date

Start

Total

Unit

8

Start and end date

18/7/2021

Start date

18/7/2021

End date

Remaining duration

7

Cancel Subscription

Device History

Serial number

State

Subscription date &

Renewal date

Renewal Reason

Warning Duration


11/07/2021


In use

8/7/2021

View All

- Device State
- Device Code, Warranty & Use By Date
- Subscription/Insertion
- Removal and Return Information

 Device 2132NV6TN	
Status	Short description
In use	Lyric (X8)
<hr/>	
Account	Warranty date
Lyric-Testfonds #3	10/7/2021
Social number	Use by
2132NV6TN	8/13/2022
Device Code	
F306	
<hr/>	
Subscription / Insertion	
Subscription	Insertion date
2020569483	9/7/2021
<hr/>	
Removal	
Removal Reason	Removal date
<hr/>	
Return	
Return reason	Return Date

Phonak Lyric™					
Home		Client	Stock	Subscriptions	Device Return
<div>  Report: Accounts with Devices Available Stock Detailed </div>					
Total Records					
18					
Device size ↑	Device: Serial number	Material number	Short description	Use by	
<input type="checkbox"/> L (6)	2132NY941	078-0064	lyricnl (L)	8/13/2022	
	2132NY943	078-0064	lyricnl (L)	8/13/2022	
	2132NY946	078-0064	lyricnl (L)	8/13/2022	
	2132NY948	078-0064	lyricnl (L)	8/13/2022	
	2132NY943	078-0064	lyricnl (L)	8/13/2022	
	2132NY94C	078-0064	lyricnl (L)	8/13/2022	
<input type="checkbox"/> M (3)	2132NY84V	078-0063	lyricnl (M)	8/13/2022	
	2132NY8C0	078-0063	lyricnl (M)	8/13/2022	
	2132NY8C4	078-0063	lyricnl (M)	8/13/2022	
<input type="checkbox"/> S (3)	2132NY77Y	078-0062	lyricnl (S)	8/6/2022	
	2132NY7KA	078-0062	lyricnl (S)	8/6/2022	
	2132NY7KC	078-0062	lyricnl (S)	8/6/2022	
<input type="checkbox"/> XL (6)	2132NY9WK	078-0065	lyricnl (XL)	8/6/2022	
	2132NY9WE	078-0065	lyricnl (XL)	8/6/2022	
	2132NY9WF	078-0065	lyricnl (XL)	8/6/2022	
Row Counts	<input checked="" type="checkbox"/>	Detail Rows	<input checked="" type="checkbox"/>	Subtotals	<input type="checkbox"/>
Grand Total		<input checked="" type="checkbox"/>			

Under [Subscriptions], all subscriptions with all states can be found. When working in a satellite office that is connected in the system, then all clients are visible but separated by account number. That list can be exported to Excel.

Phonak account number	Account Name	Subscription no.	Split ID	First name	Last name	Start date	End date	Side	Status	Subscription type	Payment model
2132NY7JY	Lyric Softwax R1 (75)	2132NY7JY	1	Lyric	Softwax	10/10/2021	10/10/2022	R	True	1 yr subscription	Nearly upfront
2132NY7JY	Lyric Softwax R1 (75)	2132NY7JY	1	Lyric	Softwax	10/10/2021	10/10/2022	R	True	1 yr subscription	Nearly upfront
2132NY7JY	Lyric Softwax R1 (75)	2132NY7JY	1	Lyric	Softwax	10/10/2021	10/10/2022	R	True	1 yr subscription	Nearly upfront
2132NY7JY	Lyric Softwax R1 (75)	2132NY7JY	1	Lyric	Softwax	10/10/2021	10/10/2022	R	True	1 yr subscription	Nearly upfront
2132NY7JY	Lyric Softwax R1 (75)	2132NY7JY	1	Lyric	Softwax	10/10/2021	10/10/2022	R	True	1 yr subscription	Nearly upfront
2132NY7JY	Lyric Softwax R1 (75)	2132NY7JY	1	Lyric	Softwax	10/10/2021	10/10/2022	R	True	1 yr subscription	Nearly upfront
2132NY7JY	Lyric Softwax R1 (75)	2132NY7JY	1	Lyric	Softwax	10/10/2021	10/10/2022	R	True	1 yr subscription	Nearly upfront
2132NY7JY	Lyric Softwax R1 (75)	2132NY7JY	1	Lyric	Softwax	10/10/2021	10/10/2022	R	True	1 yr subscription	Nearly upfront
2132NY7JY	Lyric Softwax R1 (75)	2132NY7JY	1	Lyric	Softwax	10/10/2021	10/10/2022	R	True	1 yr subscription	Nearly upfront
2132NY7JY	Lyric Softwax R1 (75)	2132NY7JY	1	Lyric	Softwax	10/10/2021	10/10/2022	R	True	1 yr subscription	Nearly upfront

The [Device Return] section allows you to return devices to Phonak that are either

- New (overstock)
- Dead on arrival
- Wrong delivery
- Expired

Choose the devices that need to be returned to Phonak. Click on the respective device and choose a return reason. Click select and then print a return form that can be sent with a device to return (if needed).

Serial number	Device size	Side	Short description	Use by
2132NY7JY	S	R/L	Lyric4 (S)	8/6/2022
2132NY7KC	S	R/L	Lyric4 (S)	8/6/2022
2132NY7KA	S	R/L	Lyric4 (S)	8/6/2022
2132NY9WK	XL	R/L	Lyric4 (XL)	8/6/2022
2132NY9WD	XL	R/L	Lyric4 (XL)	8/6/2022
2132NY9WC	XL	R/L	Lyric4 (XL)	8/6/2022
2132NY9WG	XL	R/L	Lyric4 (XL)	8/6/2022
2132NY9WF	XL	R/L	Lyric4 (XL)	8/6/2022
2132NY9WE	XL	R/L	Lyric4 (XL)	8/6/2022
2132NY941	L	R/L	Lyric4 (L)	8/13/2022
2132NY94C	L	R/L	Lyric4 (L)	8/13/2022
2132NY943	L	R/L	Lyric4 (L)	8/13/2022
2132NY948	L	R/L	Lyric4 (L)	8/13/2022
2132NY946	L	R/L	Lyric4 (L)	8/13/2022
2132NY945	L	R/L	Lyric4 (L)	8/13/2022
2132NY8AY	M	R/L	Lyric4 (M)	8/13/2022
2132NY8C0	M	R/L	Lyric4 (M)	8/13/2022
2132NY8C4	M	R/L	Lyric4 (M)	8/13/2022

Under your name, there are the following sections:

- Home
- My Profile
- My Account
- Logout

Target Tester DE

- Home
- My Profile
- My Account
- Logout

[My Profile] allows you to change your language if needed.

My Settings

Location

Language: English

Time Zone: GMT+01:00 Central European Standard Time (CET)

Cancel Save

Under [My Account] you can choose your current location (when working in different locations).

Please choose your current location:

Your current account is:

833307642
Lyric Testbunde #3
Subgroup 1
Subgroup 10124

Please choose one account. Please account number

Account Name

Address

City

Zip code

☐

833307642

Lyric Testbunde #3

Markenstraße 1

Stuttgart

70372

Testbed

In case you would like to get a weekly e-mail about the trials, subscriptions and stock about to expire, you can fill in up to five e-mail addresses. Also, choose how many months before expiration of the subscriptions you'd like to get a notification. For satellite offices, a report with a full overview over all offices can also be requested by ticking the box.

Notification Settings

Account Name

Lyric Testbunde #3

Notification E-Mail 1

Notification E-Mail 2

Notification E-Mail 3

Notification E-Mail 4

Notification E-Mail 5

Expiration Notification Time

1 month

Notification Report Satellite Offices

☐

Save

The [Dashboard] gives you an indication about how your business is doing. Explanation is found in more details in the legend within the dashboard.

Trials 64 / New Subscriptions 67 / Renewals 0 / Attivo 3 / Trial Success Rate 4 % / TSR national 25.00 %

1. Trials with starting date Jan 1st - Dec 31st

2. Renewals with starting date Jan 1st - Dec 31st

3. Total active subscriptions

4. Total active subscriptions

5. Total active subscriptions

6. Total active subscriptions

7. Total active subscriptions

8. Total active subscriptions

The [Device Analysis Form] can be downloaded.

Forms to Download


The [Search] field allows you to search data such as subscription, client names, devices etc. directly.

Search...

Enhanced log files for support

If you encounter a problem with Phonak Target, enhanced logging can help our Target software support team to find out the root cause.

Please save the log files on your desktop and/or send them by e-mail to your Lyric contact.



Information and description of symbols and System Requirements

Information and the description of symbols and an overview of system requirements can be found in the Phonak Target User Guide

This user guide is applicable for Target 11.0 and later subversions of Target 11.0 fitting software. For previous versions of the fitting software user guide contact your local Phonak representative.



Sonova AG • Laubisrütistrasse 28 •
CH-8712 Stäfa • Switzerland



Sonova Deutschland GmbH
Max-Eyth-Str. 20
70736 Fellbach-Oeffingen • Germany



0459