

# Hearing Assistance Technologies

---

- Individual and environment considerations must be analyzed before determining the best style for any individual.
- All styles improve speech reception of speaker by minimizing background noise, high reverberation levels, distance from speaker and other interfering factors.
- Generally the closer the signal is delivered to the ear the greater the enhancement.
- All systems minimally contain a remote wireless transmitter/microphone and a receiver and/or speaker.

## Transmission Modes:

<b>Roger:</b>	Digital standard remote microphone technology from Phonak that uses 2.4GHz transmission from wireless microphones to small ear-level, telecoil or universal receivers.
<b>Advantages:</b>	Several receivers can be connected to the same microphone, no frequency interference, multiple microphones are available and will work together in multi-talker network. Uses an adaptive algorithm to provide benefit at louder noise levels.
<b>Limitations:</b>	Requires a Roger microphone and a Roger receiver.
<b>FM:</b>	Remote microphone technology that uses radio frequency transmission from wireless microphone(s) to FM receivers.
<b>Advantages:</b>	Several receivers can be connected to the same microphone
<b>Limitations:</b>	Prone to transmission interference, requires frequency management

## Type of Microphones

### Omni Microphones

<b>Advantages:</b>	Sensitive to sound in all directions
<b>Limitations:</b>	May transmit more noise, especially in large crowds

### Directional Microphones

<b>Advantages:</b>	Fixed or adaptive microphone focuses on one signal minimizing amplification of surrounding noise. Some units allow user to switch between omni and directional microphone positions
<b>Limitations:</b>	If positioned incorrectly, the speaker's voice may be weak. For user selected option, multiple settings may be confusing to user. Automatic functionality is sometimes available.

## Styles of Remote Microphones

<b>Lavalier:</b>	Microphone and transmitter in one unit
<b>Advantages:</b>	Easy to hang around the neck
<b>Limitations:</b>	May pick up clothing noise. Signal may be reduced when head turned from microphone (especially if using directional microphone). Some models are cumbersome.
<b>Lapel Microphone:</b>	Connects to body-worn transmitter
<b>Advantages:</b>	Microphone clips easily onto clothes
<b>Limitations:</b>	May pick up clothing noise. Signal may be reduced when head turned from microphone (especially if using directional microphone)
<b>Conference:</b>	Microphone and transmitter in one unit; sits on tabletop
<b>Advantages:</b>	Omni-directional microphone picks up multiple speakers around the table
<b>Limitations:</b>	Speaker may be difficult to hear due to distance from microphone
<b>Boom Microphone:</b>	Connects to body-worn transmitter
<b>Advantages:</b>	Microphone stays close to speaker's mouth providing the most consistent signal.
<b>Limitations:</b>	May be uncomfortable or distracting to user.
<b>Ear Level Receivers</b>	
<b>Design Integrated:</b>	Design integrated system containing a Roger or FM receiver for signal from a compatible wireless remote microphone.
<b>Advantages:</b>	Small, discreet and always attached.
<b>Limitations:</b>	Receivers are often not compatible with multiple hearing aids.
<b>Universal Receiver : (with Audio Shoe)</b>	Adapter that interfaces Roger to HA/CI/BAHA to receive signal from wireless remote microphone.
<b>Advantages:</b>	Universal receiver with audio shoe interfaces to a variety of personal HAs/CIs & BAHA. Small in size.
<b>Limitations:</b>	Additional adapter/audio shoe connection increases susceptibility of signal to intermittency or malfunction
<b>Stand-alone:</b>	Wireless Roger receiver, may include signal processing strategies to enhance speech reception.
<b>Advantages:</b>	Provides Roger advantage to individuals with normal to near normal hearing sensitivity who require additional speech enhancement due to problems such as hearing, listening, processing and/or attention, Small, inconspicuous, Non-occluding style permits access to sounds in addition to those from remote microphone
<b>Limitations:</b>	Cannot be used without a microphone/transmitter, Not a replacement for a hearing aid

## Telecoil Receivers

**Personal Neck Loop:** Induction loop is worn around the neck

**Advantages:** Operates through hearing aid t-coil, Inexpensive, Easy to operate,

**Limitations:** Signal interruptions can occur if not properly oriented

**Wide Area Loop:** Large room or area is looped with a special wire.

**Advantages:** Easily installed in a room, classroom, or car. Operates through hearing aid t-coil. Inexpensive, Easy to operate,

**Limitations:** Requires strong hearing aid t-coil, Signal interruptions can occur if not properly oriented

**3D Mat:** Series of loops are contained in floor that could be installed under large rugs or carpeting.

**Advantages:** Consistent signal transmission, Operates through hearing aid t-coil, Durable (less maintenance)

**Limitations:** May be too costly for home use, Permanent site if installed under carpet, Requires good hearing aid t-coil

## Classroom Audio Distribution Systems

**Desktop speaker:** Speaker's voice is delivered via wireless transmitter to speaker placed strategically near student or on student's desk.

**Advantages:** Provides slightly enhanced speech when ear level option is not possible, Portable

**Limitations:** Signal enhancement less than ear level options,

**SoundField:** Speaker's voice is delivered from remote microphone/transmitter to a loud speaker strategically located in classroom or listening space (e.g. church, community room, outdoor theater). Soundfields can use a variety of transmission signal including the Roger Dynamic Soundfield, Radio Frequency Soundfields and Infrared Soundfield systems.

**Advantages:** Easy to operate, Benefits all listeners with normal hearing, No equipment to wear

**Limitations:** Not a replacement for hearing aid/CI/Baha for a user with hearing loss