

# Using the Roger X / MLx Audio Checker

Reference Sheet

**roger**



School: \_\_\_\_\_

Teacher: \_\_\_\_\_

Roger contact: \_\_\_\_\_

Phone: \_\_\_\_\_

Date: \_\_\_\_\_

## Instructions for use

1. The MLx Audio Checker should be charged before use.
2. Plug Roger X receiver in the front plug. It will only fit one way.
3. Plug the headphones into the bottom of the MLx Audio Checker. Turn on the Audio Checker.
4. Turn on the Roger inspiro microphone. You should be able to then hear the input sound from the microphone through the headphones.
5. Should you not hear any sound, please check for possible solutions:
  - A. Press **Connect** on Roger inspiro. If it now works, the receiver was not in the network.
  - B. Check to make sure the microphone is not muted.
  - C. If none of the above work, please contact your Roger contact named above.

