

Using Roger to hold a private conversation

Reference Sheet



If you are the only teacher in a class and you need to speak privately with a child who is wearing a Roger receiver, without other students hearing, you can quickly create a new network to talk with this child. This can be especially useful when the classroom network includes a soundfield system.

Create a new network – to cut off your microphone from other students and only transmit to specific child/children

1. Press the **OK** button on the Roger inspiro.
2. Scroll until **Network** is highlighted. Press **OK**.
3. Scroll until **Create NewNet** is highlighted. Press **OK**.
4. Hold Roger inspiro 4 inches (10 cm) from child's Roger receiver and press **Connect**.
5. Repeat with any other receiver you would like to connect to this new network.

Rejoin the class – to put the Roger inspiro and children back into the main classroom network

1. Press the **OK** button.
2. Scroll until **Network** is highlighted. Press **OK**.
3. Hold Roger inspiro 4 inches (10 cm) from the top of the DigiMaster loudspeaker or the receiver of any child that was not put in the new network.
4. Scroll until **JoinNet** is highlighted. Press **OK**. The Roger inspiro will prompt for next steps to join the network.
5. Reconnect the children who were in private conversation to the classroom network by holding the Roger inspiro close to a child's receiver and pressing **Connect**. Repeat for all receivers that were in the private conversation.
6. If a Roger WallPilot is installed in the classroom, children who were in the temporary network can alternatively **Connect** by walking within 3-6 feet (1-3 m) of the WallPilot.