

# Using the Roger DigiMaster X

## Reference Sheet





School:	Teacher:
Roger contact:	Phone:
Date:	Serial number:

## Connecting

### Initial set-up








1. Connect the Roger DigiMaster X to a power source.
2. Connect the Roger DigiMaster X to the room amplification system.
3. Hold Roger inspiro less than 10 cm (4 inches) from the Roger DigiMaster X and press **Connect**.

### Daily use

1. Press the power button  on the Roger DigiMaster X for 2 seconds to power on. Turn on Roger inspiro. Speak into the Roger inspiro microphone.
2. Press the power button  on the Roger DigiMaster X for 2 seconds to power off when finished.



## Roger DigiMaster X indicator light

-  Solid red – Roger DigiMaster X is off
-  Solid blue – Standard operation mode
-  Blinking blue – Roger microphone is off or out of range
-  Solid green – Microphone is muted
-  Green to red – Volume is adjusting +/-
-  Blinking yellow – Software update
-  Blinking red – Call your Roger contact.