Phonak Insight

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Quellenangaben

1. World Health Organisation guidelines on physical activity and sedentary behaviour (2020). Genf: Weltgesundheitsorganisation (WHO). Lizenz: CC BY-NC-SA 3.0 IGO
2. Kuo, P. L., Di, J., Ferrucci, L., & Lin, F. R. (2021). Analysis of hearing loss and physical activity among US adults aged 60-69 years. *JAMA Network Open*, 4(4), e215484-e215484.
3. Tsimpida, D., Kontopantelis, E., Ashcroft, D., & Panagioti, M. (2019). Socioeconomic and lifestyle factors associated with hearing loss in older adults: a cross-sectional study of the English Longitudinal Study of Ageing (ELSA). *BMJ Open*, 9(9), e031030.
4. Wells, T. S., Nickels, L. D., Rush, S. R., Musich, S. A., Wu, L., Bhattarai, G. R., & Yeh, C. S. (2020). Characteristics and Health Outcomes Associated With Hearing Loss and Hearing Aid Use Among Older Adults. *J Aging Health*, 32(7-8), 724-734.
5. Lin, T. C., Yen, M. & Liao, Y. C. (2019). Hearing loss is a risk factor of disability in older adults: A systematic review. *Arch Gerontol Geriatr*, 85, 103907. Aufrufbar unter: <https://doi.org/10.1016/j.archger.2019.103907>, Stand 1. Dezember 2021.
6. Appleton, J. (2020). AutoSense OS™ 4.0 provides better speech intelligibility and reduced listening effort Phonak Field Study News. Quelle: www.phonakpro.com/evidence, Stand Dezember 2020.
7. World Health Organization (WHO). (1948). Constitution of the World Health Organization. Genf, Schweiz: WHO.
8. Caduff, A., Feldman, Y., Ishai, P. B., & Launer, S. (2020). Physiological monitoring and hearing loss: toward a more integrated and ecologically validated health mapping. *Ear and Hearing*, 41, 120S-130S.
9. Holman J. A., Hornsby B. W. Y., Bess F. H., Naylor G. (2021). Can listening-related fatigue influence well-being? Examining associations between hearing loss, fatigue, activity levels and well-being. *Int J Audiol*. 4:1-13.
10. Buyl, R., Beogo, I., Fobelets, M., Deletroz, C., van Landuyt, P., Dequanter, S., Gorus, E., Bourbonnais, A., Bourbonnais, A., Giguère, A., Giguère, A., Giguère, A., Lechasseur, K., Gagnon, M. P., Gagnon, M. P., & Gagnon, M. P. (2020). E-Health interventions for healthy aging: A systematic review. *Systematic Reviews*, 9(1), 1-16.
11. Marcolino, M. S., Oliveira, J. A. Q., D'Agostino, M., Ribeiro, A. L., Alkmim, M. B. M., & Novillo-Ortiz, D. (2018). The impact of

mHealth interventions: Systematic review of systematic reviews. JMIR MHealth and UHealth, 6 (1).

12. 2018 Physical Activity Guidelines Advisory Committee. 2018 Physical Activity Guidelines Advisory Committee Scientific Report. Washington, DC: U.S. Department of Health and Human Services, 2018. Aufrufbar unter: <https://health.gov/our-work/nutrition-physical-activity/physical-activity-guidelines/current-guidelines/scientific-report#:~:text=The%202018%20Physical%20Activity%20Guidelines%20Advisory%20Committee%20submitted,edition%20of%20the%20Physical%20Activity%20Guidelines%20for%20Americans>, Stand 1. Dezember 2021.

13. Stockwell, S., Schofield, P., Fisher, A., Firth, J., Jackson, S., Stubbs, B., and Smith, L. (2019). Digital behavior change interventions to promote physical activity and/or reduce sedentary behaviour in older adults: A systematic review and meta-analysis. *Experimental Gerontology*, 120, 68-87.

14. Simblett, S., Greer B., Matcham F., Curtis, H., Polhemus, A., Ferrao, J., Gamble, P., Wykes, T. (2018). Barriers to and facilitators of engagement with remote measurement technology for managing health: systematic review and content analysis of findings. *Journal of Medical Internet Research*, 20 (7):e10480.

15. Kappen, D. L., Mirza-Babaei, P., & Nacke, L.E. (2019). Older Adults' Physical Activity and Exergames: A Systematic Review, *International Journal of Human-Computer Interaction*, 35(2), 140-167.

16. Martinho, D., Carneiro, J., Corchado, J. M., & Marreiros, G. (2020). A systematic review of gamification techniques applied to elderly care. *Artificial Intelligence Review*, 53(7), 4863-4901.

17. He, D., Winokur, E. S., Heldt, T. & Sodini, C. G. (2010). The ear as a location for wearable vital signs monitoring. Annual International Conference of the IEEE Engineering in Medicine and Biology Society. Conference Proceedings. 6389-6392. Aufrufbar unter: https://www.researchgate.net/publication/49627895_The_Ear_as_a_Location_for_Wearable_Vital_Signs_Monitoring, Stand 1. Dezember 2021.

18. Blazek, V., Venema, B., Leonhardt, S. and Blazek, P. (2018). Customized optoelectronic in-ear sensor approaches for unobtrusive continuous monitoring of cardiorespiratory vital signs. *International Journal of Industrial Engineering and Management*, 9(4), 197-203.

19. Budidha, K. & Kyriacou, P. A. (2014). The human ear canal: investigation of its suitability for monitoring photoplethysmographs and arterial oxygen saturation. *Physiological Measurement*, 35(2), 111-128.

20. Bunn, J., Wells, E., Manor, J., & Webster, M. (2019). Evaluation of Earbud and Wristwatch Heart Rate Monitors during Aerobic and Resistance Training. *Int J Exerc Sci*, 12(4), 374-384.

21. Hartmann V, Liu H, Chen F, Qiu Q, Hughes S and Zheng D (2019). Quantitative Comparison of Photoplethysmographic Waveform Characteristics: Effect of Measurement Site. *Front. Physiol.* 10:198, 1-8.

22. Budidha, K. & Kyriacou, P. A. (2019). Photoplethysmography for Quantitative Assessment of Sympathetic Nerve Activity (SNA) During Cold Stress. *Front Physiol*, 9(1863), 1-10.

23. Fallow, B. A., Tarumi, T., and Tanaka, H. (2013) Influence of skin type and wavelength on light wave reflectance. *International Journal of Clinical Monitoring and Computing*. 27: 313-317.

24. Lister, S., Hofland, J., & Grafton, H. (Eds.). (2020). *The Royal Marsden manual of clinical nursing procedures*. John Wiley & Sons.

25. Fox, K., Borer, J., Camm, A., Danchin, N., Ferrari, R., Sendon, J., Steg, P., Tardif, J.C., Tvasz, L., & Tendera, M. (2007). Heart Rate Working Group. Resting heart rate in cardiovascular disease. *Journal of the American College of Cardiology*, 50(9), 823-830.

26. Larsson, S. C., Drca, N., Mason, A. M., & Burgess, S. (2019). Resting Heart Rate and Cardiovascular Disease: Mendelian Randomization Analysis. *Circulation: Genomic and Precision Medicine*, 12(3), e002459.

27. Liguori, G., & American College of Sports Medicine. (2020). ACSM's guidelines for exercise testing and prescription. Lippincott Williams & Wilkins.

28. Tudor-Locke C., Craig C.L., Brown W.J., Clemes, S. A., De Cocker, K., Giles-Corti, B., Hatano, Y., Inoue, S., Matsudo, S., Mutrie, N., Oppert, J-M., Rowe, D., Schmidt, M., Schofield, G., Spence, J., Teixeira, P., Tully, M. & Blair, S. (2011). How many steps/day are enough for adults? *Int J Behav Nutr Phys Act*; 8:79.

29. Friedman, D, & Madden, K. (2019). Older people. Brukner P, Khan K. Brukner P, & Khan K(Eds.), Eds. Peter Brukner, and Karim Khan. Brukner & Khan's Clinical Sports Medicine: The Medicine of Exercise, Volume 2, 5e. McGraw-Hill. Aufrufbar unter: <https://csm.mhmedical.com/content.aspx?bookid=2667§ionid=224525454>, Stand 1. Dezember 2021.

30. World Health Organisation, Food and Agriculture Organization of the United Nations, United Nations University (2004). Human energy requirements. Report of a joint FAO/WHO/UNU Expert Consultation, Rome, Italy, 17-24 October 2001.

31. Sarant, J., Harris, D., Busby, P., Maruff, P., Schembri, A., Lemke, U., & Launer, S. (2020). The effect of hearing aid use on cognition in older adults: Can we delay decline or even improve cognitive function? *Journal of Clinical Medicine*, 9(1), 254.

Autor

Jane Woodward, MSc
Audiologie Managerin, Phonak HQ, Switzerland

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Elizabeth Stewart, AuD, PhD
Research Audiologist, Phonak Audiology Research Center (PARC), USA

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