
Power supply

Battery:	Lithium-Ion Polymer
Capacity:	700 mAh
Operating time:	Stand-by 240 hours (version > 2.17.0.0)
Charging unit:	Type: DC, stabilized Secondary voltage: 5 VDC Max. current: 650 mA
Fast-charge:	80% charged within an hour Fully charged in 1.5 hours

The battery cannot be over-charged. It may become warm during charging.

Important notice: It is recommended to use only the original Phonak charging unit.

Radio information

Antenna type:	Integrated ferrite antenna
Operation frequency:	10.6 MHz
Occupied bandwidth (99% BW):	512.0 kHz
Field Strength:	- 11.5 dB μ V/m in 10 m

Operating range:

Remote control operating range is about 100 cm (40") to hearing instrument. Status readout operating range is about 50 cm (20") from hearing instrument.

Standards

Europe

Radio:	EN 300 330
EMC:	EN 60601-1-2
El. Safety:	EN/IEC 60601-1

USA

FCC ID:	KWC-MYPILOT1
	FCC Part 15 Subpart C

Canada

IC:	2262A-MYPILOT1
	RSS-210 and ICES-003