

Hear today, hear tomorrow



Our long-standing commitment to children with hearing loss has inspired us to create Nios micro: a dedicated, state-of-the-art and attractive solution for children with a mild to moderately-severe hearing loss. Its unique benefits give invaluable support to children on the way to reaching their full potential.

More Hearing

All the audibility and clarity needed to achieve the best speech and language outcomes from day one.

Unlimited Connectivity

FM and wireless connectivity in all essential environments, for now and later.

Full Confidence

Safety, dependability and confidence in an attractively small design for a future without limits.

Nios micro is available in two price categories and in 13 colors. Try Nios micro and let your child experience the difference. www.phonak.com

Life is on

We are sensitive to the needs of everyone who depends on our knowledge, ideas and care. And by creatively challenging the limits of technology, we develop innovations that help people hear, understand and experience more of life's rich soundscapes.

Interact freely. Communicate with confidence. Live without limit. Life is on.

Your hearing care professional:

www.phonak.com

028-0473-02/V1.00/2009-03/mwlp Printed in Switzerland © Phonak AG All rights reserved

Nios micro

Dedicated to young lives

More Hearing. Unlimited Connectivity. Full Confidence.



PHONAK life is on

More Hearing



In tune with children's needs

Nios micro offers children with a mild to moderately-severe hearing loss the sound awareness and speech understanding they need to discover all the subtle, enriching sounds that make up their world.

Exceptional audibility and clarity

Packed with unique and intelligent features, such as SoundRecover and WhistleBlock Technology, Nios micro offers exceptional audibility and clarity.

▪ Breakthrough in speech and language development

The unique SoundRecover feature ensures audibility of high pitched sounds, beyond what is available with conventional amplification. It plays a crucial role in the development of speech and language, giving children with hearing loss a greater chance to keep up with their hearing friends.

▪ No feedback, more sound

WhistleBlock Technology effectively eliminates feedback (whistling), so children with a hearing loss can benefit from increased amplification, wherever they are, whatever they do.



Unlimited Connectivity

Access to the world of wireless communication

Interaction with modern communication and entertainment systems is increasingly important for children as they grow up. With Nios micro, they will benefit from the unique ability to connect to all wireless means of communication.

Design-integrated Dynamic FM solution¹

Dynamic FM provides exceptional clarity when children need to communicate across distance, in noisy environments and when echo is present. FM is ideal at home, in the car, in day-care, classroom situations and more. The elegant and convenient Dynamic FM solution helps children to reach their full potential.

Opening up the world of wireless communication²

iCom provides seamless connectivity to mobile phones, computers, laptops, MP3 players and TV. It is the ideal companion for children in the modern world of communication.



¹with design-integrated FM receiver ML12i



²iCom

Full Confidence



Safe and easy to operate

Nios micro offers unique options for parents' and caregivers' peace of mind.

Tamperproof battery lock

Prevents small children from accessing the battery.

Stick 'n Stay

Designed especially with babies and small children in mind, Stick'n Stay is an invisible, double-sided, skin-friendly tape that holds the hearing instruments securely in place.

Reassurance at the touch of a button³

Thanks to iView, it has never been easier for parents and other caregivers to check the functionality of their children's hearing systems. Battery status, volume level and current program on both hearing instruments are easily monitored on the color LCD display.



³iView