

junior

The Junior Collection

Great Solutions for little ears

Based on the CORE collection



Life is on

We are sensitive to the needs of everyone who depends on our knowledge, ideas and care. And by creatively challenging the limits of technology, we develop innovations that help people hear, understand and experience more of life's rich soundscapes.

**Interact freely. Communicate with confidence.
Live without limit. Life is on.**

www.phonak.com

V1.00/2008-09/ML © Phonak AG. All rights reserved



When developing the newest technologies, Phonak also always has little ears in mind. This results in real benefits for children, parents, caregivers, and the professionals working with them. CORE offers sophisticated pediatric solutions for children of all ages and with all degrees of hearing loss.

junior

Hearing instruments and FM systems

An attractive range of hearing instruments and FM solutions as well as a wide selection of useful accessories offer great solutions for little ears.

Hearing loss \ Segment	Mid-range	Value
Mild to moderately-severe	Versáta micro	Certéna micro
Moderately-severe to profound	Naida Junior SP V	Naida Junior SP III
Severe to profound	Naida Junior UP V	Naida Junior UP III

New – ideally suited for children with a mild to moderately-severe hearing loss, the Versáta and Certéna micro models offer a perfect set of features and accessories.

All instruments, regardless of size, can be easily coupled to sophisticated FM solutions either with a design-integrated or universal receiver. The groundbreaking Dynamic FM with the new inspiro transmitter is the new benchmark for pediatric applications.

Safety, security and fun

- A unique tamperproof battery solution to counteract inquisitive little hands
- Flexible volume/program changing options to suit the age and abilities of the child
- The mini hook ensures a comfortable and secure fit behind small ears
- The flexible Volume Control Block option means parents and caregivers can easily change the volume while the child cannot
- Child friendly colors are complemented by new animal patterns and totally transparent housings
- iCom communication interface, ideally suited for connecting kids to computers, MP3 players, mobile phones and telephone, iCom ensures that no one is left behind

iView – reassurance at the touch of a button

Battery life, volume level and current program on both hearing instruments can be easily monitored on a color LCD display. Clear information at the touch of a button provides the long-awaited reassurance professionals and parents have been asking for.



iPFG Junior Mode

From iPFG 2.0 and above Phonak is proud to support all products with the unique pediatric fitting mode. Each new release of iPFG brings with it exciting additions to the Junior Mode.

Highlights

- Evidence-based software defaults for children aged 0–4 and 5–8 years
- FM+M as the start-up program
- DSL v5a as the default fitting formula
- Individual fitting philosophies supported with customizable defaults
- Extensive help section
- Printable, tailored materials for parents, caregivers and teachers including:
 - Hearing System Passport for parents and teachers, providing personalized information on the child's hearing loss, hearing instruments and FM system
 - Communication Tips and Strategies
 - SoundRecover for Naida models (non-linear frequency compression)
 - Diary of Observations enables parents to easily monitor their child's progress
 - Information on Speech, Language and Listening Development, offering information as well as practical tips

CableFree Fitting

Fitting children can be particularly challenging due to their movement and activity. iCube allows for wireless fitting creating a more relaxed fitting environment. It allows professionals to continue programming even if the child is moving around – at a distance of up to 10 meters (30 feet).



CORE offers the most comprehensive collection of hearing instruments and FM solutions. The flexible range of accessories and the dedicated junior fitting mode offer a perfect portfolio for little ears.