


Custom Earpiece Order Form

Receiver-In-Canal Models – Sound Delivery System (SDS) for Lumity, Paradise & Marvel

Tube Length

- Measure the tube length to select the correct Receiver or SlimTube
- Measure both ears for binaural fitting



Vent Styles

- We offer two vent solutions for SlimTip (not available for SlimTip Titanium AV):



Regular

- Classic vent solution
- Parallel vent (not modifiable in clinic)



Cavity


- Optimized for best fit rate
- Easy to enlarge by the fitter

	S Receiver Gain 47(2cc) • MPO 111(2cc)	M Receiver Gain 51(2cc) • MPO 114(2cc)	MAV Receiver Gain 50(2cc) • MPO 119(2cc)	P Receiver Gain 59(2cc) • MPO 122(2cc)	UP Receiver Gain 71(2cc) • MPO 132(2cc)
Custom Earpieces					
cShell Acryl (hard)	•	•		•	•
cShell Titanium (hard)	•	•		•	•
SlimTip Acryl (hollow hard)	•	•		•	
SlimTip Silicone (solid soft)	•	•		•	
SlimTip Titanium AV (hollow hard)			•		
Instant Fit – Domes					
Cap Dome	•	•		•	
Open Dome	•	•		•	
Vented Dome	•	•		•	
Power Dome	•	•		•	
Universal SlimTip			•		


Hearing Instrument Compatibility

The SDS earpiece portfolio is compatible with Lumity, Paradise and Marvel hearing aids. It is not backwards compatible to Belong or earlier generation hearing aids. ActiveVent (MAV) is only compatible with Lumity and Paradise rechargeable RICs.


S Receiver
Mild to moderate hearing loss




M & MAV Receiver
Mild to moderately-severe hearing loss



P Receiver
Moderate to severe hearing loss



UP Receiver
Moderately severe to profound hearing loss



cShell Acryl (hard)



cShell Titanium (hard)



SlimTip Acryl (hollow hard)



SlimTip Silicone (solid soft)



SlimTip Titanium AV (hollow hard)



ActiveVent Receiver (MAV)

Custom Earpiece Order Form

Date: ____ / ____ / ____
DD MM Year

Rush Order in 48hr
(\$48 charge) N/A for Titanium

PO#: _____

Attention to _____
(at Phonak): _____

STEP 1: General Information

Ship to Account Number:

Bill to Account Number:

Phone #:

Company Name:

Address:

Contact Name (Audiologist/Dispenser):

First Name:

Last Name:

Claim Type:

- AADL BCEHP CNESST Eastern Health Greenshield
 Manitoba Health Nisga'a ODSP REGIE Supp. Health
 VAC, DND, RCMP, ISC (NIHB) WCB _____ WS-BC WSIB
 Workplace NL, WS-NB, WCB-PEI Other: _____
(Please Specify)

Claim # (Required): _____

L	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
R	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
	250	500	1K	2K	3K	4K

Please include audiogram with all orders. Audiogram required for Acoustically Optimized Venting (AOV).

STEP 2: Order Information

A) I will be fitting the earpiece(s) to hearing instrument(s): L Model: _____ R Model: _____

B) Please send earpiece(s) & hearing instrument(s): L Model: _____ Colour: _____ R Model: _____ Colour: _____

Please also send me the following accessories: _____

STEP 3: Earpiece Order

S = Standard O = Optional

00 only available in SlimTube	Please indicate Tube length (00/0/1/2/3)		Receivers					BTE	Wax System					
	L	R	S Receiver	M Receiver	P Receiver	UP Receiver ¹	ActiveVent Receiver ²	SlimTube	CeruStop	Spring Wax Guard	Receiver Tube (Ext.)	HF3	HF4	No Wax Guard
<input type="checkbox"/> L <input type="checkbox"/> R SlimTip Acryl (hollow hard)	L	R	O	S	O			O	S	O	O			O
<input type="checkbox"/> L <input type="checkbox"/> R SlimTip Silicone (solid soft)	L	R	O	S	O			O						Required
<input type="checkbox"/> L <input type="checkbox"/> R SlimTip Titanium AV (hollow hard)	L	R					Required					Required		
<input type="checkbox"/> L <input type="checkbox"/> R cShell Acryl (hard)	L	R	O	S	O	O			S	O	O	O	O	O
<input type="checkbox"/> L <input type="checkbox"/> R cShell Titanium (hard)	L	R	O	S	O	O			Required					
<input type="checkbox"/> L <input type="checkbox"/> R cShell Acryl - XL (hard)	L	R	O	O	O	S			O	O	O	O	S	O

All SlimTip receivers, excluding ActiveVent, will be delivered with a Cerustop Wax Guard in the receiver spout. The SlimTip Acryl will be built with Cerustop as an additional wax protection. ActiveVent will be delivered with an HF3 wax guard in the receiver spout.

¹UP Receiver not available with Slim hearing aid
²ActiveVent Receiver is only available for Lumity & Paradise rechargeable

STEP 4: Product Options

Faceplate colour options:

- L R Beige
L R Clear
L R Cocoa
L R Brown
L R Black

cShell Titanium only available with black faceplate

Shell colour options:

- L R Beige (Std. on cShell)
L R Clear (Std. on SlimTip)
L R Blue (L) & Red (R) (Acryl 4.0 only)
L R Brown (Silicone only)
L R Rose Tint (Silicone only)
L R Black (Sensor cShell only)

cShell Titanium & SlimTip Titanium AV only available in Titanium Gray

Vent options:

- L R Select-a-Vent (SAV)
L R Open Vent (SlimTip only)
L R AOV (Agm. required)
L R Pressure Vent
L R No Vent*

*No Vent is only available for SlimTip Titanium AV

Other options:

- L R Removal Filament
L R Extended Removal Filament
L R Soft Coat
L R Canal Lock*
L R Helix Lock (Acrylic only)*
L R Skeleton Lock*

*Same colour as shell

Additional Information: