

**Step 1: General Information & Audiogram**

Rush order in 48h (\$40.00 charge)  
(N/A for Titanium)

Date: \_\_\_/\_\_\_/\_\_\_  
          DD        MM        Year

PO#: \_\_\_\_\_

Attention to (at Phonak): \_\_\_\_\_

**Ship to Account Number:**

Phone#:  
Company Name:  
Address:  
Contact Name (Audiologist/Dispenser):

SHIP TO

**Bill to Account Number:**

Phone#:  
Company Name:  
Address:

BILL TO

Last Name: [ ] [ ] [ ] [ ] [ ] [ ] [ ] [ ] [ ] [ ] [ ] [ ]

First Name: [ ] [ ] [ ] [ ] [ ] [ ] [ ] [ ] [ ] [ ] [ ] [ ]

- AADL BCEHP CNESST Eastern Health Greenshield
- Manitoba Health Nisga'a ODSP REGIE Supp.Health
- VAC,DND,RCMP,ISC(NIHB) WCB-AB WCB-MB WCB-NS
- WCB-SASK WCB-Yukon Workplace NL,WS NB, WCB PEI
- WS-BC WSIB Other: \_\_\_\_\_

Claim # (Required): \_\_\_\_\_

**Audiogram (Please include with ALL orders)**

|   |     |     |    |    |    |    |
|---|-----|-----|----|----|----|----|
| L |     |     |    |    |    |    |
| R |     |     |    |    |    |    |
|   | 250 | 500 | 1K | 2K | 3K | 4K |

**Speech Scores**

dB HL     dB SPL

|     | Right | Left  | Binaural |
|-----|-------|-------|----------|
| SRT | _____ | _____ | _____    |
| MCL | _____ | _____ | _____    |
| UCL | _____ | _____ | _____    |

Ear Texture:  Soft  Avg.  Firm

**Step 2: Product**

(S) = Standard

|                  | Side   | IIC                        | CIC                        | MC                         | ITC                      | HS                         | FS                       |
|------------------|--|----------------------------|----------------------------|----------------------------|--------------------------|----------------------------|--------------------------|
| Virto M-Titanium | <input type="checkbox"/> L<br><input type="checkbox"/> R | <input type="checkbox"/> S | <input type="checkbox"/>   | N/A                        | N/A                      | N/A                        | N/A                      |
| Virto M-10 NW 0  | <input type="checkbox"/> L<br><input type="checkbox"/> R | <input type="checkbox"/>   | <input type="checkbox"/> S | <input type="checkbox"/>   | <input type="checkbox"/> | <input type="checkbox"/>   | <input type="checkbox"/> |
| Virto M-312 NW 0 | <input type="checkbox"/> L<br><input type="checkbox"/> R | N/A                        | N/A                        | <input type="checkbox"/> S | <input type="checkbox"/> | <input type="checkbox"/>   | <input type="checkbox"/> |
| Virto M-312      | <input type="checkbox"/> L<br><input type="checkbox"/> R | N/A                        | N/A                        | N/A                        | N/A                      | <input type="checkbox"/> S | <input type="checkbox"/> |

**Step 3: Performance Level**

Virto M-Titanium

M90  M70

Virto M Acryl

M90  M70  M50  M30

**Step 4: Power Level**

Left:  M S P SP UP<sup>1</sup>

Right:  M S P SP UP<sup>1</sup>

**Step 5: Colours**

**Shell**

- Black (06)<sup>2</sup>
- Transparent (21)
- Pink (26)
- Red/Blue (R/B) Standard for IIC
- Titanium gray (U0)<sup>5</sup>

**Faceplate**

- Black (06) N/A for Virto M-312 NW 0
- Cocoa (22) Standard for IIC
- Pink (26)
- Brown (28)

**Step 6: Options**

[ ] / (S) = Standard

[ ] = Optional

For selection options and codes, see combination matrix on the other side.

Fast order ►► Build product as standard version

|                           |  |                             |
|---------------------------|--|-----------------------------|
| Vent size                 | AOV (enter audiogram)  | <input type="checkbox"/> S  |
|                           | AOV-0 <sup>5</sup>   | <input type="checkbox"/>    |
|                           | Manually selected  | Left: _____<br>Right: _____ |
| User controls             | Standard controls  | <input type="checkbox"/> S  |
|                           | Push Button <sup>4</sup>   | <input type="checkbox"/>    |
|                           | Volume Control (VC) N/A for Titanium                               | <input type="checkbox"/>    |
|                           | T-coil (TP) Public T-coil program only <sup>3,8</sup>              | <input type="checkbox"/>    |
| Wax system                | EasyPhone (available for acoustic phone program only) <sup>3</sup> | <input type="checkbox"/>    |
|                           | Cerustop (CS)  | <input type="checkbox"/> S  |
|                           | Manually selected  | Type: _____                 |
| Removal line (RF)         | Transparent (21) <sup>6</sup>                                      | <input type="checkbox"/>    |
| Titanium surface          | Extra retention (ERS)/Matte  | <input type="checkbox"/>    |
| Acrylic/retention options | See matrix for options   |                             |
| Accessories               | TV Connector <sup>2</sup>  | <input type="checkbox"/>    |
|                           | Phonak PartnerMic <sup>TM2</sup>                                   | <input type="checkbox"/>    |
|                           | Phonak RemoteControl <sup>2</sup>                                  | <input type="checkbox"/>    |
|                           | MiniControl <sup>5</sup>   | <input type="checkbox"/>    |

Special Instructions  Call me if changes are needed

027-1190/Phonak Virto M Order Form/112/feb2023/V4

<sup>1</sup> Virto M-312 & Virto M-312 NW 0 only

<sup>2</sup> Virto M-312 only

<sup>3</sup> N/A for Virto M-312

<sup>4</sup> Standard for Virto M-312

<sup>5</sup> Virto M-Titanium & Virto M-10 NW 0 only

<sup>6</sup> Standard for IIC & CIC

<sup>7</sup> Colour will match shell colour

<sup>8</sup> Not available in combination with VC

## Combination Matrix

|                           |          | Virto M-Titanium<br> | Virto M-10 NW O<br> | Virto M-312 NW O<br> | Virto M-312<br>  |
|---------------------------|----------|---|--|---|---|
| <b>SHELL STYLE</b>        |          | <b>IIC, CIC</b>   | <b>IIC, CIC, MC, ITC, HS, FS</b>   | <b>MC, ITC, HS, FS</b>  | <b>HS, FS</b>   |
| Push button               | Standard | No push button (PX)   | Push button (PM1) <sup>8</sup>   | No push button (PX)   | Push button (PM1) – required  |
|                           | Optional | Push button (PM1) <sup>10</sup>   | No push button (PX)  | Push button (PM1) <sup>9</sup>  | -   |
| Volume Control            | Standard | -   | No VC (VCX)  | No VC (VCX)   | No VC (VCX)   |
|                           | Optional | -   | Digital VC (VC3)   | Digital VC (VC3)  | Digital VC (VC3)  |
| T-Coil                    | Standard | No T-coil (TX)  | No T-coil (TX)   | No T-coil (TX)  | -   |
|                           | Optional | T-coil (TP)   | T-coil (TP) <sup>8</sup>   | T-coil (TP) <sup>8</sup>  | -   |
| Magnetic switch function  | Standard | MiniControl (DT5) <sup>9,10</sup>   | EasyPhone (DT4) <sup>8</sup>   | EasyPhone (DT4) <sup>8</sup>  | -   |
|                           | Optional | EasyPhone (DT6) <sup>9</sup>  | No EasyPhone (RX)<br>MiniControl (DT5) <sup>9</sup>  | No EasyPhone (RX)   | -   |
| Wax system                | Standard | Cerustop (CS) – required  | Cerustop (CS)  | Cerustop (CS)   | Cerustop (CS)   |
|                           | Optional | -   | Extended rec. tube (ER)<br>Wax Spring (WS)<br>Ext. Tube + Spring (WE)<br>Wax guard not required (X)  | Extended rec. tube (ER)<br>Wax Spring (WS)<br>Ext. Tube + Spring (WE)<br>Wax guard not required (X)     | Extended rec. tube (ER)<br>Wax Spring (WS)<br>Ext. Tube + Spring (WE)<br>Wax guard not required (X) |
| Removal line              | Standard | Removal line transparent (RF) <sup>6</sup>  | Removal line transparent (RF) <sup>6</sup>   | -   | -   |
|                           | Optional | -   | Removal line transparent (RF)  | Removal line transparent (RF)   | Removal line transparent (RF)   |
| Acrylic/retention options | Optional | -   | Canal lock <sup>7</sup><br>Helix lock <sup>7</sup><br>Raised VC<br>Canal bell                        | Canal lock <sup>7</sup><br>Helix lock <sup>7</sup><br>Raised VC<br>Canal bell                           | Canal lock <sup>7</sup><br>Helix lock <sup>7</sup><br>Raised VC<br>Canal bell                       |

<sup>6</sup> Standard for IIC and CIC

<sup>7</sup> Colour will match shell colour

<sup>8</sup> Not available in combination with VC

<sup>9</sup> Choose between MiniControl or EasyPhone

<sup>10</sup> Choose between Push button or MiniControl

## Titanium Information

### Vent Options

**AOV** = Acoustically Optimized Vent. Phonak will set the correct vent diameter and style in relation to Audiogram/shell geometry/feedback canceler and required low frequency gain. This is the recommended vent type.

**AOV-O** = Acoustically Optimized Vent-Open. This specific vent is more open and recommended for first time users and clients who are sensitive to occlusion. Only available for Virto M-Titanium.

### Titanium FitGuide


Measure the possible insertion depth. The device will be built in accordance with the indicated FitGuide value. If the measurement falls between two numbers, please use the smaller number.

- Side M – for devices with M-receiver
- Side P – for devices with P-receiver
- Models with SP receiver will be built without FitGuide information

### Impression Taking

Recommendations for taking impressions:

- EasyView Otoblock placed past the second bend **1**
- Low viscosity material **2**
- Always fill the complete concha with impression material **3**
- Intertragal notch must be visible **4**
- When fitting Titanium, impression length must cover the evaluated FitGuide value **5**



Check out our video on how to take ear impressions!



All of our products, including custom products and spare parts, can be ordered online at the Phonak store.