

Phonak EON cShell 6.0 order form.

Customer information

Ship to account number:

Address:

City: State: Zip:

Bill to account number:

Third party patient number:

Fitting date:

Purchase order number:

Contact name:

Phone:

Email:

Hearing instrument

Qty	Instrument type	Technology level			
		90	70	50	30
	Phonak Audéo™ E-Sphere	☐	☐	☐	-
	Phonak Audéo E-R	☐	☐	☐	☐
	Phonak CROS E-R	-	-	-	-

Specify color choice:

See page 2 for color options

cShell style



cShell Acrylic (062-0056-01)

☐ L ☐ R

Shell colors: Black (06), Cocoa (22), Pink (26), Tan (14), Brown (28), Red (R), Blue (B), Transparent (21)

Faceplate colors: Black (06), Cocoa (22), Pink (26), Tan (14), Brown (28)



cShell Titanium (062-0057-01)

☐ L ☐ R

Shell color: Titanium gray (U0)

Faceplate color: Black (06)

Accessories and Roger™

Accessories	TV Connector (076-3006-0611)	☐
	PartnerMic (076-4001-P511)	☐
	RemoteControl (076-0065-P5)	☐
Roger	Roger On™ 3 (056-3032-xx011) Champagne (P5), Graphite Gray (P7)	☒ P5 ☐ P7

Preferences

If necessary, may we change the following: ☐ Please call

Yes No

☐

☐

UP to P OR P to M OR M to S if audiogram permits

☐

☐

Build larger if required

☐

☐

I would like to receive a paper user guide

Patient information

Last name:

First name:

Age:

Gender:

Is this patient order being fully or partially paid for by state or federal government funds? ☐ Yes ☐ No

Audiogram (Required for AOV):

	Hz	250	500	1K	2K	3K	4K
Left:	AC						
Right:	AC						

HI warranty

☐ 4th year ☐ 5th year (Patient age required at time of order)

Shipping and handling options[†]

☐ 2 Business Days (\$20.99) ☐ 3–5 Business Days (\$19.99)

[†] Based on location and courier

Options

Shell color	Standard color (Transparent 21)	☒ S
	Custom color	Code:
Faceplate color	Standard color (Transparent 21)	☒ S
	Custom color	Code:
Receiver	S Receiver	L ☐ R ☐
	M Receiver (standard)	L ☒ S R ☒ S
	P Receiver	L ☐ R ☐
	UP Receiver (may result in cShell XL)	L ☐ R ☐
Length	0/1/2/3 Size 2 is standard	Left: Right:
Vent size	AOV (enter audiogram)	☒ S
	Customer specific 12/20/25/30/3L/none	Left: mm Right: mm
Wax system	Standard wax system**	☒ S
	Custom wax system	
Retention*	Helix Lock (HL)	☐
	Canal Lock (CL)	☐
	Skeleton Lock (SL)	☐
	Removal line (always included)	☒
Removal line (RF)	Line extension + customer specific	
Titanium surface	Extra retention (ERS)	☐

Special instructions

Additional information

Alternative cShell build style

We offer an alternative cShell build style with an anterior tube outlet. This build style provides high cosmetic acceptance in spite of a slightly larger size and can be used in case of build problems due to anatomical limitations. This build style is available for acrylic shells only.



Fig.1 cShell - standard tube



Fig.2 cShell - anterior tube outlet

Acoustically Optimized Vent (AOV)

We recommend to select AOV. We will set the correct vent diameter and style consider audiogram/shell geometry/feedback canceler and required low frequency gain. The client's audiogram is mandatory for the calculation of the AOV.

Extra Retention Surface (ERS)

We offer a rougher surface for cShell titanium. The extra retention surface can improve retention in the ear.

Wax systems


cShell Acrylic	Description
X	No wax guard
CB	CeruStop
H4	HF4 Filter, default with UP receivers
ER	Extended receiver tube

Color codes

Code	Color
06	Black
13/21	Transparent
14	Tan
22	Cocoa
26	Pink
28	Brown
R/B	Red/Blue
U0	Titanium Gray

Receiver tube length

Measure the tube length to select the correct receiver. Measure both ears for binaural fitting.



EON hearing instrument color options



8A
Sienna Red



P1
Sand Beige



P3
Sandalwood



P4
Chestnut



P5
Champagne



P6
Silver Gray



P7
Graphite Gray



P8
Velvet Black