

Phonak RIC Coupling Guide.

EON



Audéo Sphere Audéo R


Infinio




Audéo R

Receivers (SDS 6.0)


Lengths: 00, 0, 1, 2, 3




Standard (S)
47/112 (2cc)



Medium (M)
51/114 (2cc)



Power (P)
59/122 (2cc)



UltraPower (UP)
71/132 (2cc)
cShell only



CROS E-R



CROS I-R



CROS Tube 6.0
(sizes 0-3)



CROS Tip 6.0


Lumity




Audéo 312 Audéo R/RT

Receivers (SDS 5)


Lengths: 0, 1, 2, 3




Standard (S)
47/111 (2cc)



Medium (M)
51/114 (2cc)



Power (P)
59/122 (2cc)



UltraPower (UP)
71/132 (2cc)
cShell only



CROS R



CROS Tube
(sizes 0-3)



CROS Tip

EasyGuard™ Domes

Now compatible with SDS 5 and 6.0



Open
S, M, L



Vented
S, M, L



Power
S, S/M, M, M/L, L

Domes



Cap Dome



Open Dome



Vented Dome



Power Dome

Custom earpieces



SlimTip Acrylic



SlimTip Silicone




cShell Acrylic




cShell Titanium

ActiveVent (MAV)

4.0 (Lumity)
50/119 (2cc)



6.0 (Infinio)
50/119 (2cc)



Universal SlimTip AV
Sizes A, B, C and D



SlimTip
Titanium AV

Custom Product Guide.

Infinio



Shell colors



Faceplate colors




*Virto I-R and Virto I-10 NW O only


Paradise



Titanium shell color




Fitting ranges



Mild to moderate (M) 40/109
Mild to moderately-severe (P) 50/115
Mild to severe (SP) 60/119
Moderate to severe/profound (UP) 70/127 R and 312 models only

Wax protection



CeruStop wax filter


SlimTip custom earpiece (RIC/BTE)



cShell custom earpiece (RIC)



Shell colors



Faceplate colors



Wax protection options

