



There is a Phonak for every ear.

Receiver-in-the-canal



Audéo EON



Sphere



R

Audéo Infinio



R

Audéo Lumity



312



RT

CROS



EON R



Infinio R



Lumity R

Behind-the-ear



Naída Lumity



M
312



PR



SP
13



UP
675

Sky Lumity



M
312



PR



SP
13



UP
675

Bimodal solutions



Naída
Link M
13



Sky
Link M
13

In-the-ear



Virto Infinio



R



Titanium
10



10 NW 0
10

Virto Paradise



312

Roger™



Roger On™



Roger Select™



Roger Table Mic

Wireless accessories



PartnerMic



RemoteControl



TV Connector

Apps



myPhonak



myPhonak Junior



myRogerMic

Hearing aids and accessories are not shown to scale or actual size

Phonak Charger Compatibility Guide.

	ChargerGo RIC EON	Charger RIC EON	ChargerGo RIC Sphere Infinio	ChargerGo RIC Infinio	Charger RIC Infinio	ChargerGo ITE	Charger Ease	Charger Combi BTE 2*
								
	075-0108-11	075-0073-11	075-0057-11	075-0055-11	075-0054-11	075-0101-11	075-0080-11	075-0065-11
Audéo E-Sphere	•	•						
Audéo E-R	•	•						
CROS E-R	•	•						
Audéo I-Sphere			•					
Audéo I-R				•	•			
CROS I-R				•	•			
Virto I-R						•		
Audéo L-R/RT							•	
CROS L-R							•	
Naída L-PR								•
Sky L-PR								•
Roger Focus II								•

* Same compatibility applies to Phonak Charger Combi BTE (075-0062-xx)