

PHONAK



Audéo EON R and CROS EON R

Clarity you hear.
Connection you feel.



EON

Clarity you hear. Connection you feel.

Phonak Audéo™ EON R hearing aids are designed to help you hear what matters. Powered by Phonak AI, EON automatically adapts to changing sound environments to help deliver a more natural listening experience and keep you connected to the people and moments in your life.

For people with single-sided deafness, CROS EON R pairs with Phonak Audéo EON and transmits sound from the non-hearing side to the better hearing side, helping make sounds from that side audible through the better-hearing ear.




Clarity that speaks for itself

Phonak hearing aids are powered by AI and made for living, adapting seamlessly to your surroundings to keep you connected all day.

Phonak EON is hearing technology that works for you, so you can focus on life.

- **Clarity that moves with you automatically.¹**
- **Easier conversations start now.** Enjoy being in the moment.
- **Our most natural sound available.²** Outdoors or indoors, EON automatically adjusts to different sound environments thanks to Phonak AI.



¹ Preuss, M., Ulloa Sanchez, D. (2025). Phonak Audéo R Infinio outperforms competitors for speech in loud noise performance. Phonak Insight retrieved from www.phonak.com/evidence

² Vaisberg, J., Brunner, J., & Ulloa-Sanchez, D. (2026). Hearing confidently in wind: WindBlock 2.0 delivers comfort, natural sound quality and clearer speech. Phonak Field Study News retrieved from www.phonak.com/evidence



Unlimited connections. Seamless listening.

Connect to the sounds that matter

Phonak EON Made For All connectivity lets you stream music, take calls, and connect to Auracast™ broadcast audio directly through your hearing aids. Easily connect your smartphone or other Bluetooth® enabled devices with a few taps and enjoy seamless listening.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc., and any use of such marks by Sonova AG is under license.

The Auracast™ word mark and logos are trademarks owned by Bluetooth SIG, Inc., and any use of such marks by Sonova AG is under license.

Experience the front row with Auracast™ broadcast audio

Auracast™ broadcast audio is found in places like theaters, museums, lecture halls and gyms. It lets you stream performances, tour guidance and announcements directly to your hearing aids.

Joining a broadcast is easy

Scan a QR code with the myPhonak app or use an Auracast™ assistant on a compatible smartphone to stream audio directly to your hearing aids.





Designed to remove limits

Crafted to deliver all-day comfort and durability, ready for wherever life takes you.

- **Designed to fit your life:** with a listening experience tailored just for you
- **Trusted by experts, backed by science, chosen by millions worldwide:** enjoy peace of mind for years to come
- **Sleek, lightweight design:** with durable materials and a comfortable fit
- **Automatic as you move:** adapts to your listening environment as you transition throughout your day
- **Stay powered on the go:** with compact, easy-to-use charging

From effortless connectivity to all-day comfort, Phonak Audéo E-R and CROS E-R are designed to keep pace with your life.



Meet the new Phonak EON

Audéo EON R

The Audéo R is a rechargeable receiver-in-canal solution for people with mild to profound hearing loss. It offers clarity that moves with you automatically,¹ and Made for All connectivity, including Auracast™ broadcast audio.

CROS EON R

Specifically designed for people with unaidable hearing loss on one side, CROS technology transmits sound wirelessly from the unaidable ear to the better-hearing ear, helping you stay engaged in the conversation.

Compatibility note: For those experiencing unaidable hearing loss on one side, Audéo EON R pairs seamlessly with CROS EON R as part of an integrated solution.

Battery

Power that keeps up with you. Reliable all-day battery life means you can focus on living without interruption.

- Audéo EON R: Up to 38-hour battery life*
- CROS EON R: Up to 20-hour battery life*

* The maximum battery runtime is determined by battery condition, hearing loss, environmental sound classification, streaming usage and activated feature set.

Colors

- | | |
|---|---|
|  Sienna Red  |  Champagne |
|  Sand Beige |  Silver Gray |
|  Sandalwood |  Graphite Gray |
|  Chestnut |  Velvet Black |

Chargers



ChargerGo RIC EON

Includes an integrated battery for charging on the go



Charger RIC E

Compact charger with magnetic lid and device retention



Maximize hearing potential with Roger

Experience clearer conversations even in noisy environments with Roger On™ — your personal third ear. This wireless microphone connects directly to Phonak hearing aids, transmitting voices during one-to-one chats, group conversations, hearing from a distance, or from TV or video calls for easier listening.

Roger On supports better speech understanding³ and helps you stay connected in challenging situations, so you stay in life's special moments.



Discover how Roger On can empower your hearing journey at phonak.com/Roger.

³ Thibodeau L. M. (2020). Benefits in Speech Recognition in Noise with Remote Wireless Microphones in Group Settings. *Journal of the American Academy of Audiology*, 31(6), 404–411.

Managing your hearing just got easier

Take control of your hearing experience right from your smartphone with the myPhonak app.

- Adjust volume and save personalized programs
- Easily connect to Auracast™ broadcast audio
- Locate your hearing aids easily with the Find My Hearing Aids feature
- Designed for simplicity and ease of use with an accessible, user-friendly experience for everyone


How to get the app
Scan to download



App Store®



Google Play™



App Store is a service mark of Apple Inc.
Google Play is a trademark of Google LLC.

At Phonak, we see technology as an enabler of life, not the focus of it.

Since 1947, we have advanced hearing technology to help people stay connected to what matters most: the people, moments and experiences that shape our lives.

Today, we bring together Swiss-engineered performance, a seamless connected hearing experience and purposeful design to create hearing solutions built for where life happens.

So, you can stay confident, connected and fully in it.

Life is on.
Phonak is made for it.

phonak.com

028-2916-03/V1.00/2026-08/hc © 2026 Sonova AG. All rights reserved.