

OGA Infinio cShell 6.0 Order Form

(36C79125D0007)



Government Services | Phone 888-561-7555 | Fax 630-836-9770

Customer Information

Payer:

ROES (DALC PH5120018626)

Other

RACHAP (Enclose payment and if paying by cc, RACHAP authorization form)

Ship To:

Address:

City: _____ **State:** _____ **Zip:** _____

PO #:

Contact Name:

Phone Number:

Email Address:

Patient's Last Name:

Patient's 4 Digit SSN:

Age: _____ **Gender:** _____

cShell style (6.0 is only compatible with Infinio RIC devices)



cShell Acrylic [062-0056-01]

L R

Shell colors: Black [06], Cocoa [22], Pink [26], Tan [14], Brown [28], Red [R], Blue [B], Transparent [21]

Faceplate colors: Black [06], Cocoa [22], Pink [26], Tan [14], Brown [28]



cShell Titanium [062-0057-01]

L R

Shell color: Titanium gray [U0]***

Faceplate color: Black [06]

Audiogram (Required for AOV):

	Hz	250	500	1K	2K	3K	4K
Left:	AC						
Right:	AC						

Preferences

If necessary, may we change the following: Please call

Yes No

- UP to P **OR** P to M **OR** M to S if audiogram permits
- Build larger if required
- I would like to receive a paper user guide

Options

Shell Color	Standard color (Transparent 21)	<input checked="" type="checkbox"/> [S]
	Custom color	Code: _____
Faceplate Color	Standard color	<input checked="" type="checkbox"/> [S]
	Custom color	Code: _____
Receiver	S Receiver	L <input type="checkbox"/> R <input type="checkbox"/>
	M Receiver (standard)	L <input checked="" type="checkbox"/> [S] R <input checked="" type="checkbox"/> [S]
	P Receiver	L <input type="checkbox"/> R <input type="checkbox"/>
	UP Receiver (may result in cShell XL)	L <input type="checkbox"/> R <input type="checkbox"/>
Length	0/1/2/3 Size 2 is standard	Left: _____ Right: _____

Options

Vent size	AOV (enter audiogram)	<input checked="" type="checkbox"/> [S]
	Customer specific 12/20/25/3L/none	Left: _____ mm Right: _____ mm
Wax system	Standard wax system**	<input checked="" type="checkbox"/> [S]
	Custom wax system	
Retention*	Helix Lock (HL)	<input type="checkbox"/>
	Canal Lock (CL)	<input type="checkbox"/>
	Skeleton Lock (SL)	<input type="checkbox"/>
Removal line (RF)	Transparent (21)	<input type="checkbox"/>
	Line extension + customer specific	
Titanium surface	Extra retention (ERS)	<input type="checkbox"/>

Special Instructions

Additional Information

Alternative cShell Build Style

We offer an alternative cShell build style with an anterior tube outlet. This build style provides high cosmetic acceptance in spite of a slightly larger size and can be used in case of build problems due to anatomical limitations. This build style is available for acrylic shells only.



Fig.1 cShell - standard tube




Fig.2 cShell - anterior tube outlet

Acoustically Optimized Vent (AOV)

We recommend to select AOV. We will set the correct vent diameter and style consider audiogram/shell geometry/feedback canceler and required low frequency gain. The client's audiogram is mandatory for the calculation of the AOV.

Extra Retention Surface (ERS)

We offer a rougher surface for cShell titanium. The extra retention surface can improve retention in the ear.

Wax systems

cShell Acrylic	Description
X	No wax guard
CB	CeruStop
H4	HF4 Filter, default with UP receivers

Color codes

Code	Color
06	Black
13/21	Transparent
14	Tan
22	Cocoa
26	Pink
28	Brown
R/B	Red/Blue
UO	Titanium Gray***

*** Only available with Titanium material

Receiver tube length

Measure the tube length to select the correct receiver. Measure both ears for binaural fitting.

