

# OGA Infinio/Ultra SlimTip/CROS Tip 6.0 Order Form (36C79125D0007)



Government Services | Phone 888-561-7555 | Fax 630-836-9770

## Customer Information

### Payer:

- ROES (DALC PH5120018626)
- Other
- RACHAP (Enclose payment and if paying by cc, RACHAP authorization form)

Ship To: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_

### PO #:

Contact Name: \_\_\_\_\_

Phone Number: \_\_\_\_\_

Email Address: \_\_\_\_\_

Patient's Last Name: \_\_\_\_\_

Patient's 4 Digit SSN: \_\_\_\_\_

Age: \_\_\_\_\_ Gender: \_\_\_\_\_

## SlimTip style



SlimTip Acrylic [062-0063-01]

L  R

Shell style: Hollow  
 Standard color: Transparent (21)  
 Custom colors: Tan (14), Cocoa (22), Pink (26), Brown (28), Red/Blue (R/B)



SlimTip Silicone [062-0064-01]

L  R

Shell style: Solid  
 Standard color: Transparent (21)



SlimTip Titanium AV [062-0065-01]

Only compatible with ActiveVent Receiver MAV

L  R

Shell style: Hollow Color:  
 Titanium gray (U0)\*\*\*\*



CROS Tip [062-0059-01]

L  R

Shell style: Hollow  
 Standard color: Transparent (21)  
 Custom colors: Tan (14), Cocoa (22), Pink (26), Brown (28), Red/Blue (R/B)

## Audiogram (Required for AOV):

	Hz	250	500	1K	2K	3K	4K
Left:	AC						
Right:	AC						

## Options

Shell Color	Standard color	<input checked="" type="checkbox"/>
	Custom color	Code: _____
Connection	S Receiver	<input type="checkbox"/> L <input type="checkbox"/> R
	M Receiver (standard)	<input checked="" type="checkbox"/> L <input checked="" type="checkbox"/> R
	P Receiver	<input type="checkbox"/> L <input type="checkbox"/> R
Length	Length 00*/0/1/2/3 Size 2 is standard	Left: _____ Right: _____
	AOV (enter audiogram)**	<input checked="" type="checkbox"/>
Vent size	Customer specific**	Left: _____ mm Right: _____ mm
	Cavity (Y)	<input checked="" type="checkbox"/>
Vent style	Regular (P)	<input type="checkbox"/>
	Helix Lock*** (HL)	<input type="checkbox"/>
	Retention**	Canal Lock (CL) <input type="checkbox"/> Skeleton Lock (SL) <input type="checkbox"/>
Removal line	Removal line (always included)	<input checked="" type="checkbox"/>
	Line extension + customer specific length	
Titanium surface	Extra retention (ERS)	<input type="checkbox"/>

## Preferences

If necessary, may we change the following:  Please call

- |                          |                          |  |
|--------------------------|--------------------------|--|
| Yes                      | No                       |  |
| <input type="checkbox"/> | <input type="checkbox"/> | P to M <b>OR</b> M to S if audiogram permits |
| <input type="checkbox"/> | <input type="checkbox"/> | Build larger if required                     |
| <input type="checkbox"/> | <input type="checkbox"/> | I would like to receive a paper user guide   |

## Special instructions


## Additional Information

### Acoustically Optimized Vent (AOV)

We recommend to select AOV. We will set the correct vent diameter and style consider audiogram/shell geometry/feedback canceler and required low frequency gain. The patient's audiogram is mandatory for the calculation of the AOV.

### Vent styles

We offer two vent options for SlimTips.



Regular:

- Classic vent solution
- Fixed diameter (not enlargeable)




Cavity:

- Optimized for best fitrate
- Easy to enlarge by the fitter

### Tube length

Measure the tube length to select the correct SlimTube. Measure both ears for binaural fitting.



### Wax Systems

The following wax systems are used for different types of SlimTip/CROS Tip.

SlimTip Acrylic	CeruStop mounted in SlimTip <sup>†</sup>
SlimTip Silicone	CeruStop in receiver
SlimTip Titanium AV	HF3 mounted in ActiveVent receiver
CROS Tip	No wax filter needed

<sup>†</sup> CeruStop from receiver must be removed

### SlimTip Titanium AV

This variant of the SlimTip is designed for ActiveVent™ (AV) and is compatible to the MAV receiver (MAV = Medium ActiveVent) only. There is no additional vent option needed. It is offered in durable Titanium to maximize fitrate and increase comfort.

Our SlimTip Link AV connector makes the assembly of the MAV Receiver simple and easy.

ActiveVent is compatible to Infinio rechargeable RIC hearing aids.

### Extra Retention Surface (ERS)

We offer a rougher surface for Titanium. The extra retention surface can improve retention in the ear.

### Hearing instrument compatibility

The SlimTip 6.0 portfolio is compatible to Infinio hearing aids or newer. It is not backward compatible to Lumity or earlier generation hearing aids.

### Color codes

Code	Color
14	Tan
21	Transparent
22	Cocoa
26	Pink
28	Brown
R/B	Red/Blue
U0	Titanium Gray****

\*\*\*\*Only available with Titanium material

All of our products, including custom products and spare parts, can be ordered online in the Phonak store.

Sonova USA Inc. is not responsible and assumes no liability for any device or accessory not manufactured by Sonova which is sent by you to Sonova. Please ensure that you only include Phonak branded devices and accessories herein.

Products, options and accessories are subject to change without notice.