

Phonak

Charge and Care

Technical Data

Dimensions (L x W x H)	102 x 78 x 52 mm (4 x 3 x 2")
Weight	85 g (3 oz)
Color	Champaign
Operating conditions	5° to +40° Celsius (41° to 104° Fahrenheit)
Transport and storage conditions	-20° to +60° Celsius (-4° to +140° Fahrenheit)
Humidity during transportation	Up to 90% (non condensing)
Humidity during storage	0% to 70% (non condensing)
Humidity during operation	0% to 85% (non condensing)
Atmospheric pressure during operation	800 hPa to 1060 hPa
Atmospheric pressure during storage and transportation	500 hPa to 1060 hPa
Included Parts	Power supply micro-USB cable (for charging only)
External Interfaces	micro-USB port
Battery	None
Input	Type: DC, stabilized 5 V +/- 10% Max. current: 500 mA

