

Phonak ChargerGo ITE

Technical data

Dimensions (L x W x H)	99 x 57 x 54 mm (3.9 x 2.2 x 2.1") +/- 5%
Weight	84 g (3 oz) +/- 10%
Color	Graphite gray
Operating temperature	+5 °C to +40 °C (+41 °F to +104 °F)
Transport and storage temperature	-20 °C to +60 °C (-4 °F to +140 °F)
Operating humidity	0% to 90% (non-condensing)
Transport and storage humidity	0% to 93% (non-condensing)
Atmospheric pressure during use, storage and transport	500 hPa to 1060 hPa
Included parts	Power supply USB-C cable
External interfaces	USB-C port
Battery	Rechargeable lithium-ion
Charging time	See user guide
Input	5 V DC +/-5%, current 500 mA
Rechargeable lithium-ion battery capacity	1.184 Wh



The battery cannot be overcharged. It may become warm during charging.

V2.00/2026-01/ds © 2026 Sonova AG All rights reserved