

Phonak Charger Virto 1

Technical data.

Dimensions (L x W x H)	99 x 59 x 49 mm (3.9 x 2.3 x 1.9") +/- 5%
Weight	68 g (2.4 oz) +/- 20%
Color	Champagne
Operating conditions	+5° to +40° Celsius (+41° to +104° Fahrenheit)
Transport and storage conditions	-20° to +60° Celsius (-4° to +140° Fahrenheit)
Humidity during transportation and storage	0% to 93% (non condensing)
Humidity during operation	0% to 90% (non condensing)
Atmospheric pressure during use	500 hPa to 1060 hPa
Atmospheric pressure during storage and transportation	500 hPa to 1060 hPa
Included Parts	Power supply USB-C cable (for charging only)
External Interfaces	USB-C port
Input	Type: DC, stabilized 5 V +/-10% Max. current: 600 mA



V1.00/2024-03/na © 2024 Sonova AG All rights reserved