

Phonak Infinio Ultra.

The very best. For you.





Phonak AI hearing aids are designed to help people hear better in the real world. Infinio Ultra delivers exceptional sound and speech understanding, no matter the situation. Powered by adaptive AI precision to help your patients hear the world more clearly and confidently. Every step of the way.



Adaptive

The best sound for your patients. No matter the situation. Anytime. Anywhere.

- Adaptive Phonak Digital 3.0 delivers instant comfort & instant love with 93% patient satisfaction from the first fit1
- AutoSense OS™ 7.0 harnesses the power of AI to accurately identify and adjust now trained with 18 times more sound environments than before, boosting speech understanding in noise²
- SmartSpeech Technology delivers up to 8 dB SNR benefit in group conversations, outperforming four other brands, automatically with AutoSense OS 7.0³

Easier maintenance for everyone. Phonak does the work, so you don't have to.



EasyGuard™ revolutionizes wax management with its innovative, sound-transmitting, sealed membrane

- Designed to protect the receiver from wax buildup, helping to extend the lifespan of hearing aids
- Designed to reduce service appointments by up to 38%⁴







- Stewart, E., Adler, M., & Seitz-Paquette, K., (2024). Adaptive Phonak Digital (APD) 3.0 is the preferred first fit compared to a key competitor device. Phonak Field Study News retrieved from www.phonak.com/evidence
- Wright, A. et al. (2025). AutoSense OS 7.0 improves speech understanding with highly rated sound quality for challenging listening environments. Phonak Field Study News retrieved from www.phonak.com/evidence
- Preuss, M., & Ulloa Sanchez, D. (2025). Phonak Audéo R Infinio outperforms competitors for speech in loud noise performance. Phonak Insight retrieved from www.phonak.com/evidence
- ⁴ Based on Phonak internal data. Contact claims@sonova.com if you are interested in more information



With Infinio Ultra complete the feedback and real-ear test in under 13 seconds: now faster and available via Remote Support!




All the benefits of the latest technology for your patients with 40% faster binaural firmware updates



Al Dome Proposer leverages Al to recommend the optimal dome, helping you and your patients achieve the perfect fit from the very first moment

Connected



Pairs with up to 8 Bluetooth® devices, with the flexibility to have 2 simultaneously connected





Simpler, faster, better

Ultra-responsive ERA chip provides uninterrupted, hands-free calls and



streaming across up to 2x the distance, thanks to up to 6x more transmission power



One-step pairing to the myPhonak app, Roger™ microphones, smartphones and tablets for seamless connection to all your sound sources

Tested. Reliable. Proven. Recommend it with confidence.


- Our most tested hearing aid with over 10,000 hours and 135 individual stringent stress tests*
- Electric components are protected with a durable coating ensuring a prolonged lifespan and consistent performance

^{*} Sonova applies 135 individual stringent tests to all our products

Spheric Speech Clarity like never before

Any voice. Any direction. Always clear.

Spheric Speech Clarity 2.0 extracts, enhances and integrates speech from unlimited voices in real-time to deliver unprecedented speech clarity in noise from any direction.⁵ DEEPSONIC™, the world's first dedicated AI chip for all-around, speech from noise separation, is now 30% more efficient delivering a new level of clear speech⁴ throughout the day.



Easier listening, more energy for what matters**

^{*} When compared to two key competitors against StereoZoom 2.0; adults with moderate to severe hearing loss

^{**}Relative performance change compared to StereoZoom 2.0, subjects with moderate-severe hearing loss



- ⁵ Raufer, S., Kohlhauer, P., Jehle, F., Kühnel, V., Preuss, M., Hobi, S. (2024). Spheric Speech Clarity proven to outperform three key competitors for clear speech in noise. Phonak Field Study News retrieved from www.phonak.com/evidence
- 6 Wright, A., Kuehnel, V., Keller, M., Seitz-Paquette, K., Latzel, M. (2024). Spheric Speech Clarity applies DNN signal processing to significantly improve speech understanding from any direction and reduce the listening effort. Phonak Field Study News retrieved from www.phonak.com/evidence
- ⁷ Tian, E. et al. (2025) Spheric Speech Clarity significantly improves speech understanding in a multi-talker scenario and reduces spatial listening effort with Mandarin speakers. Phonak Field Study News. Retrieved from www.phonak.com/evidence
- ⁸ Latzel, M et al. (2025) Spheric Speech Clarity improves speech understanding and reduces both listening effort and fatigue. Phonak Field Study News. Retrieved from www.phonak.com/evidence

The Infinio Ultra portfolio

Phonak Audéo™ R Infinio



Performance levels: 90, 70, 50, 30



Highlights:

- AutoSense OS 7.0 adapts 24% more precisely to listening situations* so that wearers can enjoy exceptional sound quality in every moment
- Let your patients immerse themselves in their favorite tunes with AutoSense OS 7.0, our most precise sound classification system, recognizing music styles with exceptional accuracy
- Offer your patients the freedom of instant compatibility with all Bluetooth® devices
- Compatible with Phonak patented EasyGuard dome
- Up to 31-hour battery life**, ensuring all-day peace of mind

Colors:



Compatible chargers:



- * Compared to previous generations
- ** The maximum battery runtime is determined by battery condition, hearing loss, environmental sound classification, streaming usage and activated feature set

Phonak Audéo Sphere™ Infinio



Performance levels: 90, 70





Highlights:

- Featuring all the benefits of Audéo R Infinio
- Spheric Speech Clarity 2.0, works in real time to catch conversations from multiple speakers and any angle, simultaneously, even in the noisiest environments
- Sphere integrates the ultra-responsive ERA chip with the DEEPSONIC chip, to effectively separate speech from noise
- Offer your patients the freedom of instant compatibility with all Bluetooth® devices
- Compatible with Phonak patented EasyGuard dome
- Up to 56 hours of battery life*

Colors:



Compatible charger:



The maximum battery runtime is determined by battery condition, hearing loss, environmental sound classification, streaming usage and activated feature set

Key feature comparison



Experience the clarity of Phonak Audéo Sphere™ Infinio

90 Sphere

Spheric Speech Clarity 2.0 (full strength)

70 Sphere

Spheric Speech Clarity 2.0 (medium strength) SpeechSensor







Phonak Audéo I-R

90 Premium

Speech in loud noise Music

Speech in car
Comfort in noise
Speech in noise
Calm situation
Streaming speech

Streaming music SpeechSensor Speech Enhancer

StereoZoom 2.0 Dynamic Noise Cancellation Motion Sensor Hearing ActiveVent Receiver*

Roger Technology Speech in loud noise (manual) Speech in 360° (manual) Tap Control SoundRelax

DuoPhone Real Ear Sound UltraZoom WindBlock SoundRecover2 WhistleBlock Tinnitus Balance 20 channels

70 Advanced Music Speech in car Comfort in noise Speech in noise Calm situation Streaming speech Streaming music

StereoZoom 2.0 Dynamic Noise Cancellation Motion Sensor Hearing ActiveVent Receiver*

Roger Technology Speech in loud noise (manual) Speech in 360° (manual) Tap Control

Real Ear Sound UltraZoom WindBlock SoundRecover2 WhistleBlock

SoundRelax

DuoPhone

* Requires ActiveVent Receiver with Infinio RIC hearing aids

50 Standard

16 channels

Speech in car	
Comfort in noise	
Speech in noise	
Calm situation	
Streaming speech	30 E
Streaming music	Spee
Motion Sensor Hearing	Calm
ActiveVent Receiver*	Strea
Roger Technology	Strea
SoundRelax	Activ
DuoPhone	Roge
Real Ear Sound	Real I
UltraZoom	Ultra
WindBlock	Wind
SoundRecover2	Soun
WhistleBlock	Whist
Tinnitus Balance	Tinnit

ssential

30 Essential
Speech in noise
Calm situation
Streaming speech
Streaming music
ActiveVent Receiver*
Roger Technology
Real Ear Sound
UltraZoom
WindBlock
SoundRecover2
WhistleBlock
Tinnitus Balance
12 channels

Tinnitus Balance

20 channels

Refer to the product information document for details on feature availability at different performance levels

Maximize hearing potential with Roger™ microphones



Roger is an intelligent, wireless technology that transmits speech directly to hearing aids via RogerDirect™. Functioning like a third ear, Roger gives an extra boost in noisy environments and over distance.



Up to 62% more understanding of speech in noise and over distance with Roger, compared to people without hearing loss⁹



Reduce the effects of delay in remote microphones with Roger, having the fastest latency in the industry*

Roger microphones



Roger On™



Roger Select™



Roger Table Mic


- When tested against three broadcast audio competitors
- ⁹ Thibodeau, L. (2014). Comparison of speech recognition with adaptive digital and FM wireless technology by listeners who use hearing aids. American Journal of Audiology, 23(2), 201-210



myPhonak 7.2 & myPhonak Junior 2.0

Simplified experience and new features

Phonak innovation extends to more than just devices, with a host of improvements being made for myPhonak and myPhonak Junior to ensure a more user-friendly experience and greater flexibility.





Target 11.1 updates to benefit you and your patients

Phonak Target fitting software is designed to support your fitting success. Offer your patients a positive and efficient fitting experience.

- Feedback test completed in under 13 seconds
- 40% faster simultaneous firmware updates
- Simultaneous Trial performance upgrades
- Hands-free call disablement for compliance with high security facilities, such as SCIF (Sensitive Compartmented Information Facility)
- Empowering you to personalize features in the workflow



Practice Forward: Empowering your practice



Your journey is one of continuous learning and adaptation in the dynamic world of hearing care. 93% of patients reported higher trust in their audiologist when the professional demonstrated up-to-date knowledge and expertise. ¹⁰ Crafted just for you, our advanced technologies and best practices amplify the impact you can make in your community. Let's navigate the future together.

Contact your Phonak representative to learn more about Practice Forward.



life is on

Since 1947, Phonak is dedicated to preserving life quality by opening new acoustic worlds. Already back then, in the days of our foundation, our company was driven by a formative conviction: We believe that well-hearing equates to well-being and thus, is essential for living life to the fullest. In fact, the sense of hearing is directly linked to social, emotional, cognitive and physical well-being. Today as in future, we thrive to offer the broadest portfolio of innovative hearing solutions. And, together with our hearing care professionals, we keep on focusing on what matters most: improving speech understanding, changing people's lives and having a positive effect on society as a whole.

phonak.com



