Phonak cShell 6.0 order form.



Customer information Patient information Ship to account number: Last name: Address: First name: City: State: Zip: Age: Gender: Bill to account number: Is this patient order being fully or partially paid for by state or federal government funds? Yes No Third party patient number: Audiogram (Required for AOV): Fitting date: Purchase order number: Hz 250 500 1K 2K 3K 4K Contact name: AC Left: Riaht: AC. Phone: Email: HI warranty Hearing instrument 4th year 5th year (Patient age required at time of order) Qty Instrument type Technology level Shipping and handling options[†] 50 70 90 30 2 Business Days (\$20.99) C.O.D. (\$21.99 + shipping) Phonak Audéo™ I Sphere _ 3-5 Business Days (\$19.99) Phonak Audéo I-R _ [†] Based on location and courier Phonak CROS I-R **Options** Specify color choice:_ Standard color (Transparent 21) (S) See page 2 for color options **Shell Color** Custom color Code: cShell style Standard color (5) **Faceplate Color** cShell Acrylic (062-0056-01) Custom color Code: $\bigcap L \bigcap R$ S Receiver $L \cap$ $R \bigcirc$ Shell colors: Black (06), Cocoa (22), Pink (26), Tan (14), Brown (28), L(S) R(S)M Receiver (standard) Red (R), Blue (B), Transparent (21) Receiver P Receiver R (Faceplate colors: Black (06), Cocoa (22), Pink (26), Tan (14), Brown (28) UP Receiver (may result in cShell XL) $R \bigcirc$ cShell Titanium (062-0057) Left: 0/1/2/3 Length \square L \square R Size 2 is standard Right: Shell color: Titanium gray (U0) AOV (enter audiogram) (S) Faceplate color: Black (06) Vent size Left: mm Customer specific 12/20/25/3L/none Accessories and Roger™ Right: mm Standard wax system** (S)TV Connector (076-3006-0611) Wax system Custom wax system Accessories PartnerMic (076-4001-P511) Helix Lock (HL) RemoteControl (076-0065-P5) Retention* Canal Lock (CL) Roger On™ 3 (056-3032-xx011) Champagne (P5), Graphite Gray (P7) (P5) (P7) Skeleton Lock (SL) (P7) Roger Roger Select[™] 3 (056-3033-xx011) Champagne (P5), Graphite Gray (P7) (P5) Transparent (21) Roger Table Mic (056-4012-P5011) Removal line (RF) Line extension + customer specific **Preferences** Titanium surface Extra retention (ERS) Please call If necessary, may we change the following: **Special Instructions** Yes No UP to P OR P to M OR M to S if audiogram permits Build larger if required I would like to receive a paper user guide

/(S) = Standard

= Optional

Legend

^{*} Not available for cShell Titanium ** HF4 for UP receiver, CeruStop for all other receivers

Additional Information

Alternative cShell Build Style

We offer an alternative cShell build style with an anterior tube outlet. This build style provides high cosmetic acceptance in spite of a slightly larger size and can be used in case of build problems due to anatomical limitations. This build style is available for acrylic shells only.





Fig.1 cShell - standard tube

Fig.2 cShell - anterior tube outlet

Acoustically Optimized Vent (AOV)

We recommend to select AOV. We will set the correct vent diameter and style consider audiogram/shell geometry/feedback canceler and required low frequency gain. The client's audiogram is mandatory for the calculation of the AOV.

Extra Retention Surface (ERS)

We offer a rougher surface for cShell titanium. The extra retention surface can improve retention in the ear.

Wax systems

cShell Acrylic	Description
X	No wax guard
СВ	CeruStop
H4	HF4 Filter, default with UP receivers

Color codes

Code	Color
06	Black
13/21	Transparent
14	Tan
22	Cocoa
26	Pink
28	Brown
R/B	Red/Blue
U0	Titanium Gray

Receiver tube length

Measure the tube length to select the correct receiver. Measure both ears for binaural fitting.



Infinio hearing instrument color options





Champagne



Silver Gray



Graphite Gray