

Phonak SlimTip/CROS Tip 6.0 order form.

PHONAK
life is on

Customer information

Ship to account number:

Address:

City: State: Zip:

Bill to account number:

Third party patient number:

Fitting date:

Purchase order number:

Contact name:

Phone:

Email:

Hearing instrument

Qty	Instrument type	Technology level			
		90	70	50	30
	Phonak Audéo™ I Sphere	<input type="checkbox"/>	<input type="checkbox"/>	-	-
	Phonak Audéo I-R	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Phonak CROS I-R	-	-	-	-

Specify color choice:

See page 2 for color options

SlimTip style



SlimTip Acrylic (062-0063-01)

☐ L ☐ R

Shell style: Hollow

Standard color: Transparent (21)

Custom colors: Tan (14), Cocoa (22), Pink (26), Brown (28), Red/Blue (R/B)



SlimTip Silicone (062-0064-01)

☐ L ☐ R

Shell style: Solid

Standard color: Transparent (21)



SlimTip Titanium AV (062-0065-01)

Only compatible with ActiveVent Receiver MAV

☐ L ☐ R

Shell style: Hollow

Color: Titanium gray (U0)



CROS Tip (062-0059)

☐ L ☐ R

Shell style: Hollow

Standard color: Transparent (21)

Custom colors: Tan (14), Cocoa (22), Pink (26), Brown (28), Red/Blue (R/B)

Preferences

If necessary, may we change the following: ☐ Please call

Yes No

☐
☐
☐

P to M OR M to S if audiogram permits

Build larger if required

I would like to receive a paper user guide

Patient information

Last name:

First name:

Age:

Gender:

Is this patient order being fully or partially paid for by state or federal government funds? ☐ Yes ☐ No

Audiogram (Required for AOV):

	Hz	250	500	1K	2K	3K	4K
Left:	AC						
Right:	AC						

HI warranty

☐ 4th year ☐ 5th year (Patient age required at time of order)

Shipping and handling options[†]

☐ 2 Business Days (\$20.99) ☐ C.O.D. (\$21.99 + shipping)
☐ 3-5 Business Days (\$19.99)

[†] Based on location and courier

Options

Shell Color	Standard color	<input checked="" type="checkbox"/> S
	Custom color	Code:
Connection	S Receiver	<input type="checkbox"/> L <input type="checkbox"/> R
	M Receiver (standard)	<input checked="" type="checkbox"/> S L <input checked="" type="checkbox"/> S R
	P Receiver	<input type="checkbox"/> L <input type="checkbox"/> R
Length	Length 00*/0/1/2/3	Left:
	Size 2 is standard	Right:
Vent size	AOV (enter audiogram)**	<input checked="" type="checkbox"/> S
	Customer specific**	Left: mm Right: mm
Vent style	Cavity (Y)	<input checked="" type="checkbox"/> S
	Regular (P)	<input type="checkbox"/>
Retention**	Helix Lock*** (HL)	<input type="checkbox"/>
	Canal Lock (CL)	<input type="checkbox"/>
	Skeleton Lock (SL)	<input type="checkbox"/>
Removal line	Removal line (always included)	<input checked="" type="checkbox"/>
	Line extension + customer specific length	
Titanium surface	Extra retention (ERS)	<input type="checkbox"/>
	TV Connector (076-3006-0611)	<input type="checkbox"/>
Accessories	PartnerMic (076-4001-P511)	<input type="checkbox"/>
	Remote Control (076-0065-P5)	<input type="checkbox"/>
	Roger On™ 3 (056-3032-xx011) Champagne (P5), Graphite Gray (P7)	<input type="checkbox"/> P5 <input type="checkbox"/> P7
Roger	Roger Select™ 3 (056-3033-xx011) Champagne (P5), Graphite Gray (P7)	<input type="checkbox"/> P5 <input type="checkbox"/> P7
	Roger Table Mic 3 (056-4012-P5011)	<input type="checkbox"/>

Special instructions

*Not available for SlimTip Titanium AV and CROS Tip **Not available for SlimTip Titanium AV ***Not available for SlimTip Silicone

Legend ☐ / ☒ S = Standard ☐ = Optional


Additional Information

Acoustically Optimized Vent (AOV)

We recommend to select AOV. We will set the correct vent diameter and style consider audiogram/shell geometry/feedback canceler and required low frequency gain. The patient's audiogram is mandatory for the calculation of the AOV.


Vent styles

We offer two vent options for SlimTips.



Regular:

- Classic vent solution
- Fixed diameter (not enlargeable)




Cavity:

- Optimized for best fitrate
- Easy to enlarge by the fitter

Tube length

Measure the tube length to select the correct SlimTube. Measure both ears for binaural fitting.



Wax Systems

The following wax systems are used for different types of SlimTip/CROS Tip.

SlimTip Acrylic	CeruStop mounted in SlimTip*
SlimTip Silicone	CeruStop in receiver
SlimTip Titanium AV	HF3 mounted in ActiveVent receiver
CROS Tip	No wax filter needed

* CeruStop from receiver must be removed

SlimTip Titanium AV

This variant of the SlimTip is designed for ActiveVent™ (AV) and is compatible to the MAV receiver (MAV = Medium ActiveVent) only. There is no additional vent option needed. It is offered in durable Titanium to maximize fitrate and increase comfort.

Our SlimTip Link AV connector makes the assembly of the MAV Receiver simple and easy.

ActiveVent is compatible to Infinio rechargeable RIC hearing aids.

Extra Retention Surface (ERS)

We offer a rougher surface for Titanium. The extra retention surface can improve retention in the ear.

Hearing instrument compatibility

The SlimTip 6.0 portfolio is compatible to Infinio hearing aids or newer. It is not backward compatible to Lumity or earlier generation hearing aids.

Color codes

Code	Color
14	Tan
21	Transparent
22	Cocoa
26	Pink
28	Brown
R/B	Red/Blue
U0	Titanium Gray

Infinio hearing instrument color options



P1
Sand Beige



P5
Champagne



P6
Silver Gray



P7
Graphite Gray

All of our products, including custom products and spare parts, can be ordered online in the Phonak store.

Sonova USA Inc. is not responsible and assumes no liability for any device or accessory not manufactured by Sonova which is sent by you to Sonova. Please ensure that you only include Phonak branded devices and accessories herein.

Products, options and accessories are subject to change without notice.