

Phonak Virto™ Infinio order form.



Virto 10 NW O and Titanium

Custom hearing instrument ordering available on the **Phonak Store**

Customer information

Ship to account:
 Address: _____
 City: _____
 State: _____ Zip: _____

Bill to account:
 Third party patient number: _____
 Date: _____
 Purchase order number: _____
 Contact name: _____
 Phone number: _____
 Email address: _____

HI warranty options

4th year

Product



Virto I-10 NW O
 L R



Virto I-Titanium
 L R

FitGuide values
 L: _____ R: _____

Performance level

Virto I-10 NW O
 I90 I70 I50 I30

Virto I-Titanium
 I90 I70

Power level

Left: M (40/109) P (50/115) SP (60/119)
 Right: M (40/109) P (50/115) SP (60/119)

Color

Shell

06 Black
 14 Tan
 22 Cocoa
 26 Pink
 28 Brown
 R/B Red/Blue¹
 21 Transparent
 U0 Titanium gray (I-Titanium only)

Faceplate

06 Black²
 14 Tan
 22 Cocoa²
 26 Pink²
 28 Brown

Patient information

Last name: _____
 First name: _____
 Age: _____ Gender: _____

Is this patient order being fully or partially paid for by state or federal government funds? Yes No

Audiogram (Required for AOV/AOV-0):

	HZ	250	500	1K	2K	3K	4K
Left:	AC						
Right:	AC						

Shipping and handling options*

2 Business Days (\$20.99)** 3-5 Business Days (\$19.99)

* Based on location and courier
 ** Not applicable for Titanium

Shell style, options and accessories

Fast order **▶▶** This **(S)** indicates our standard options. We will use these options unless you indicate otherwise.

Shell size	Smallest possible	<input checked="" type="checkbox"/> (S)
	Mini canal	<input type="checkbox"/>
Vent size	AOV (requires audiogram)	<input checked="" type="checkbox"/> (S)
	AOV-0 (Virto I-Titanium)	<input type="checkbox"/>
	Manually selected (12/20/25/3L/SAV/No vent)	Left: _____ mm Right: _____ mm
User controls	Standard controls	<input checked="" type="checkbox"/> (S)
	EasyPhone & push button (DT6, PM1)	<input type="checkbox"/>
	EasyPhone without push button (DT6, PX)	<input type="checkbox"/>
	MiniControl without push button (DT5, PX)	<input type="checkbox"/>
	Volume control (VC)	<input type="checkbox"/>
	T-coil	<input type="checkbox"/>
Wax system	Cerustop (CS) Titanium is only available with CS	<input checked="" type="checkbox"/> (S)
	Manually selected (X-No Waxguard, ER-Extended Tube, WS-Wax Spring, WE-WaxSpring with extended Tube)	Type: _____
Removal line (RF)	Transparent (21)	<input type="checkbox"/>
Titanium surface	Extra retention (ERS)/Matte	<input type="checkbox"/>
Acrylic/retention options	Other	Type: _____
Accessories	Phonak MiniControl	<input type="checkbox"/>
Special Instructions	<input type="checkbox"/> Call me if changes are needed to this order	
	<input type="checkbox"/> Must call to drop receiver size	
	<input type="checkbox"/> I would like to receive a paper user guide	

1 Standard for IIC and CIC shell styles in Virto I-10 NW O
 2 Available for Virto I-Titanium

For selection options and codes see combination matrix on the back side.
 Legend / (S) = Standard = Optional

Combination matrix: Virto 10 NW O and Titanium

		Virto I-10 NW O	Virto I-Titanium
Push button	Standard	Push button (PM1) ⁵	No push button (PX)
	Optional	No push button (PX)	Push button (PM1)
Volume Control	Standard	No VC (VCX)	-
	Optional	Digital VC (VC3) ⁶	-
T-Coil	Standard	No T-coil (TX)	No T-coil (TX)
	Optional	T-coil (TP) ⁵	T-coil (TP)
Magnetic switch function	Standard	EasyPhone (DT6) ⁵	MiniControl (DT5)
	Optional	MiniControl (DT5) No magnetic switch (RX)	EasyPhone (DT6)
Wax system	Standard	Cerustop (CS)	Cerustop (CS)
	Optional	Extended rec. tube (ER) Wax Spring (WS) Ext. Tube + Spring (WE) Wax guard not required (X)	-
Acrylic/retention options	Optional	Canal lock ^{3,4} Helix lock ^{3,4} Raised VC	-

3 Chargeable option

4 Color will match shell color

5 Not available in combination with volume control

6 Not available in combination with push button

Titanium Information

Vent Options

AOV = Acoustically Optimized Vent. Phonak will set the correct vent diameter and style in relation to Audiogram/shell geometry/feedback canceler and required low frequency gain. This is the recommended vent type.

AOV-O = Acoustically Optimized Vent-Open. This specific vent is more open and recommended for first time users and clients who are sensitive to occlusion. Only available for Virto I-Titanium.

Titanium FitGuide

Measure the possible insertion depth. The device will be built in accordance with the indicated FitGuide value. If the measurement falls between two numbers, please use the smaller number.

- Side M - For devices with M-receiver
- Side P - For devices with P-receiver
- Models with SP receiver will be built without FitGuide information

Impression Taking

Recommendation for taking deep impressions:

- Use EasyView Otoblock **1**
- Use Dreve Otoform A softX 25 shore A (#2956) **2**
- Always fill the complete concha with impression material **3**
- Intertragal notch must be visible **4**
- Impression length must cover the evaluated FitGuide value **5**

