

Phonak cShell 6.0 order form.

PHONAK
life is on

Customer information

Ship to account number:

Address:

City: State: Zip:

Bill to account number:

Third party patient number:

Fitting date:

Purchase order number:

Contact name:

Phone:

Email:

Hearing instrument

Qty	Instrument type	Technology level			
		90	70	50	30
	Phonak Audéo™ I Sphere	<input type="checkbox"/>	<input type="checkbox"/>	-	-
	Phonak Audéo I-R	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Phonak CROS I-R	-	-	-	-

Specify color choice:

See page 2 for color options

cShell style



cShell Acrylic (062-0056-01)

☐ L ☐ R

Shell colors: Black (06), Cocoa (22), Pink (26), Tan (14), Brown (28), Red (R), Blue (B), Transparent (21)

Faceplate colors: Black (06), Cocoa (22), Pink (26), Tan (14), Brown (28)



cShell Titanium (062-0057)

☐ L ☐ R

Shell color: Titanium gray (U0)

Faceplate color: Black (06)

Accessories and Roger™

Accessories	TV Connector (076-3006-0611)	<input type="checkbox"/>
	PartnerMic (076-4001-P511)	<input type="checkbox"/>
	RemoteControl (076-0065-P5)	<input type="checkbox"/>
Roger	Roger On™ 3 (056-3032-xx011) Champagne (P5), Graphite Gray (P7)	<input type="checkbox"/> P5 <input type="checkbox"/> P7
	Roger Select™ 3 (056-3033-xx011) Champagne (P5), Graphite Gray (P7)	<input type="checkbox"/> P5 <input type="checkbox"/> P7
	Roger Table Mic (056-4012-P5011)	<input type="checkbox"/>

Preferences

If necessary, may we change the following: ☐ Please call

- | | | |
|--------------------------|--------------------------|--|
| <input type="checkbox"/> | <input type="checkbox"/> | UP to P OR P to M OR M to S if audiogram permits |
| <input type="checkbox"/> | <input type="checkbox"/> | Build larger if required |
| <input type="checkbox"/> | <input type="checkbox"/> | I would like to receive a paper user guide |

* Not available for cShell Titanium

** HF4 for UP receiver, CeruStop for all other receivers

Patient information

Last name:

First name:

Age:

Gender:

Is this patient order being fully or partially paid for by state or federal government funds? ☐ Yes ☐ No

Audiogram (Required for AOV):

	Hz	250	500	1K	2K	3K	4K
Left:	AC						
Right:	AC						

HI warranty

☐ 4th year ☐ 5th year (Patient age required at time of order)

Shipping and handling options[†]

☐ 2 Business Days (\$20.99) ☐ C.O.D. (\$21.99 + shipping)
☐ 3-5 Business Days (\$19.99)

[†] Based on location and courier

Options

Shell Color	Standard color (Transparent 21)	<input checked="" type="checkbox"/> S
	Custom color	Code:
Faceplate Color	Standard color	<input checked="" type="checkbox"/> S
	Custom color	Code:
Receiver	S Receiver	L <input type="checkbox"/> R <input type="checkbox"/>
	M Receiver (standard)	L <input checked="" type="checkbox"/> S R <input checked="" type="checkbox"/> S
	P Receiver	L <input type="checkbox"/> R <input type="checkbox"/>
	UP Receiver (may result in cShell XL)	L <input type="checkbox"/> R <input type="checkbox"/>
Length	0/1/2/3 Size 2 is standard	Left: Right:
	AOV (enter audiogram)	<input checked="" type="checkbox"/> S
Vent size	Customer specific 12/20/25/3L/none	Left: mm
		Right: mm
Wax system	Standard wax system**	<input checked="" type="checkbox"/> S
	Custom wax system	
Retention*	Helix Lock (HL)	<input type="checkbox"/>
	Canal Lock (CL)	<input type="checkbox"/>
	Skeleton Lock (SL)	<input type="checkbox"/>
Removal line (RF)	Transparent (21)	<input type="checkbox"/>
	Line extension + customer specific	
Titanium surface	Extra retention (ERS)	<input type="checkbox"/>

Special Instructions

Legend ☐ / ☒ S = Standard ☐ = Optional

Additional Information

Alternative cShell Build Style

We offer an alternative cShell build style with an anterior tube outlet. This build style provides high cosmetic acceptance in spite of a slightly larger size and can be used in case of build problems due to anatomical limitations. This build style is available for acrylic shells only.



Fig.1 cShell - standard tube




Fig.2 cShell - anterior tube outlet

Acoustically Optimized Vent (AOV)

We recommend to select AOV. We will set the correct vent diameter and style consider audiogram/shell geometry/feedback canceler and required low frequency gain. The client's audiogram is mandatory for the calculation of the AOV.

Extra Retention Surface (ERS)

We offer a rougher surface for cShell titanium. The extra retention surface can improve retention in the ear.

Wax systems


cShell Acrylic	Description
X	No wax guard
CB	CeruStop
H4	HF4 Filter, default with UP receivers

Color codes

Code	Color
06	Black
13/21	Transparent
14	Tan
22	Cocoa
26	Pink
28	Brown
R/B	Red/Blue
U0	Titanium Gray

Receiver tube length

Measure the tube length to select the correct receiver. Measure both ears for binaural fitting.



Infinio hearing instrument color options



P1
Sand Beige



P5
Champagne



P6
Silver Gray



P7
Graphite Gray