

# Phonak RIC coupling guide

## Infinio




Phonak Audéo™  
R Infinio




Phonak Audéo  
Sphere™ Infinio

## Receivers (SDS 6.0)


Lengths: 00, 0, 1, 2, 3  
CeruStop wax protection




Standard (S)  
47/112 (2cc)



Medium (M)  
51/114 (2cc)



Power (P)  
59/122 (2cc)



UltraPower (UP)  
71/132 (2cc)  
cShell only



Phonak CROS  
R Infinio



CROS Tube 6.0  
(sizes 0-3)



CROS Tip 6.0

## Lumity



Phonak  
Audéo L-312




Phonak  
Audéo L-R/RT




Phonak  
Slim L-R

## Receivers (SDS 5)


Lengths: 0, 1, 2, 3  
CeruStop wax protection




Standard (S)  
47/111 (2cc)



Medium (M)  
51/114 (2cc)



Power (P)  
59/122 (2cc)



UltraPower (UP)  
71/132 (2cc)  
cShell only



Phonak CROS  
Lumity-R



CROS Tube  
(sizes 0-3)



CROS Tip

## Domes



Cap Dome



Open Dome



Vented Dome




Power Dome

## Custom ear pieces




SlimTip  
Acrylic / Silicone




cShell  
Acrylic / Titanium

## ActiveVent (MAV)

4.0 (Lumity)  
50/119 (2cc)



6.0 (Infinio)  
50/119 (2cc)



Universal SlimTip AV  
Sizes A, B, C and D



SlimTip  
Titanium AV

## Phonak Audéo Lumity Fit



Sensor Receiver (M)  
51/114 (2cc)



Sensor Receiver (P)  
59/122 (2cc)  
cShell only



Cap Dome



Open Dome



Clear Vented  
Dome



Sensor  
cShell

# Custom product guide

## Custom product models



Phonak Virto™  
R Infinio



Phonak Virto  
10 NW O Infinio



Phonak Virto  
Titanium Infinio



Phonak Virto  
P-312 NW O

## Shell colors



26  
Pink



14  
Tan



22  
Cocoa



28  
Brown



06  
Black



R  
Red  
Transparent



B  
Blue  
Transparent



21  
Transparent

## Faceplate colors



26  
Pink



14\*  
Tan



22  
Cocoa



28\*  
Brown



06  
Black


## Titanium shell color




Titanium  
gray

\*Virto I-R and Virto I-10 NW O only


## Fitting ranges




Mild to moderate  
(M) 40/109



Mild to moderately-severe  
(P) 50/115



Mild to severe  
(SP) 60/119



Moderate to severe/  
profound (UP) 70/127  
R and 312 models only

## Wax protection options




CeruStop  
wax filter

## SlimTip custom earpiece (RIC/BTE)



SlimTip  
Acrylic / Silicone

## cShell custom earpiece (RIC)



cShell  
Acrylic / Titanium

## Shell colors



26  
Pink



14  
Tan



22  
Cocoa



28  
Brown



06  
Black



R  
Red  
Transparent



B  
Blue  
Transparent



21  
Transparent

## Faceplate colors



26  
Pink



14  
Tan



22  
Cocoa



28  
Brown



06  
Black



21  
Transparent

## Wax protection options



CeruStop  
wax filter



H4 wax filter  
Default for UP receivers