

Phonak RIC coupling guide

Infinio Ultra




Audéo R




Audéo Sphere

Receivers (SDS 6.0)


Lengths: 00, 0, 1, 2, 3
CeruStop wax protection




Standard (S)
47/112 (2cc)



Medium (M)
51/114 (2cc)



Power (P)
59/122 (2cc)



UltraPower (UP)
71/132 (2cc)
cShell only



CROS R



CROS Tube 6.0
(sizes 0-3)



CROS Tip 6.0

Lumity




Audéo 312




Audéo R/RT

Receivers (SDS 5)


Lengths: 0, 1, 2, 3
CeruStop wax protection




Standard (S)
47/111 (2cc)



Medium (M)
51/114 (2cc)



Power (P)
59/122 (2cc)



UltraPower (UP)
71/132 (2cc)
cShell only



CROS R



CROS Tube
(sizes 0-3)



CROS Tip

EasyGuard™ Vented Domes (Infinio Ultra)



S



M



L

Domes



Cap Dome



Open Dome



Vented Dome



Power Dome

Custom ear pieces



SlimTip Acrylic



SlimTip Silicone




cShell Acrylic




cShell Titanium

ActiveVent (MAV)

4.0 (Lumity)
50/119 (2cc)



6.0 (Infinio)
50/119 (2cc)



Universal SlimTip AV
Sizes A, B, C and D



SlimTip
Titanium AV

Custom product guide

Infinio



New

Virtro R



Virtro 10 NW O



Virtro Titanium

Paradise



Virtro 312



Virtro 312 NW O

Shell colors



26
Pink



14
Tan



22
Cocoa



28
Brown



06
Black



R
Red
Transparent



B
Blue
Transparent



21
Transparent

Faceplate colors



26
Pink



14*
Tan



22
Cocoa



28*
Brown



06
Black


Titanium shell color




Titanium
gray

*Virtro I-R and Virtro I-10 NW O only


Fitting ranges




Mild to moderate
(M) 40/109



Mild to moderately-severe
(P) 50/115



Mild to severe
(SP) 60/119



Moderate to severe/
profound (UP) 70/127
R and 312 models only

Wax protection



CeruStop
wax filter

SlimTip custom earpiece (RIC/BTE)



SlimTip Acrylic



SlimTip Silicone

cShell custom earpiece (RIC)



cShell Acrylic



cShell Titanium

Shell colors



26
Pink



14
Tan



22
Cocoa



28
Brown



06
Black



R
Red
Transparent



B
Blue
Transparent



21
Transparent

Faceplate colors



26
Pink



14
Tan



22
Cocoa



28
Brown



06
Black



21
Transparent

Wax protection options



CeruStop
wax filter



H4 wax filter
Default for UP receivers