

Phonak

Charger Combi BTE

Technical data

Dimensions (L x W x H)	101.3 x 76.1 x 49 mm (4 x 3 x 1.9 ")
Weight	90 g (3.1 oz)
Color	Champagne
Operating conditions	5° to +40° Celsius (41° to 104° Fahrenheit)
Transport and storage conditions	-20° to +60° Celsius (-4° to +140° Fahrenheit)
Humidity during transportation	Up to 90% (non condensing)
Humidity during storage	0% to 70% (non condensing)
Humidity during operation	0% to 85% (non condensing)
Atmospheric pressure during use	800 hPa to 1060 hPa
Atmospheric pressure during storage and transportation	500 hPa to 1060 hPa
Included parts	Power supply Micro-USB cable (for charging only)
External interfaces	Micro-USB port
Battery	None
Charging time	See user guide
Input	Type: DC, stabilized 5 V +/- 5% Max. current: 150 mA

