

Phonak Sky™ Link L.

Product information



Product description



Sky Link L

Battery size	675	
LED pediatric patterns	•	
Telecoil	•	
Bluetooth®	•	
Hydrophobic coating	•	
IP rating	IP68 ¹	
Dimensions	L x W x D	42.6 x 18.8 x 8.6 mm (1.68 x 0.74 x 0.34")
	Weight	4.60 g / 0.162 oz

		HE11	HE11 680
Max. Power Output (dB SPL)	2 cc coupler	141	134
Max. gain (dB)	2 cc coupler	84	78
Frequency range (Hz)		<100 - > 4700	<100 - >4700
Working current (mA)		3.3	2.8

Fitting range



Housing and earhook colors

Housing



P1

Sand Beige



P4

Chestnut



P8

Velvet Black



Q3

Caribbean
Pirate



T3

Precious Pink



M7

Ocean Blue

Earhook



R1

Green



R2

Yellow



R3

Pink



R4

Orange



Q6

Blue



13

Transparent



Q7

Purple

¹ IP68 indicates that the hearing aid is water and dust resistant. It survived continuous immersion in 1 meter of fresh water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard.

² When using DSL fitting formula.

Performance levels

AutoSense Sky OS 5.0

Calm situation	•
Speech in noise	•
Comfort in noise	•
Music	•
Speech in loud noise	•
Speech in car	•
Comfort in echo	•
Media speech + mic	•
Media music + mic	•

Additional programs

Max. additional programs	4
Comfort in Echo	•
Speech in Loud Noise	•
Speech in 360°	•
Speech in Noise	•
Calm Situation	•
Comfort in Noise	•
Music	•
Acoustic phone	•
Custom program	•
Phone via T-coil ± mic	•
Public T-coil ± mic	•

Streaming programs

RogerDirect + mic	•
Phonak PartnerMic™ + mic	•
Phone call + mic	•

Features

SpeechSensor	•
Speech Enhancer	•
StereoZoom 2.0	•
Dynamic Noise Cancellation	•
RogerDirect	•
EchoBlock	20
SoundRelax	•
DuoPhone	•
Real Ear Sound	•
Roger™ and Directional setting	•
UltraZoom	•
Broadband Booster	•
WindBlock	•
User Preference Tuning	•
Auto Acclimatization	•
Environmental balance	•
SoundRecover2	•
SoundFinder	•
WhistleBlock	•
NoiseBlock	•
QuickSync	•
AOV	•
APD 2.0	•
Fine tuning channels	20

Compatibility

Sky Link L can be used:

- as a monaural fitting
- as a binaural fitting with two Sky Link L hearing aids
- as a bimodal fitting with an Advanced Bionics Sky CI L90 sound processor

Fitting

Phonak Target	11.2 or higher
Noahlink Wireless / Noahlink Wireless 2	•
NOAHLink	•
HI-PRO, HI-PRO2	•

Acoustic coupling

Sound outlet

Standard earhook	Damped (HE11 680), transparent and 6 additional pediatric colors Undamped (HE 11) only in transparent color
Tamper-proof ear-hook (mini)	Damped (HE11 680), transparent and 6 additional pediatric colors Undamped (HE 11) only in transparent color

Earpiece for ear-hook

Earmold	
---------	--

Wireless communication portfolio

Sky Link-L	
Wireless accessories	
TV Connector	•
Phonak PartnerMic	•
Phonak RemoteControl	•
Roger	
Roger microphones	•
Apps	
myPhonak Junior	•
AB Remote	•

myPhonak Junior is compatible with Sky Link L only if used in a monaural configuration or paired to another Sky Link L.

life is on

Since 1947, Phonak is dedicated to preserving life quality by opening new acoustic worlds. Already back then, in the days of our foundation, our company was driven by a formative conviction: We believe that well-hearing equates to well-being and thus is essential for living life to the fullest. In fact, the sense of hearing is directly linked to social, emotional, cognitive and physical well-being. Today as in future, we thrive to offer the broadest portfolio of innovative hearing solutions. And, together with our hearing care professionals, we keep on focusing on what matters most: improving speech understanding, changing people's lives and having a positive effect on society as a whole.

phonak.com



028-2844-02/V1.00/2025-11/jb © 2025 Sonova AG. All rights reserved.