

Phonak Sky™ Link M.

Product information



Product description




Sky Link M

Battery size	13	
Multi-function button	•	
Telecoil	•	
Bluetooth® 4.2 wireless technology	•	
Hydrophobic coating	•	
IP rating	IP68 ¹	
Dimensions	L x W x D	40.0 x 16.9 x 8.5 mm (1.57 x 0.67 x 0.33")
	Weight	4.40 g / 0.155 oz

		HE11	HE11 680
Max. Power Output (dB SPL)	2 cc coupler	139	133
Max. gain (dB)	2 cc coupler	81	75
Frequency range (Hz)		< 100 - > 5000	< 100 - > 5000
Working current (mA)		2.1	2.0

Fitting range



Mild to severe hearing loss, all audiometric configurations

■ HE11 680

Housing and earhook colors

Housing



P1

Sand Beige



P4

Chestnut



P8

Velvet Black



Q3

Caribbean
Pirate



T3

Precious Pink



M7

Ocean Blue

Earhook



R1

Green



R2

Yellow



R3

Pink



R4

Orange



Q6

Blue



13

Transparent



Q7

Purple

¹ IP68 indicates that the hearing aid is water and dust resistant. It survived continuous immersion in 1 meter of fresh water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Sonova AG is under license.

Performance levels

AutoSense Sky OS 3.0

Calm situation	•
Speech in noise	•
Comfort in noise	•
Music	•
Speech in loud noise	•
Speech in car	•
Comfort in echo	•
Media speech + mic	•
Media music + mic	•

Additional programs

Max. additional programs	4
Comfort in echo	•
Speech in loud noise	•
Speech in 360°	•
Speech in noise	•
Calm situation	•
Comfort in noise	•
Off-ear compatible*	•
Music	•
Acoustic phone	•
Custom program	•

Streaming programs

RogerDirect + mic	•
Phonak PartnerMic™ + mic	•
Phone call + mic	•

* Used in bimodal fittings only

Features

UltraZoom	•
SNR-Boost	•
DuoPhone	•
SoundRecover2	•
User Preference Tuning	•
Real Ear Sound	•
Fine tuning channels	20
WhistleBlock	•
NoiseBlock	•
WindBlock	•
EchoBlock	•
SoundRelax	•
QuickSync	•
Acoustically Optimized Vent	•
auto Acclimatization	•
Environmental balance	•
RogerDirect™	•
Roger and directional setting	•

Fitting

Phonak Target	7.0 or higher
Advanced Bionics Target CI	•
Noahlink Wireless / Noahlink Wireless 2	•
NOAHLink	•
HI-PRO, HI-PRO2	•

Acoustic coupling

Sound outlet	Standard earhook	Damped (HE 11 680) transparent and 6 additional pediatric colors (see page 2)
		Undamped (HE 11) transparent and 6 additional pediatric colors
	Tamperproof earhook	Mini damped (HE 11 680) transparent and 6 additional pediatric colors (see page 2)
		Mini undamped (HE 11) transparent and 6 additional pediatric colors
	SlimTube 4.0	Available in left and right and five different lengths (00, 0, 1, 2, 3)
Earpiece for SlimTube	SDS 4.0	Power dome and all other domes SDS 4.0 available in three sizes (S, M, L)
	SlimTip	Acrylic
		Silicone
Earpiece for earhook	Earmold	
	Universal eartip	

Phonak Trial™

If a trial device is required, Phonak Sky M-SP Trial device is recommended. This trial device is not compatible with the AB sound processor.

Compatibility

Sky Link M can be used:

- as a monaural fitting
- as a binaural fitting with two Sky Link M hearing aids
- as a bimodal fitting with an Advanced Bionics Sky CI M90 sound processor

Wireless communication portfolio

Sky Link-M	
Wireless accessories	
TV Connector	•
Phonak PartnerMic	•
Phonak RemoteControl	•
Roger	
Roger microphones	•
Apps	
myPhonak Junior	•
AB Remote	•

myPhonak Junior is compatible with Sky Link M only if used in a monaural configuration or paired to another Sky Link M.

life is on

Since 1947, Phonak is dedicated to preserving life quality by opening new acoustic worlds. Already back then, in the days of our foundation, our company was driven by a formative conviction: We believe that well-hearing equates to well-being and thus is essential for living life to the fullest. In fact, the sense of hearing is directly linked to social, emotional, cognitive and physical well-being. Today as in future, we strive to offer the broadest portfolio of innovative hearing solutions. And, together with our hearing care professionals, we keep on focusing on what matters most: improving speech understanding, changing people's lives and having a positive effect on society as a whole.

phonak.com



027-0606-02/V3.00/2024-08/jb © 2020-2024 Sonova AG All rights reserved

