

Phonak Paradise Receiver-in-the-Canal hearing aids

Unleash the wonder within

All with Phonak Paradise's unrivaled* sound quality¹, universal connectivity, brilliant speech understanding^{1,2,3} and crisp natural sound⁴

Phonak Audéo™



- Favorite dinner and date companion
- Batteries, t-coil, colors -- more options than a Swiss Army knife
- Power to perform

Phonak Audéo Fit™



- I ♥ you
- Supercharged app!
- 10,000 step champion
- Wear it with your fine Swiss timepiece

Phonak Audéo Life™



- No umbrella needed
- Peace of mind
- Splash in the pool like a kid again
- Charging made easy

Phonak Slim™



- Selfie ready
- A perfect couple
- Like your favorite pair of shoes
- Not your traditional hearing aid

*Compared to Phonak Audéo™ Marvel

1. Appleton, J. (2020) AutoSense OS 4.0 - significantly less listening effort and preferred for speech intelligibility. Phonak Field Study News retrieved from www.phonakpro.com/evidence, accessed November 2020.
2. Appleton, J., & Voss S.C., (2020) Motion-based beamformer steering leads to better speech understanding and overall listening experience. Phonak Field Study News. retrieved from www.phonakpro.com/evidence, accessed August 2020.
3. Wright, A. (2020). Adaptive Phonak Digital 2.0: Next-level fitting formula with adaptive compression for reduced listening effort. Phonak Field Study News, retrieved from www.phonakpro.com/evidence, accessed August 19th 2020.
4. Taphuntsang, D. (2020). Market research ID 4387. Please contact marketinsight@phonak.com if you are interested in further information.

A Sonova brand

PHONAK
life is on

Tailored to your lifestyle



	Audéo	Audéo Fit	Audéo Life	Slim
Automatically adapt to multiple environments with AutoSense OS™ 4.0	•	•	•	•
Reduced listening effort in noise for those with mild to severe hearing loss ^{1,3}	•	•	•	•
Automatically distinguish between streamed speech and music	•	•	•	•
Tap Control (available only in 90 & 70 rechargeable models)	•	•	•	•
Universal Bluetooth® connectivity (iOS® & Android™)	•	•	•	•
RogerDirect™ compatibility	•	•	•	•
Rechargeable option	•	•	•	•
Zinc-air battery option	•			
CROS compatibility	•			
UP receiver compatibility	•		•	
ActiveVent™ Receiver compatibility	•		•	•
Tracks average wearing time in different sound environments (rechargeable models)*	•	•	•	•
Tracks steps and activity levels with optional goal setting (rechargeable models)*	•	•	•	•
Tracks heart rate and distance walked or ran*		•		
Tested beyond IP68. Waterproof (pool water, salt water, rain) up to 50 cm and sweatproof			•	
Unique 7° angle which follows the natural curve of the head				•

*with the myPhonak app. Some features may require user registration.

iOS is a trademark of Cisco Technology, Inc. Android is a trademark of Google LLC.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Sonova AG is under license.



 Sonova AG
Laubisrütistrasse 28
CH-8712 Stäfa · Switzerland

 Sonova Deutschland GmbH
Max-Eyth-Str. 20
70736 Fellbach-Oeffingen · Germany



028-2358-02/V2.00/2022-01/dr © 2022 Sonova AG All rights reserved