

Phonak ITE maintenance kit.

The kit contains:

- 10 wipes
- 1 cleaning spray 30ml
- 2 drying capsules
- 1 drying cup



Cleaning wipe:

The cleaning wipe is a cleaning product for the maintenance of hearing aids and earmolds. It removes earwax responsible for the decrease in sound quality, increases the life of hearing systems and preserves their effectiveness.



Cleaning spray:

The cleaning spray is a cleaning product for the maintenance of hearing aids and earmolds. It removes earwax responsible for the decrease in sound quality, increases the life of hearing systems and preserves their effectiveness.



Drying capsule:

The drying capsule is indicated for dehumidifying hearing aids and earmolds, in order to avoid breakdown due to moisture.



Drying cup:

The drying cup is an essential accessory to be used with drying capsules. It ensures a hermetic environment in order to guarantee efficient dehumidification of hearing aids thanks to the action of dehydrating capsules.



029-1355-02 | V1.00 | 2023-12 | MMW © 2023 Sonova AG All rights reserved