

OGA Phonak Infinio and Paradise Custom Digital Hearing Aids (36C79125D0007)

Government Services | Phone 888-561-7555 | Fax 630-836-9770

Customer Information

Payer:

- ☐ ROES (DALC PH5120018626)
- ☐ Other
- ☐ RACHAP (Enclose payment and if paying by cc, RACHAP authorization form)

Ship To:

Address:

City: State: Zip:

PO #:

Contact Name:

Phone Number:

Email Address:

Patient's Last Name:

Patient's 4 Digit SSN:

Age:

Gender:

Step 1: Select Product

Paradise Wireless (312 battery)



Virto P90 FS

☐ L ☐ R



Virto P90 HS

☐ L ☐ R



Virto P90 CA

☐ L ☐ R

Infinio Non-Wireless (10A battery)



Virto I90 CIC

☐ L ☐ R



Virto I90Titanium IIC

☐ L ☐ R

Step 2: Select Power Level

M ☐ L ☐ R P ☐ L ☐ R SP ☐ L ☐ R UP¹ ☐ L ☐ R

Step 3: Colors

Shell

- ☐ Black (06)¹
- ☐ Transparent (21)
- ☐ Cocoa (22)
- ☐ Pink (26)
- ☐ Tan (14)
- ☐ Brown (28)
- ☐ Red/Blue (R/B)³
- ☐ Titanium gray (U0)²

Faceplate

- ☐ Black (06)
- ☐ Cocoa (22)
- ☐ Pink (26)
- ☐ Tan (14)
- ☐ Brown (28)

Step 4: Options

For selection options and codes see combination matrix on the back side.

Legend

☐ / ☐ S = Standard

☐ = Optional

Vent size	AOV (enter audiogram)	☐ S
	AOV-O (Virto I-Titanium only)	☐
	Manually selected	Left: Right:
User controls	Standard controls	☐ S
	Push button	☐
	Volume control (VC)	☐
	T-coil (IIC/CIC only)	☐
	Other	
Wax system	Cerustop (CS)	☐ S
	Manually selected	Type:
Removal line (RF)	Transparent (21) ³	☐
Titanium surface	Extra retention (ERS)/Matte	☐
Acrylic/retention options	See matrix for options	
Accessories (Virto P FS, HS, CA only)	TV Connector 2.0	☐
	Phonak PartnerMic™	☐
	Phonak RemoteControl	☐
	Roger On™ V3	☐

1 Only available for Virto P FS, HS, CA

2 Only available for Virto I90Titanium IIC

3 Standard for IIC and CIC

Combination matrix

		Virto I90Titanium IIC	VirtoI90 CIC	Virto P90 CA	Virto P90 HS	Virto P90 FS
						
Battery Size		10	10	312	312	312
Wireless		No	No	Yes	Yes	Yes
Push button	Standard	No push button (PX)	Push button (PM1) ⁶	Push button (PM1)	Push button (PM1)	Push button (PM1)
	Optional	Push button (PM1)	No push button (PX)	-	-	-
Volume Control	Standard	-	No VC (VCX)	No VC (VCX)	No VC (VCX)	No VC (VCX)
	Optional	-	Digital VC (VC3)	Digital VC (VC3)	Digital VC (VC3)	Digital VC (VC3)
T-Coil ⁵	Standard	No T-coil (TX)	No T-coil (TX)	-	-	-
	Optional	T-coil (TP)	T-coil (TP) ⁶	-	-	-
Magnetic switch function	Standard	MiniControl (DT5)	EasyPhone (DT4) ⁶	-	-	-
	Optional	EasyPhone (DT6)	MiniControl (DT3) No EasyPhone (RX)	-	-	-
Wax system	Standard	Cerustop (CS) – required	Cerustop (CS)	Cerustop (CS)	Cerustop (CS)	Cerustop (CS)
	Optional	-	Extended rec. tube (ER) Wax Spring (WS) Ext. Tube + Spring (WE) Wax guard not required (X)	Extended rec. tube (ER) Wax Spring (WS) Ext. Tube + Spring (WE) Wax guard not required (X)	Extended rec. tube (ER) Wax Spring (WS) Ext. Tube + Spring (WE) Wax guard not required (X)	Extended rec. tube (ER) Wax Spring (WS) Ext. Tube + Spring (WE) Wax guard not required (X)
Removal line	Standard	Removal line transparent (RF) ⁴	Removal line transparent (RF) ⁴	-	-	-
	Optional	-	Removal line transparent (RF)	Removal line transparent (RF)	Removal line transparent (RF)	Removal line transparent (RF)
Acrylic / retention options	Optional	-	Canal lock ⁵ Helix lock ⁵ Raised VC	Canal lock ⁵ Helix lock ⁵ Raised VC	Canal lock ⁵ Helix lock ⁵ Raised VC	Helix lock ⁵ Raised VC

4 Standard for IIC and CIC

5 Color will match shell color

6 Not available in combination with VC

Step 5: Audiometric Data and Preprogramming

Audiometric Information

Please include with all preprogramming requests.
Enter 250-4000Hz for AOV Venting

L						
R						
	250	500	1000	2000	3000	4000

Preprogramming Options

- ☐ ☐ Long Term User [EXL04]
☐ ☐ NAL NL1 Fitting Formula [FIT03]
☐ ☐ Add Speech in Loud Noise Program [STZP0]
☐ ☐ Add Speech in 360° [ZMCP0]
☐ ☐ T+Mic in PRG1 [TCPP0]
CIC/IIC only
☐ ☐ WhistleBlock on [WHB23]

- ☐ ☐ NAL NL2 Fitting Formula [FIT04]
☐ ☐ Add Tinnitus Balance Program [TPMP0]
☐ ☐ Add Acoustic DuoPhone Program [ADMP0]
☐ ☐ Adaptive Phonak Digital Contrast [FIT05]
☐ ☐ Change Soft Noise Reduction to Weak [SNR04]
☐ ☐ Disable Sound Recover [SRB00]

Step 6: Special Instructions and Notes

- ☐ Call if changes are required
☐ No phone call required to make changes