

Charger compatibility guide.

Phonak Charger
Virto Infinio



075-0075-11

Phonak ChargerGo
RIC Sphere Infinio



075-0057-11

Phonak ChargerGo
RIC Infinio



075-0055-11

Phonak Charger
Ease



075-0080-11

Phonak Charger
Case Go™



075-3017-11

Phonak Charger
Case Combi 2



075-0063-11

Phonak Charger
Combi BTE 2



075-0065-11

Phonak Virto™ I-R	•					
Phonak Audéo™ I-R Sphere		•				
Phonak Audéo I-R			•			
Phonak CROS I-R			•			
Phonak Audéo L-R/RT				•		
Phonak Audéo L-R Life					•	
Phonak Audéo L-R Fit				•		
Phonak CROS L-R				•		
Phonak Naída™ L-PR						•
Phonak Audéo P-R/RT					•	
Phonak Audéo P-R Life					•	
Phonak Audéo P-R Fit					•	
Phonak CROS P-R					•	
Phonak Naída P-PR						•
Phonak Audéo M-R/RT					•	
Phonak Bolero M-PR						•

• 071-0003 Power Pack compatible

Power Pack is not recommended to be charged with third-party battery chargers (power banks) and should only be charged with the combi case power supply. Power Pack always ships ground.



Phonak
Government Services

Charger compatibility guide.

	Phonak Charger Virto Infinio	Phonak ChargerGo RIC Sphere Infinio	Phonak ChargerGo RIC Infinio	Phonak Charger Ease	Phonak Charger Case Go™	Phonak Charger Case Combi 2	Phonak Charger Combi BTE 2
							
	075-0075-11	075-0057-11	075-0055-11	075-0080-11	075-3017-11	075-0063-11	075-0065-11
External Power Supply (EPS) 5V US – no cable included 075-0030-11	•	•	•	•	•	•	•
USB A to USB C cable 043-3116	•	•	•	•			
USB microUSB 1m 28AWG No data 043-0002						•	•
USB A to USB C cable with ferrite bead 043-3116-1					•		

Power Supply

External Power Supply (EPS) 5V 1A US	075-0030-11	
---	-------------	---

USB A to USB C with ferrite	043-3116-1	
-----------------------------	------------	---

Cable USB Micro USB 1m 28AWG No data	043-0002	
---	----------	---

USB A to USB C	043-3116	
----------------	----------	---