

# Phonak RIC options.

## Audéo I90-R Sphere

The go-to RIC for your patients



- Exceptional sound quality thanks to APD 3.0 & AutoSense OS™ 6.0
- World's first hearing aid with dedicated deep neural network (DNN) chip
- Instant speech separation from noise all around — with an unprecedented 10dB SNR improvement<sup>1</sup>
- Unprecedented speech clarity from any direction in background noise<sup>2</sup>
- 34% less listening effort with speech from a distance<sup>3</sup>

## Specialty Solutions

### Audéo I90-R\*

Standard Rechargeable




### Audéo L90-RT

T-coil



### Audéo L90-312

Disposable battery



## Universal Connectivity

- iPhone®, Android™, flip phones, etc.
- Computers, tablets, etc.
- Widget feature
- 2 Active connections
  - i.e. 2 phones, phone & tablet, phone & computer, etc.

\* Compatible with Infinio CROS

<sup>1</sup> Rauber, S., Kohlhauer, P., Jehle, F., Kühnel, V., Preuss, M., Hobi, S. (2024). Spheric Speech Clarity proven to outperform three key competitors for clear speech in noise. Phonak Field Study News retrieved from [www.phonak.com/evidence](http://www.phonak.com/evidence)

<sup>2</sup> Wright, A., et al "Spheric Speech Clarity applies DNN signal processing to significantly improve speech understanding from any direction and reduce the listening effort." Phonak Field Study News retrieved from [www.phonak.com/evidence](http://www.phonak.com/evidence)

<sup>3</sup> Latzel, M, Heeren, J and Lesimple, C. (2024) „Speech Enhancer reduces listening effort and fatigue." Phonak Field Study News retrieved from [www.phonak.com/evidence](http://www.phonak.com/evidence)

iPhone is a trademark of Apple Inc., registered in the U.S. and other countries and regions. Android is a trademark of Google LLC.

