

Phonak Government Services Selection Guide | November 2025

Hearing Aids

RIC

Phonak Audéo™ I90




Phonak Audéo L90




Fitting Ranges


S Receiver




M/M Sensor Receiver



P Receiver




UP Receiver
N/A on Slim



Custom


Phonak Virto™ I90




Phonak Virto P90




All styles
40/109 (M)




All styles
50/115 (P)



All styles
60/119 (SP)






FS, HS, CA (312 battery)
70/127 (UP)




*Wireless connectivity

BTE

Phonak Naída™ L90



Phonak Naída Link M



Phonak Government Services Selection Guide | November 2025

CROS



CROS I-R
Compatible with Audéo I90-R



CROS L-R
Compatible with Audéo L90-R/RT
and Naïda L90-PR

Wireless Accessories



TV Connector 2.0



Phonak
PartnerMic™



Phonak
RemoteControl

Roger On™ V3

All-in-one microphone
and media streaming device



Colors

Infinio/Lumity



H0
Beige



P1
Sand Beige



P3*
Sandalwood



P4
Chestnut

* N/A for Audéo
or CROS Lumity



P5
Champagne



P6
Silver Gray



P7
Graphite Gray



P8
Velvet Black

Naïda Link M



P1
Sand Beige



P4
Chestnut



P6
Silver Gray



P8
Velvet Black

Custom shells



26
Pink



14
Tan



22
Cocoa



28
Brown



06
Black



R
Red
Transparent



B
Blue
Transparent



21
Transparent

Custom faceplates



26
Pink



14*
Tan



22
Cocoa



28*
Brown



06
Black

* N/A on Virto Titanium

Virto Titanium – Shell color U0 Titanium Gray Et faceplate colors



26
Pink



22
Cocoa



06
Black

Coupling/Earmolds



Cap
Dome



Open
Dome



Vented
Dome



Power
Dome



SlimTip
(Acrylic)



SlimTip
(Silicone)



cShell



cShell
Titanium



MAV
TITNM
SlimTip



Phonak
CROS



Phonak
Skeleton



Phonak
Semi-Skeleton



Phonak
Canal



Phonak
Canal Lock



Phonak
Helix Lock



Phonak
Half Shell



Phonak
Full Shell
Carved



Phonak
Standard
Full Shell