



Phonak Virto™ 1.

Product information for non-wireless models

Powered by the innovative technology of the Phonak Infinio platform, Virto I hearing aids are the solution for those clients that prefer a discreet, custom-made solution.



Product description

	 Virtuo I-10 NW O 190/170/150/130	 Virtuo I-Titanium 190/170
Power level		
M	•	•
P	•	•
SP	•	•
Hardware configuration / options / technology		
Battery size	10	10
Shell material	Acrylic	Titanium ¹
IP rating	IP68 ²	IP68 ²
Microphone technology	Omni	Omni
Telecoil	Optional	Optional
Push button (Program / Volume)	Optional	Optional
Volume control (wheel)	Optional	
EasyPhone (for acoustic phone)	•	Optional
MiniControl (for program or volume change)	•	•

Colors

Faceplate colors



26
Pink



14³
Tan



22
Cocoa



28³
Brown



06
Black

Titanium shell color



Titanium
gray

Acrylic shell colors



26
Pink



14
Tan



22
Cocoa



28
Brown



06
Black



R
Red Transparent



B
Blue Transparent



21
Transparent

¹ Available in two different surface roughnesses: standard or extra retention (ERS)

² IP68 indicates that the hearing aid is water and dust resistant. It survived continuous immersion in 1 meter of fresh water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard.

³ Virtuo I-10 NW O only

Performance levels

AutoSense OS™ 6.0	Premium I90	Advanced I70	Standard I50	Essential I30
Music	•	•		
Speech in Car	•	•	•	
Comfort in Noise	•	•	•	
Speech in Noise	•	•	•	•
Calm Situation	•	•	•	•
Additional programs				
Max. additional programs	4	4	4	4
Speech in car	•	•	•	
Speech in Noise	•	•	•	•
Calm Situation	•	•	•	•
Comfort in Noise	•	•	•	•
Music	•	•	•	•
Acoustic phone	•	•	•	•
Custom program	•	•	•	•
Phone via T-coil ± mic	•	•	•	•
Public T-coil ± mic	•	•	•	•




Features	Premium I90	Advanced I70	Standard I50	Essential I30
Speech Enhancer	•			
SoundRelax	•	•	•	
Roger compatibility ¹	•	•	•	•
EasyPhone	•	•	•	•
WindBlock	•	•	•	•
User Preference Tuning	•	•	•	•
Tinnitus Balance	•	•	•	•
Auto Acclimatization	•	•	•	•
Environmental balance	•	•	•	•
SoundRecover2	•	•	•	•
WhistleBlock	•	•	•	•
NoiseBlock	•	•	•	•
AOV	•	•	•	•
Fine tuning channels	20	20	16	12

* Available via Roger Neckloop when optional telecoil is installed



Fitting

	Virto I-10 NW O / Virto I-Titanium
Phonak Target	10.0 or higher
NOAHlink	•
HI-PRO, HI-PRO2	•

Acoustic coupling

		M	P	SP
Maximum Power Output (dB SPL)	2 cc coupler	109	115	119
Maximum gain (dB)	2 cc coupler	40	50	60
Frequency range (Hz)	2 cc coupler	<100 - 7000	<100 - 6700	<100 - 7000
Working current (mA)		1.1	1.2	1.1
Fitting range		 <p>For mild to moderate hearing loss, all audiometric configurations</p>	 <p>For mild to moderately-severe hearing loss, all audiometric configurations</p>	 <p>For moderate to severe hearing loss, all audiometric configurations</p>

Shell styles

	 Virto I-10 NW O	 Virto I-Titanium
IIC	•	•
CIC	•	•
MC	•	
ITC	•	
HS	•	
FS	•	

Wireless communication

	Virto I-10 NW O I90/I70/I50/I30	Virto I-Titanium I90/I70
Roger		
Roger Neckloop	•*	•*

* Available if telecoil is installed

life is on

Since 1947, Phonak is dedicated to preserving life quality by opening new acoustic worlds. Already back then, in the days of our foundation, our company was driven by a formative conviction: We believe that well-hearing equates to well-being and thus is essential for living life to the fullest. In fact, the sense of hearing is directly linked to social, emotional, cognitive and physical well-being. Today as in future, we thrive to offer the broadest portfolio of innovative hearing solutions. And, together with our hearing care professionals, we keep on focusing on what matters most: improving speech understanding, changing people's lives and having a positive effect on society as a whole.

phonak.com



027-0710-02/N3.00/2025-01/na © 2025 Sonova AG All rights reserved

