

Playful design

With an exclusive Mix & Match concept incorporating seven colors for earhooks and 14 colors for hearing aid housings, we give every child the opportunity to design a hearing aid that truly reflects their personality.



Life is on

At Phonak, we believe that hearing well is essential to living life to the fullest. For more than 70 years, we have remained true to our mission by developing pioneering hearing solutions that change people's lives to thrive socially and emotionally. Life is on.

www.phonak.com/kids

Phonak Sky™ B-PR



When a child can enjoy 24 hours of hearing with one simple charge, life is on

- ¹ Flexer, C. (2018). The ears are doorways to the brain. *Phonak Insight*, retrieved from www.phonakpro.com/evidence, accessed February 19th, 2018.
- ² Rakita, L. (2016). AutoSense OS: Hearing well in every listening environment has never been easier. *Phonak Insight*. Retrieved from www.phonakpro.com/evidence, accessed February 19th, 2018.
- ³ Latzel, M. (2016). AutoSense OS - Benefit of the next generation of technology automation. *Phonak Field Study News*, retrieved from www.phonakpro.com/evidence, accessed February 19th, 2018.
- ⁴ Stephenson, B., Omisore, D., Pedersen, H. M., & Taphuntsang (2016). The Phonak rechargeable solution. *Field Study News*, retrieved from www.phonakpro.com/evidence, accessed February 19th, 2018.

sonova
HEAR THE WORLD



028-1767-02V1.00|2018-02|na Printed in XXXX © 2018 Sonova AG All rights reserved

A Sonova brand

PHONAK
life is on

Phonak pediatric solutions:

The key to unlocking a child's full potential

We tend to think that we hear with our ears, but science confirms that we really hear with the brain.¹ Our ears help sounds reach the brain and as such, can be considered 'doorways' to the brain.

Hearing loss can be thought of as a 'doorway problem', an obstruction that prevents vital auditory information from reaching the brain.

Phonak has developed unique pediatric evidence-based solutions that open these obstructed doorways and provide clear and clean signals, thereby giving your child access to the millions of words and thousands of hours of listening required for spoken language, literacy, and learning in general.

Maximized hearing performance made simple

At Phonak, we believe in simplifying lives. We understand your child's hearing aids must automatically provide excellent sound quality, seamlessly adjust to every environment, and fit in with your child's active life.^{2,3}

Keeping kids in charge all day long

Sky B-PR is the first rechargeable hearing aid designed specifically with children in mind. Equipped with a built-in lithium-ion battery, it simplifies life by offering first-of-its-kind benefits in a pediatric hearing aid.⁴

- **Quick-charging** – Fully charged in just three hours. A quick 30-minute charge will provide up to six hours of battery life
- **Long-lasting** – A fully charged battery will last the whole day,* even with up to 10 hours of Roger™** use during the day
- **Easy to use** – Automatically adjusts to your child's listening environment
- **Hassle-free** – No batteries to purchase or change

Always on

Sky B-PR comes with a choice of easy-to-use charging options:

- Charger Case
- Mini Charger
- Power Pack



Performance where kids need it the most

Phonak Sky B-PR offers the dedicated features of Sky B hearing aids, including AutoSense Sky OS and adaptive SoundRecover2, so you can rest assured your child is getting outstanding performance wherever and whenever needed. When combined with Roger technology, children can master the complexity of everyday listening situations.

Safe and strong

- **Water resistant** (IP68 rating)*** – No need to remove your child's hearing aids when life gets a bit wet
- **Robust** – Tough enough to keep up with the most active children
- **Safety** – Tamperproof earhooks safe for curious hands, plus the built-in battery cannot be removed
- **Status notification** – Easy-to-understand indicator light confirms that the hearing aids are powered on, the Roger program is active, and warns of a low battery

* 16 hours of hearing with 10 hours of Roger use or audio streaming

** Roger technology can be accessed via the ComPilot II with a Roger X or via Roger MyLink

*** IP68 indicates that the hearing aid is water resistant and dust tight. It survived continuous immersion in 1 meter of water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard, no traces of dust were evident within the housing