


# Your hearing. Your solution.



Name

Date

# Your audiogram



## Hearing aid styles

### ITE (In-the-Ear)



### RIC (Receiver-In-Canal)



### BTE (Behind-the-Ear)



Available in a host of hair, skin, fashion and fun colors.

# Select your hearing goals


## Premium (P)



- Follow discussions in the car



- Enjoy conversations at noisy social events



- Engage in outdoor sports activities



- Enjoy a visit to the museum

## Advanced (A)



- Hear voices from any direction



- Enjoy fuller and richer music



- Enjoy conversations at restaurants

## Standard (S)



- Converse with colleagues



- Chat on the phone

## Essential (E)



- Enjoy personal discussions



Your recommended performance level \_\_\_\_\_

Notes \_\_\_\_\_

---


---

---

---

# Phonak programs and features overview

**AutoSense OS**, the operating system within Phonak Venture hearing aids, has been developed to automatically adapt to your sound environment. Whether you're in a noisy restaurant, car, concert hall or simply sitting at home, AutoSense OS makes listening effortless and doesn't require any manual interaction. Improved speech understanding in various situations will keep you connected throughout the day, making interacting and engaging a real pleasure.

		P	A	S	E
AutoSense OS programs	 <b>Comfort in Echo</b> Reduces distracting echo, making listening in echoey environments such as large public buildings more comfortable.	•			
	 <b>Speech in Car</b> Listening effort is reduced so you can enjoy the journey and the conversation.	•			
	 <b>Speech in Loud Noise</b> Provides understanding of a single voice in front of you, for example while gathering with your family in a restaurant.	•			
	 <b>Music</b> Offers a full and rich music experience so you can enjoy listening to your favorite music.	•	•		
	 <b>Comfort in Noise</b> Uncomfortable loud sounds are reduced enhancing your listening comfort.	•	•	•	
	 <b>Speech in Noise</b> Noise is reduced from the back allowing you to focus on the conversation in front of you.	•	•	•	•
Additional programs	 <b>Calm Situation</b> Creates a comfortable listening environment to simply relax or enjoy a one to one conversation.	•	•	•	•
	 <b>Speech in Wind</b> Allows you to enjoy conversations even on a windy day, for example while playing golf.	•			
	 <b>Speech in 360°</b> Automatically selects the direction of speech, improving speech intelligibility without the need for you to face the speaker. Follow conversations for example in a business meeting around a large table or while driving.	•	•		
	 <b>Speech in Loud Noise</b> Provides understanding of a single person speaking directly in front of you in a noisy crowd, for example while gathering with your family in a restaurant.	•	•		
Signature features	 <b>Acoustic phone</b> Makes understanding on the phone easy by streaming the caller's voice to both ears, any time, any where. Talk on your home or cell phone.	•	•	•	
	 <b>SoundRecover</b> Makes high-pitched sounds easier to hear in all situations, for example birds singing, children's voices or the ring of a doorbell.	•	•	•	•
	 <b>WhistleBlock</b> Eliminates annoying whistling from hearing aids so you can keep focused on the conversation without worry of interruption.	•	•	•	•
	 <b>NoiseBlock</b> Reduces uncomfortable background noise such as traffic or the constant hum of a ceiling fan without affecting your ability to understand.	•	•	•	•

For a full list of Phonak program and feature benefits, please visit [www.phonak.com](http://www.phonak.com) or speak to your hearing care professional.



# Phonak Wireless Communication Portfolio

---

## Control

Easy to use remote control for your hearing aids and an App for discreet access.



PilotOne  
PilotOne II



RemoteControl App

## Phone

Talk on the phone at home and on your cell phone with ease.



DECT  
DECT II



EasyCall  
EasyCall II

## Connectivity

Wirelessly stream TV, music, and conversations in stereo quality.



ComPilot Air II



ComPilot  
ComPilot II



TVLink II

## Noise

Understand speech over distance in groups or one-on-one in the most challenging, noisy environments.



Roger Pen  
Roger EasyPen



Roger Clip-On Mic



Roger receiver

## Distance

Enhance one-to-one conversations over distance.



RemoteMic

# Phonak has the right hearing solution for you – regardless of your type of hearing loss, lifestyle, personal preferences, age or budget.

## Premium

Developed to satisfy those who desire and expect nothing but the best. Thanks to a complete set of trailblazing solutions that uniquely address all hearing challenges, these hearing aids support understanding in even the most difficult listening situations, for example at a busy social event.

## Advanced

All the performance levels have the perfect balance between performance and price. An impressive offering of features that enhances hearing ability in a variety of listening environments, including those with moderate levels of background noise, for example at a restaurant.

## Standard

Ideally suited for people who expect value in state-of-the-art technology. A competitive feature set that has the capability to support the hearing demands of daily life with minimal to moderate levels of background noise, for example using the telephone in a public place.

## Essential

The Essential products are designed for people looking for an affordable yet effective solution to enhance daily communication abilities, mostly in quiet or minimal background noise, for example one to one conversations.

## Choose to enjoy better hearing that matches your lifestyle!

### Next appointment

---

Date

---

Time

This is your personal record of your solution.  
Please bring it with you to every appointment.  
Feel free to add notes of your own.

Your hearing care professional: