





# C&C Line





Hearing system cleansing & care












## Cleaning and Care Kits

Product	Description	Content	Part.No.	
	<b>Kit K1 BTE</b> with standard tube/earmold	Cleansing, disinfection and care kit for behind-the-ear hearing aids with earmold	Phonak Hygiene bag with Dry Cup with sieve, 2 C&C Dry Cap, Tool-Set, 10 C&C Wipe, Air blower, 6 C&C Clean Tab, C&C Spray 18 ml, 1 Phonak microfiber cloth, Manual drying, Manual Cleansing	098-0041
	<b>Kit K2 ITE and BTE</b> with Slim/Microtube	Cleansing, disinfection and care kit for in-the-ear and behind-the-ear hearing aids with slim tube and dome	Phonak Hygiene bag with C&C Dry Cup with sieve, 2 C&C Dry Cap, Tool-Set, 10 C&C Wipe, 1 C&C Spray 18 ml, 1 Phonak microfiber cloth, Manual Drying, Manual Cleansing	098-0042
	<b>Kit K3 ITE and BTE</b> non disinfectant, drying only	Care kit for in-the-ear and behind-the-ear hearing systems	Phonak Hygiene bag with C&C Dry Cup with sieve, 4 C&C Dry Cap, Tool-Set, 1 Phonak microfiber cloth, Manual Drying	098-0396
	<b>Kit K4 BTE</b> earmold	Non disinfectant, drying and cleansing only Cleansing and Care kit for behind-the-ear hearing aids and earmolds	Phonak Hygiene bag with C&C Dry Cup with sieve, 4 C&C Dry Cap, 8 C&C Clean Tab, Air Blower, Tool-Set, Manual Drying, Manual Cleansing	098-0397

## Refills

Product	Recommended application	Part.No.	
	<b>C&amp;C Wipe</b> (25 tissues) Disinfectant wipes for hearing systems and earmolds	BTE: every 2 days - lasts for 50 days ITE: every day - lasts for 25 days	098-0012
	<b>C&amp;C Clean Tab</b> (20 tablets) Special cleansing tablets for earmolds	BTE: every 3 days to clean the earmold - lasts for 60 days ITE: not for ITE use!	098-0013
	<b>C&amp;C Dry Cap</b> (6 capsules) Capsules for the drying of hearing systems and earmolds	BTE: every day - a capsule lasts for approx. 2 weeks ITE: every day - a capsule lasts for approx. 2 weeks	098-0014
	<b>C&amp;C Spray</b> (18 ml) Disinfectant spray for hearing systems and earmolds	BTE: every 2 days - lasts for approx. 80 applications ITE: every day - lasts for approx. 80 applications	098-0032

# Accessories

Product		Part.No.
	Phonak Hygiene Bag	017-0372
	Phonak MiniKit Neoprene Pouch	098-0033
	Microfiber cloth red (12 x 17cm)	098-0034
	Air Blower color: dark blue, with small metal tip	098-0018
	C&C Dry Cup, with sieve	098-0194
	C&C Dry Bag	098-0017
	C&C Clean Cup, with sieve	098-0108
	C&C Wipe Dispenser (90) Disinfectant wipes for hearing systems and earmolds	098-0357
	Tool Set	098-0016