

Phonak Ladegeräte – Kompatibilitätsübersicht.

	ChargerGo RIC I EON	Charger RIC EON	ChargerGo RIC Sphere Infinio	ChargerGo RIC Infinio	Charger RIC Infinio	Charger Virto I	Charger Ease	Charger Combi BTE 2*
								
	075-0108-12	075-0073-12	075-0057-12	075-0055-12	075-0054-12	075-0075-12	075-0080-10	075-0065-12
Audéo E-Sphere	•	•						
Audéo E-R	•	•						
CROS E-R	•	•						
Audéo I-Sphere			•					
Audéo I-R				•	•			
CROS I-R				•	•			
Virto I-R						•		
Audéo L-R/RT							•	
CROS L-R							•	
Naída L-PR								•
Sky L-PR								•
Roger Focus II								•

* Die gleiche Kompatibilität gilt für Phonak Charger Combi BTE (075-0062-12)

 Kompatibel mit Phonak Power Pack (071-0003)