

Phonak Audéo P RIC Order Form

For mild to severe hearing losses

Rechargeable and zinc air RIC models capable of connecting to a broad range of devices via universal Bluetooth® and Phonak's Airstream technology.

PHONAK
life is on

Step 1: Order details

Account Number : _____

Hearing Clinic: _____

Delivery Address: _____

Contact Name: _____ Phone Number: _____

Email Address: _____

Date Ordered: _____ Date Required: _____

Step 2: Patient information

Name: _____

Audiogram: (Required for Intellivent option)

	kHz	0.25	0.5	1	2	4
LEFT AC (dB HL)						
RIGHT AC (dB HL)						

Step 3: Warranty/VIP service

4th year (\$160 per device)

VIP 24 hour service (\$50 ex GST) (VIP not available for Titanium products)

Step 4: Choose hearing instrument (see following page for mould order, HSP codes and colour options)

L	R	Model	Part No.	Colour
		P90-312	050-0769	
		P70-312	050-0788	
		P50-312	050-0787	
		P30-312	050-0786	

L	R	Model	Part No.	Colour
		P90-13T	050-0770	
		P70-13T	050-0791	
		P50-13T	050-0790	
		P30-13T	050-0789	

L	R	Model	Part No.	Colour
		P90-R	050-0771	
		P70-R	050-0794	
		P50-R	050-0793	
		P30-R	050-0792	

L	R	Model	Part No.	Colour
		P90-RT	050-0772	
		P70-RT	050-0797	
		P50-RT	050-0796	
		P30-RT	050-0795	

L	R	Model	Part No.	Colour
		P90-RL	050-0773	
		P70-RL	050-0800	
		P50-RL	050-0799	
		P30-RL	050-0798	

L	R	Model	Part No.	Colour
		P90-R Fit	050-0870	
		P70-R Fit	050-0869	

Step 5(a): Choose receiver option (Sound Delivery System version 4.0 only)

Standard	Medium	Power	UltraPower – cShell only HF3 filter recommended
S 0-L S 0-R	M 0-L M 0-R	P 0-L P 0-R	UP 0-L UP 0-R
S 1-L S 1-R	M 1-L M 1-R	P 1-L P 1-R	UP 1-L UP 1-R
S 2-L S 2-R	M 2-L M 2-R	P 2-L P 2-R	UP 2-L UP 2-R
S 3-L S 3-R	M 3-L M 3-R	P 3-L P 3-R	UP 3-L UP 3-R

Step 5(b): Choose sensor receiver option for Audéo P-R Fit only (Sound Delivery System version 6.0 only)

Medium	Power – Sensor cShell only
Sensor Receiver 0L Sensor Receiver 0R	Sensor Receiver 0L Sensor Receiver 0R
Sensor Receiver 1L Sensor Receiver 1R	Sensor Receiver 1L Sensor Receiver 1R
Sensor Receiver 2L Sensor Receiver 2R	Sensor Receiver 2L Sensor Receiver 2R
Sensor Receiver 3L Sensor Receiver 3R	Sensor Receiver 3L Sensor Receiver 3R

Step 5(c): Choose ActiveVent Receiver option (compatible with Audéo P-R and Audéo P-RT only)

Audéo P-RL is compatible with ActiveVent Receiver, however due to the moving parts within ActiveVent and the increased potential for water intrusion into the receiver, it is not recommended for use with Audéo P-RL.

Moderate	(a) Choose earpieces – SlimTip Titanium AV (impression required)	(b) Choose Earpieces – Universal Silicone SlimTip AV (impression not required)
MAV Receiver 0L MAV Receiver 0R	L R SlimTip Titanium AV	Universal SlimTip AV A-L Universal SlimTip AV C-L
MAV Receiver 1L MAV Receiver 1R		Universal SlimTip AV A-R Universal SlimTip AV C-R
MAV Receiver 2L MAV Receiver 2R		Universal SlimTip AV B-L Universal SlimTip AV D-L
MAV Receiver 3L MAV Receiver 3R		Universal SlimTip AV B-R Universal SlimTip AV D-R

Step 6 (a): Choose dome coupling option (custom ear piece option below)

Dome style			Dome size		
Cap [#]	Left	Right			
Open [#]	Left	Right	Small	Medium	Large
Vented	Left	Right	Small	Medium	Large
Power	Left	Right	Small	Medium	Large
Vented Clear [#]	Left	Right	Small	Medium	Large



CapOpenVentedPowerVented Clear

Please note these are the only options available to order with Audéo P-R Fit

Step 6 (b): Choose custom ear piece option (impression required)

SlimTip (removeable earpiece)		cShell Acrylic		cShell Titanium**		Sensor cShell (062-0051)* (Only custom ear piece option for Audéo P-R Fit)		Impressions	
L	R	L	R	L	R	L	R		
Material		Shell colour		Shell colour		Shell colour		New impressions enclosed	
Hard Acrylic	Hard Acrylic								
Soft Silicone	Soft Silicone								
Shell colour		Faceplate colour		Faceplate colour		Faceplate colour		 	
 Clear (21)	 Clear (21)	 Beige/Pink (26)	 Beige/Pink (26)	 Pink (26)	 Pink (26)	 Pink (26)	 Pink (26)		
Rosa Tint (T)	Rosa Tint (T)	Cocoa (22)	Cocoa (22)	Cocoa (22)	Cocoa (22)	Cocoa (22)	Cocoa (22)		
Beige/Pink (26) (hard only)	Beige/Pink (26) (hard only)	Brown (28)	Brown (28)	Brown (28)	Brown (28)	Brown (28)	Brown (28)		
Tan (14) (hard only)	Tan (14) (hard only)	Transparent (13)	Transparent (13)	Transparent (13)	Transparent (13)	Transparent (13)	Transparent (13)		
		Tan (14)	Tan (14)	Tan (14)	Tan (14)	Tan (14)	Tan (14)	Okay to change to cShell if SlimTip not possible	
		Black Transp. (47)	Black Transp. (47)	Black Transp. (47)	Black Transp. (47)	Black Transp. (47)	Black Transp. (47)		

* Please note cShell Titanium products are not included free of charge with hearing aid order.
S = Standard, will default if no other selection ticked
^ Please note cShell Titanium is only available in Titanium shell colour and Black faceplate colour.
Please note Sensor cShell is only available to order with Audéo P-R Fit

Step 7: Choose venting, wax management and handling options

Venting (mm)*		Wax management		Handling	
L	R	L	R	L	R
	AOV Smart vent size calculation. Requires AC thresholds.		Ext. receiver tube + wax spring (WE)		Removal Line (RF)
	Occluded		None (X)		Removal Line + 5mm extension (RFE)
	1.0mm				Canal Lock (CL)
	1.1 - 1.6mm				Skeleton Lock (SL)
	1.7 - 2.2mm				
	2.3 - 3.0mm				
	Open				

*Cerustop default option for SlimTip HA and only option for SlimTip Titanium.

Please note: A coupling code representing actual vent size is printed on custom earpieces. Please enter this code into Phonak Target when programming your client's hearing aids.

**Not applicable for cShell Titanium.

Step 8: CROS Fitting?

CROS P transmitter required. Please complete the Phonak CROS P Order Form
Do not require a CROS P transmitter

Special requests:

Step 9: Choose accessories

Cleaning kits		Complimentary charging accessory (for Audéo P-R/P-RT/P-RL/P-R Fit)		Other accessories	
Choose one complimentary option		Phonak Charging Case Combi (075-0040-14)		TV Connector (076-3006-06)	
C&C Cleaning Kit (098-0041)		Phonak Charger BTE RIC (075-3016-10)		Remote Control (076-0065-P5)	
D-Dry+ Electronic Kit (098-0511-00014)		Phonak Charger Case Et Go (075-3017-10)		PartnerMic (076-4001-P510)	
		Phonak Life Charger (075-3023-10)		Roger Accessories - please complete the Phonak Roger Order Form	
		Additional charging accessory (for Audéo P-R/RT/P-R Fit)			
		Phonak Power Pack (071-0003)			

Additional information

HSP codes

Model	HSP codes	Category
P90-312	B679PHO	C2
P70-312	B680PHO	C2
P50-312	B681PHO	C2
P30-312	B682PHO	C2

Model	HSP codes	Category
P90-13T	B683PHO	C2
P70-13T	B684PHO	C2
P50-13T	B681PHO	C2
P30-13T	B686PHO	C2

Model	HSP codes	Category
P90-R	B687PHO	C2
P70-R	B688PHO	C2
P50-R	B689PHO	C2
P30-R	B690PHO	C2

Model	HSP codes	Category
P90-RT	B691PHO	C2
P70-RT	B692PHO	C2
P50-RT	B693PHO	C2
P30-RT	B694PHO	C2

Model	HSP codes	Category
P90-RL	B750PHO	C2
P70-RL	B751PHO	C2
P50-RL	B752PHO	C2
P30-RL	B753PHO	C2

Model	HSP codes	Category
P90-R Fit	B759PHO	C2
P70-R Fit	B760PHO	C2

Housing colours

P1 Sand Beige*	P3 Sandalwood	P4 Chestnut*	P5 Champagne
P6 Silver Gray*	P7 Graphite Gray	P8 Velvet Black	H0 Beige

* Only available for Audeo Fit

Models



Audéo P-R/P-RT/P-RL/P-R Fit charging accessories

 Phonak Charger Case Combi 075-0040-14 Compatible with P-R, P-RT & P-R Fit	 Phonak Power Pack 071-0003 Compatible with Charger Case Combi only	 Phonak Charger BTE RIC 075-3016-10 Compatible with P-R & P-RT	 Phonak Life Charger 075-3023-10 Compatible only with P-RL	 Phonak Charger Case & Go 075-3017-10
--	---	--	---	---

Other accessories

 TV Connector 076-3006-06	 RemoteControl 076-0065-P5	 PartnerMic™ 076-4001-P510
--	---	--