

Step 1: Order Details

Account Number: 20

Hearing Clinic:

Delivery Address:

Contact Name:Phone Number:

Email Address:

Date Ordered:Date Required:

Step 2: Patient Information

Last Name:

First Name:

Step 3: Which hearing aids are being used?

Please list hearing aid(s) and platform:

Step 4: Choose your Roger Microphone

Roger Microphones	Colour (refer to colour chart)	Qty
Roger for Adults (order 03 receiver below)		
Roger On V2		
Roger Select		
Roger Table Mic II		
Roger Table Mic II Bundle		
Roger Clip-On Mic		
Roger for Education (order 02 receiver below)		
Roger Pass-around ¹		
Roger Multimedia Hub ¹		
Roger Touchscreen Mic		

¹ Requires Touchscreen Mic

Roger iN Microphones <i>Includes 2 receiver installations</i>	Colour (refer to colour chart)	Qty
Roger On iN V2		
Roger Select iN		
Roger Table Mic II iN		
Roger Table Mic II iN Bundle		



Step 5: Choose your Roger Receiver (Not Required for Roger iN)

Roger Receivers (choose ONE) Roger Focus II also available. Please use the Roger Focus II Order Form.	Roger Receiver Options			
	Choose model (refer to list below)	Colour (refer to colour chart)	Version and Quantity	
			Education (02)	Adult (03)
Roger Design-integrated Receiver				
Roger X only				
Roger X with Audio Shoe				
Roger X with ComPilot				
Roger NeckLoop				
Roger Installer				



Design-integrated for implants



- Compatible with:

 - Sonnet
- Compatible with:

 - Nucleus 7
- Compatible with:

 - Naída CI Q

Design-integrated for Phonak hearing aids



- BTE & RIC (Size 675):

 - Belong
 - Venture
 - Vitus
- BTE & RIC (Size 13):

 - Belong
 - Venture
 - Vitus
- Compatible with:

 - Naída Q-RIC
 - Naída S-CRT
 - Bolero Q-M13
 - Cassia M H2O
- Compatible with:

 - Cassia microP
 - Cassia SP
 - Dalia microP
 - Dalia SP
- Compatible with:

 - Naída Q-SP
 - Naída S SP
 - Naída SP
- Compatible with:

 - Naída Q-UP
 - Naída S UP
 - Naída UP
 - Baseo Q15-UP

Information regarding RogerDirect:

- RogerDirect is a feature available in Paradise hearing aids that allows Roger receivers to be installed into the hearing aids.
- Roger iN microphones have two Roger receiver installations available
- It is possible to install existing Roger X receivers (if serial number is higher than 1944xxxx) using the Roger installer. This process will also allow connection to Roger devices that have already been purchased or outside of the Roger iN options (e.g. Roger Clip-On Mic).



Sand Beige	P1								
Amber Beige	P2								
Sandalwood	P3								
Chestnut	P4								
Champagne	P5								
Silver Gray	P6								
Graphite Gray	P7								
Velvet Black	P8								
Ruby	P9								
Petrol	Q1								
Electric Green	Q2								
Caribbean Pirate	Q3								
Dragon Orange	Q4								
Vanity Pink	Q5								
Blue Lagoon Transparent	Q6								
Royal Purple Transparent	Q7								
Precious Pink	T3								
Beige	01								
Brown	03								
Light Gray	04								
Black	06								
Pure Transparent	13								
Blue Transparent	37								
Purple Transparent	38								
Gray	23								
Light Pink	33								
Light Blue	34								
Taupe	35								
Chestnut Brown	50								
Palladium Black	70								
Chocolate	E5								
Cliff	E6								
Mid Gray	E8								
Glacier	E9								
Dark Cliff Gray	F0								
Sterling	F3								
Brown	L0								
White	L8								
Black	L9								
Sand	M1								
Charcoal	M2								
Brown	T1								
Grey	T2								
Alpine White	XN/T7								
Princess Pink	XP								
Lava Red	M6								
Blue Ocean	M7								
Majesty Purple	M8								
Anthracite	V1								
Beige	V2								
Black	V3								
Ebony	V4								
Nordic Grey	V5								
White	V6								

Roger On V2 Colour Chart:

P7

P5

Graphite Gray

Champagne

Roger Select Colour Chart:

P7

P5

W3

Graphite Gray

Champagne

Pearl White