

# Phonak CROS P Order Form

Phonak CROS P is a solution for unilateral hearing loss designed to help the wearer follow conversations. CROS P is worn on the unaidable ear to pick up sounds and wirelessly transmits them to a Phonak Audéo™ P hearing aid worn on the ear with better hearing. CROS P has no acoustical output.

CROS P includes two models and is compatible with Audéo P-13T, Audéo P-R and Audéo P-RT hearing aids in all performance levels.

## Step 1: Order details

Account Number: 20  
 Hearing Clinic:  
 Delivery Address:  
 Contact Name: Phone Number:  
 Email Address:  
 Date Ordered: Date Required:

## Step 2: Patient information

Name:

## Step 3: VIP service

VIP 24-hour service (\$50 ex GST)

## Step 4: Choose hearing instrument

Model	Part No.	HSP	Colour	Qty
CROS P-R	050-0860-xx	N0276CR		
CROS P 13	050-0861-xx	N0275CR		


## Colour chart



## Step 5: Choose your Acoustic Coupling Option

### CROS Slim Tube with Domes


CROS Tube is available in left/right and four different lengths (0, 1, 2, 3) with standard dome available in three sizes (S, M, L)



L R Select options in section a and b only below

### CROS Slim Tube with Slim Tip\*

CROS Tube is available in left/right and four different lengths (0, 1, 2, 3), with a custom made earpiece



L R Select options in section b and d only below

^ Impression required

### CROS Tip\*

Fully custom made earpiece for CROS, available in four different cable lengths (0, 1, 2, 3)



L R Select options in section c and d only below

^ Impression required

## Step 6: Choose from the following options

### (a) Domes 4.0

Dome style	Dome size	
	Left	Right
Cap	Left	Right
Open	Left	Right
Vented	Left	Right
Power	Left	Right


### (b) CROS Slim Tubes




### (c) CROS Tip Cables - must be paired with a CROS Tip below

054-0932	CROS Tip Cable 4.0 0L	054-0933	CROS Tip Cable 4.0 0R
054-0934	CROS Tip Cable 4.0 1L	054-0935	CROS Tip Cable 4.0 1R
054-0936	CROS Tip Cable 4.0 2L	054-0937	CROS Tip Cable 4.0 2R
054-0938	CROS Tip Cable 4.0 3L	054-0939	CROS Tip Cable 4.0 3R

### (d) Slim Tip 4.0 or CROS Tip



SlimTip  
L  
R



CROS Tip  
L  
R

	Shell colour	Clear (21)	Rose Tint (7)	Beige/Pink (26)	Tan (14)	Vent style	Cavity	Regular	Shell attribute	Hollow	Solid
SlimTip Acrylic (062-0006-01)											-
SlimTip Silicone (062-0007-01)											-
CROS Tip 4.0 Hard Acrylic (062-0016-**) - must be paired with a CROS Tip cable											-

S = Standard, will default if no other selection ticked

### Generic options (all earpieces)

Removal line Transparent (21) Removal line extension + 5mm  
 Locks Canal Skeleton

# Phonak CROS P Order Form



Phonak CROS P is a solution for unilateral hearing loss designed to help the wearer follow conversations. CROS P is worn on the unaidable ear to pick up sounds and wirelessly transmits them to a Phonak Audéo™ P hearing aid worn on the ear with better hearing. CROS P has no acoustical output.

CROS P includes two models and is compatible with Audéo P-13T, Audéo P-R and Audéo P-RT hearing aids in all performance levels.

---

## Step 7: Choose your wireless Paradise hearing instrument for the better ear

Audéo P-RT (NOTE: Only CROS P-R is compatible. Please also complete the Phonak Audéo P-RIC order form to specify correct inclusions).

Audéo P-R (NOTE: Only CROS P-R is compatible. Please also complete the Phonak Audéo P-RIC order form to specify correct inclusions).

Audéo P-13 (NOTE: Please also complete the Phonak Audéo P-RIC order form to specify correct inclusions).

---

## Step 8: Choose your accessories

### Wireless Communication Accessories

Please use the Phonak Audéo P-RIC order form to specify accessories.

### Roger Accessories

Refer to Phonak Roger Order Form

---

## Step 9: Check List

Phonak Audéo P RIC order form accompanies this order

Roger order form accompanies this order (if applicable)