

# Phonak Custom Earpiece 4.0 Order Form

Items on this order form are **only compatible with Marvel platform devices and newer**.  
For older platforms (e.g. Belong, Venture etc) please use the alternative earmould order form.

## Step 1: Order details

Account Number: 20

Hearing Clinic:

Delivery Address:

Contact Name: Phone Number:

Email Address:

Date Ordered: Date Required:

## Step 2: Patient information

Name:

Audiogram:

kHz	0.25	0.5	1	2	4
Left HL					
Right HL					

VIP 24 hour service (\$50 ex GST)  
VIP not possible for Titanium products

## Step 3: Product, style and options

	SlimTip 4.0		cShell 4.0		cShell 4.0 Titanium								
	L	R	L	R	L	R							
	Connection type		Shell colour		Vent style		Shell attribute						
	S Receiver	M Receiver	P Receiver	UP Receiver*	Clear (21)	Rosa Tint (7)	Beige/Pink (26)	Tan (14)	Cavity	Regular	Hollow	Solid	
SlimTip Acrylic (062-0006-01)				-	<input checked="" type="checkbox"/>				<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>	-	
SlimTip Silicone (062-0007-01)				-	<input checked="" type="checkbox"/>		-	-	-	<input checked="" type="checkbox"/>	-	<input checked="" type="checkbox"/>	
cShell Acrylic (062-0008-01)							<input checked="" type="checkbox"/>		-	<input checked="" type="checkbox"/>	-	-	
cShell Titanium (062-0009)*					Only available in Titanium Gray				-	<input checked="" type="checkbox"/>	-	-	

\*Please note SlimTip Titanium and cShell Titanium products are not included free of charge with hearing aid order. \*HF3 filter recommended

### SlimTip Acrylic options

Wax system  Cerustop (CS) Ext. receiver tube + wax spring (WE) None

### SlimTip Silicone options

Wax system  None

### cShell Titanium options

Cerustop (CS) HF 4 (H4) HF 3 (H3) ERS

### cShell Acrylic options (HF3 recommended for UP receivers)

Wax system Cerustop (CS) HF 4 (H4) HF 3 (H3) SmartGuard (SG) Ext. receiver tube + wax spring (WE)  
Faceplate colour  Pink (26) Cocoa (22) Brown (28) Transparent (13) Tan (14) Black Transparent (47)

### Generic options (all earpieces)

Tube length (0-3) Left: \_\_\_\_\_ Right: \_\_\_\_\_

Vent size AOV - let Phonak choose. Audiogram req'd. Occluded 1mm 1.1-1.6mm 1.7-2.2mm 2.3-3.0mm Open

Removal line Transparent (21) Removal line extension + 5mm

Locks Canal Skeleton

## Step 4: Preferences

If necessary, may we change the following:

- Step down power level - if audiogram permits (maintain 6dB Gain Reserve) Yes No
- cShell will be built with an anterior tube outlet in case of anatomical limitations. (See over for more details) Yes No

Special instructions:

# Phonak Custom Earpiece 4.0 Order Form

## Additional information

### Order details

- Please fill in client / order details
- Enter delivery date required
- Enter Audiogram (required for AOV)
- Select warranty and rush order details

### Product, Style and Options

- Select earpiece type
- Select receiver type
- Select colours
- Select tube length (mandatory for cShell)
- Select wax system\*
- Select vent
- Select removal line

\*All SlimTip receivers will be delivered with a CeruShield wax guard in the receiver spout. The SlimTip can be built with an additional wax protection.

### Preferences

- Let us know about special wishes
- Select if we may change receiver level or cShell build style under consideration of the audiogram, if necessary
- Select if BTE or xReceiver should be delivered together with the earpiece

### Alternative cShell build style

We offer an alternative cShell build style with an anterior tube outlet. This build style provides high cosmetic acceptance in spite of a slightly larger size and can be used in case of build problems due to anatomical limitations.



Fig.1 cShell - standard tube





Fig.2 cShell - anterior tube outlet


### Tube length

Measure the tube length to select the correct Receiver. Measure both ears for binaural fitting.




### Vent Styles

We offer two vent solutions for SlimTip.



Regular:

- Classic vent solution
- Fixed diameter (not enlargable)



Cavity:

- Optimised for best fitrate
- Easy to enlarge by the fitter

### Hearing Instrument Compatibility

The SDS 4.0 earpiece portfolio is compatible to Marvel hearing aids or newer. It is not backward compatible to Belong or earlier generation hearing aids.