

# Phonak Custom Earpiece 4.0 & SDS 5.0 Receiver Order Form

For Lumity, Terra+, Paradise and Marvel platform devices.

## Step 1: Order details

Account Number: 20

Hearing Clinic:

Delivery Address:

Contact Name: Phone Number:

Email Address:

Date Ordered: Date Required:

## Step 2: Patient information


Name:

Audiogram:

kHz	0.25	0.5	1	2	4
Left HL					
Right HL					

VIP 24 hour service (\$57.75 incl GST)  
VIP not possible for Titanium products

## Step 3: Product, style and options



Please note if a different configuration is required for each ear, please use two separate forms.

	Connection type				Shell colour				Vent style		Shell attribute		
	S Receiver	M Receiver	P Receiver	UP Receiver*	Clear (21)	Rosa Tint (7)	Beige/Pink (26)	Tan (14)	Cavity	Regular	Hollow	Solid	
SlimTip Acrylic (062-0006-01)				-	<input checked="" type="checkbox"/>				<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>	-	
SlimTip Silicone (062-0007-01)				-	<input checked="" type="checkbox"/>		-	-	-	<input checked="" type="checkbox"/>	-	<input checked="" type="checkbox"/>	
cShell Acrylic (062-0008-01)							<input checked="" type="checkbox"/>		-	<input checked="" type="checkbox"/>	-	-	
cShell Titanium (062-0009)*					Only available in Titanium Gray				-	<input checked="" type="checkbox"/>	-	-	

\*Please note SlimTip Titanium and cShell Titanium products are not included free of charge with hearing aid order. \*HF3 filter recommended

### SlimTip Acrylic options

Wax system  Cerustop (CS) Ext. receiver tube + wax spring (WE) None

### SlimTip Silicone options

Wax system  None

### cShell Titanium options

Cerustop (CS) HF 4 (H4) HF 3 (H3) ERS

### cShell Acrylic options (HF3 recommended for UP receivers)

Wax system Cerustop (CS) HF 4 (H4) HF 3 (H3) SmartGuard (SG) Ext. receiver tube + wax spring (WE)  
Faceplate colour  Pink (26) Cocoa (22) Brown (28) Transparent (13) Tan (14) Black Transparent (47)

### Generic options (all earpieces)

Tube length (0-3) Left: \_\_\_\_\_ Right: \_\_\_\_\_

Vent size AOV - let Phonak choose. Audiogram req'd. Occluded 1mm 1.1-1.6mm 1.7-2.2mm 2.3-3.0mm Open

Removal line Transparent (21) Removal line extension + 5mm

Locks Canal Skeleton

## Step 4: Preferences

If necessary, may we change the following:

- Step down power level - if audiogram permits (maintain 6dB Gain Reserve) Yes No
- cShell will be built with an anterior tube outlet in case of anatomical limitations. (See over for more details) Yes No

Special instructions:

# Phonak Custom Earpiece 4.0 Order Form

## Additional information

### Order details

- Please fill in client / order details
- Enter delivery date required
- Enter Audiogram (required for AOV)
- Select warranty and rush order details

### Product, Style and Options

- Select earpiece type
- Select receiver type
- Select colours
- Select tube length (mandatory for cShell)
- Select wax system\*
- Select vent
- Select removal line

\*All SlimTip receivers will be delivered with a CeruShield wax guard in the receiver spout. The SlimTip can be built with an additional wax protection.

### Preferences

- Let us know about special wishes
- Select if we may change receiver level or cShell build style under consideration of the audiogram, if necessary
- Select if BTE or xReceiver should be delivered together with the earpiece

### Alternative cShell build style

We offer an alternative cShell build style with an anterior tube outlet. This build style provides high cosmetic acceptance in spite of a slightly larger size and can be used in case of build problems due to anatomical limitations.



Fig.1 cShell - standard tube





Fig.2 cShell - anterior tube outlet


### Tube length

Measure the tube length to select the correct Receiver. Measure both ears for binaural fitting.




### Vent Styles

We offer two vent solutions for SlimTip.



Regular:

- Classic vent solution
- Fixed diameter (not enlargable)



Cavity:

- Optimised for best fitrate
- Easy to enlarge by the fitter

### Hearing Instrument Compatibility

The SDS 4.0 earpiece portfolio and SDS 5.0 receivers is compatible with Lumity, Paradise and Marvel hearing aids. It is not backward compatible to Belong or earlier generation hearing aids.