

# Phonak Earmould 4.0 Order Form

For Marvel platform BTE devices and newer.

## Step 1: Order details

Account Number: 20

Hearing Clinic: \_\_\_\_\_

Delivery Address: \_\_\_\_\_

Contact Name: \_\_\_\_\_ Phone Number: \_\_\_\_\_

Email Address: \_\_\_\_\_

Date Ordered: \_\_\_\_\_ Date Required: \_\_\_\_\_

## Step 2: Patient information

Name: \_\_\_\_\_

Audiogram: \_\_\_\_\_  
(Required for AOV vent option) VIP 24 hour service (\$50 ex GST)

	kHz	0.25	0.5	1	2	4
Left AC (dB HL)						
Right AC (dB HL)						

## Step 3: Earmould style

	Canal Mould (CU)	L	R
	Canal Lock (CL)	L	R
	Helix Lock (HL)	L	R
	Semi Skeleton (SS)	L	R
	Skeleton (SK)	L	R
	Carved Half Shell (HC)	L	R
	Carved Full Shell (SC)	L	R
	Ultra-carved Full Shell (SC)	L	R
	Full Shell (SU)	L	R
	CROS B (CB)	L	R

\* Compatible to tubing 13M and 13T

\*\* Compatible to tubing 13T and 13W (silicone only)

## Step 4: Options

Earmould material	Hard Acrylic	
	Silicone S40	
	Silicone S70	
Earmould colour (Hard Acrylic)	Refer to colour chart overleaf	
Earmould colour (Silicone)	Refer to colour chart overleaf	
Connector (optional)	Plastic Tube Lock (TRS)**	
Tubing	Medium Tube (13M) – Standard	
	Thick Wall Tube (13T)	
	Dry Wall Tube (13D)	
	Dry Thick Wall Tube (13W)	
SlimTube	SlimTube Cut (STC4)	
Vent size	AOV (enter audiogram)	
	Customer specific	Left: mm
		Right: mm
Canal length	Short (RS)	
	Medium (RM) – Standard	
	Long (RL)	
Additional options	Removal Line	
	Tapered Canal Tip	
	Extended Helix (Top Lock)	
	Concha Lock	

## Step 5: Preferences

Indicate the hearing instrument to be fitted with the earmould.

Model: \_\_\_\_\_

Colour: \_\_\_\_\_

## Special requests:

Call-back requested

# Phonak Earmould 4.0 Order Form

Additional information

## Acoustically Optimised Vent (AOV)


We recommend to select AOV. We will set the correct vent diameter and style consider audiogram / shell geometry / feedback canceler and required low frequency gain. The client's audiogram is mandatory for the calculation of the AOV.

## Canal Length

Canal length short – at 1st bend

Canal length medium – between 1st and 2nd bend

Canal length long – at 2nd bend



## Plastic Tube Lock (TRS)

For silicone ear moulds we offer a plastic tube lock. The plastic tube lock is for an adhesive free sound tube connection.

## Special Colours

We offer opaque colours upon request. Please provide a note in the special comment field.

## Preferences

Let us know if you have specific requests. Select if BTE should be delivered together with the earmould.

## Earmould Colour Chart

- |   |                                    |                                 |                                 |
|---|------------------------------------|---------------------------------|---------------------------------|
| <input type="checkbox"/> Transparent (Standard) |                                    |                                 |                                 |
| <input type="checkbox"/> White                  | <input type="checkbox"/> Rosa Tint | <input type="checkbox"/> Beige  | <input type="checkbox"/> Tan    |
| <input type="checkbox"/> Red                    | <input type="checkbox"/> Orange    | <input type="checkbox"/> Purple | <input type="checkbox"/> Yellow |
| <input type="checkbox"/> Blue                   | <input type="checkbox"/> Green     | <input type="checkbox"/> Black  | <input type="checkbox"/> Pink   |

## SlimTube Cut (STC4)

For small ears with narrow canals, we offer a Slimtube Cut.

SlimTube Cut (STC4) = SlimTube with removed connector

The SlimTube Cut is compatible to Marvel hearing aids or newer. It is not backward compatible to Belong or earlier generation hearing aids.

## Tube Length

Measure the tube length to select the correct SlimTube. Measure both ears for binaural fitting.

